


Instruction Manual for AD-1675 Tabletop Breeze Break for the AD-4212 Series

1. Features


- The tabletop breeze break protects the balance from breezes due to air conditioning and body motion, thereby reducing balance weighing error.
- The transparent panel assembly consists of antistatic plastic material, which protects the balance from static electricity.
- The hole on the back panel can be used to pass the RS-232C and AC adapter cables through.
- Sliding panels on top, right and left side enable samples to be inserted or removed from each side.

2. Specifications

- Weight : 2 kg
- Size : 237 (Width) × 371 (Depth) × 238 (Height) mm
- Material : Transparent panel: Antistatic plastic
Frame: Aluminum
- Applicable models (balances made by A&D)
 - AD-4212A series
 - AD-4212B series
 - AD-4212C series
 - AD-4212D series



3. External dimensions



Caution

- Do not use organic solvents to clean the tabletop breeze break as that may damage the transparency and antistatic properties of the panels.
- Use a soft lint free cloth dampened with warm water and a mild detergent when cleaning the tabletop breeze break.

A&D Company, Limited

3-23-14 Higashi-Ikebukuro, Toshima-ku, Tokyo 170-0013 JAPAN
Telephone: [81] (3) 5391-6132 Fax: [81] (3) 5391-6148