

Precision Health Scale

UC-324NFC

Instruction Manual

Contents

Before Use	
Precautions	1
Scale Components	2
Display Symbols	
Automatic Power-Off Function	3
Near Field Communication (NFC) Function	3
Installing (or Replacing) the Batteries	
Measuring Your Weight on a Carpet	
Operations	
Setting the Unit	5
Measuring Your Weight Using the NFC Function	6
References	
Maintaining the Scale	7
Before Requesting Service	
Caution on Using the NFC Function	
Specifications	
•	


Precautions

- Please read this booklet before using your scale.
- Place your scale on a solid level surface, such as tile, wood, or laminate flooring.
- Do not place your scale on a soft, wet or slippery floor.
 If the scale must be used on a carpet, be sure to attach the plastic add-on feet as shown on Page 4.
- Your scale contains delicate, high-precision parts. Avoid extremes in temperature, humidity, direct sunlight, shock, moisture and dust.
- Do not step on your scale with your feet wet. Be sure that your feet are dry to prevent slipping or serious injury.
- Step on the center portion of your scale gently.
- A caregiver should attend when a disabled person uses the scale.
- Consult with your physician before you start a diet or exercises.
- Do not disassemble or modify the scale.

Indications for Use

This scale is designed for personal use.

Scale Components



Display Symbols

Symbols	Name	Description
0	Ready/Stable symbol	This symbol is displayed when the measurement display is stable.
+	Battery indicator	The normal battery power indicator with charge remaining
	Low battery indicator	If this symbol is displayed, replace all four batteries with new ones.
	WEEE label	
CE	CE mark	

Automatic Power-Off Function

The scale has an "Automatic power-off" function.

- If "0.00 kg /0.0 lb" and the Ready/ Stable symbol are displayed for about forty-five seconds, the scale will automatically turn off.
- If the O symbol and the weight are displayed for about three minutes, the scale will automatically turn off.

Near Field Communication (NFC) Function

To transmit data to or from the scale, communication devices that comply with the communication specifications of the scale are required

Examples of communication devices that comply

- NFC-enabled cellular phone
- Personal computer with NFC reader/writer
- Access Point

Notes

■ For information on communication methods and applications for the communication devices, refer to the instruction manuals for the devices.


Installing (or Replacing) the Batteries

Caution

- Insert each battery with its positive (+) and negative (-) terminals properly aligned with the corresponding symbols in the battery compartment.
- Improper installation of the batteries may cause battery leakage.
- Replace all four batteries with new ones, when the **___** symbol is displayed.
- Do not mix new and used batteries. It may cause damage to the scale or shorten the battery life.
- The **___** symbol does not appear when the batteries are completely drained.
- The battery life is dependent upon the environmental temperature, battery type and usage.
- Remove the batteries from the battery compartment, when the scale is not to be used for a long time. They may leak and damage the scale.
- The batteries provided with the scale are for trial only and may have a short life.

Steps for Installing (or Replacing) Batteries

- 1. Open the battery compartment cover on the underside of the scale.
- 2. Install four new batteries with the positive (+) and negative (-) terminals properly aligned with the corresponding symbols in the battery compartment.
- 3. Close the battery compartment cover.



Measuring Your Weight on a Carpet

Place the scale on a solid level surface, such as tile, wood, or laminate flooring. If the scale must be used on a carpet, be sure to attach the plastic add-on feet.

Steps for Attaching the Add-on Feet

- Each add-on foot has two pins.
 Match those pins to the holes on the scale leg and insert firmly.
- 2. Attach add-on feet on all four legs.

Add-on foot Scale leg


Notes

- Please be careful not to fall when using the add-on feet.
- We recommend using this scale on a solid and level surface.

Setting the Unit

Setting the unit manually

- Confirm that the scale is turned off.
 Press the SET switch. The scale beeps and displays the unit currently set.
- 2. Press the SET switch to select the unit.
- 3. After three seconds of non-operation, the scale confirms the selected unit with two beeps, displays "End" and automatically turns off.




Setting the unit by the communication device

(Only for devices that rewriting of the unit information is enabled.)

- 1. Confirm that the scale is turned off.

 Hold the communication symbol of the NFC-enabled device close to the communication area ((2)) of the scale.
- 2. When the data transmission is complete, the scale displays "End" and automatically turns off.

Note




Measuring Your Weight Using the NFC Function

After measurement, the communication device receives the data.

- Confirm that the scale is turned off.
 Step on the scale gently and stand still during the measurement. The scale turns on, displays all of the display segments.
- Wait the finish of the measurement.
 The process of the measurement is displayed with segment "-".
- 3. Step off the scale when the stable symbol O is displayed with three beeps (Sound: pi pi pi). If you do not step off the scale, the scale beeps every two seconds (Sound: pi) and the scale turns off after two minutes.
- 4. The scale displays all of the display segments. Do not step on the scale.
- 5. The measurement result is displayed.
- 6. Point the communication symbol of the NFC-enabled device to the communication area () on the scale.
- 7. When the personal data transmission is complete, the scale displays "End".

Notes

- The NFC-enabled device may not be recognized when it is slanted or placed out of the communication area. If "{ / ;" or "{ /;" appears, start from step 5 again.
- After the measurement result is displayed about three minutes, the scale will automatically turn off and the result will be deleted.



Maintaining the Scale

- Clean the scale with a soft dry cloth or a cloth moistened with water and a mild detergent. Never use alcohol, benzene, thinner or other harsh chemicals to clean the scale.
- Remove the batteries from the battery compartment, when the scale is not going to be used for a long time.

Before Requesting Service

	Problems	Possible Reason / Corrective Action	
E	is displayed.		
£ -[]	is displayed	Exceeds the maximum capacity.	
Err	is displayed.	Exceeds the maximum capacity. Does not reach the minimum capacity.	
E - 10	is displayed.		
E -	is displayed.	The communication device may be slanted or placed out of the communication area of the scale.	
E - 12	is displayed.		
H	is displayed.	Fuer and the manufacture are ait.	
L	is displayed.	Exceeds the maximum capacity.	
	is displayed.		
Nothing is displayed by		Remove and reinstall the batteries.	
tapping on the scale or		Or, replace all four batteries with new ones.	
pressing the SET switch.			
		Place the scale on a solid level surface.	
The me	easurement	Remove any foreign objects that might be touching the	
results are too much		bottom of the scale.	
different from the		When the scale is placed on a carpet, be sure that the	
expected values.		add-on feet are attached or try moving the scale to a	
		level surface such as tile or wood.	

If you can not locate and fix the problem, contact your supplier, or the A&D service group. The A&D service group will support authorized suppliers with technical information, spare parts and units.

Caution on Using the NFC Function

The scale has an NFC Forum Type 3 Tag wireless interface module built in.

When the scale is used near wireless communication devices which use the same frequency as that of the scale (13.56 MHz) or use its harmonic frequencies, electromagnetic interference may occur. In that case, turn off the devices that are not in use, or place the devices 1 meter or more away from the scale.

Specifications

Type UC-324NFC

Maximum capacity 150 kg / 330 lb Minimum capacity 10 kg / 22 lb

Minimum display 0 to 100 kg: 0.05 kg, 100 to 150 kg: 0.1 kg

0 to 330 lb: 0.1 lb

Wireless communication NFC wireless interface module

Setting items Unit: kg or lb

Power source Four type AA batteries (R6P or LR6, 1.5 volt x 4)

Battery life Approximately 1 years (with two measurements a day

using Alkaline batteries)

Dimensions 350 (D) X 350 (W) X 39 (H) mm

Mass 2.3 kg including the batteries

Accessories Instruction manual, Batteries (4), Add-on feet (4)

Specifications are subject to change without prior notice or any obligation on the part of the manufacturer.

Compliance with European Directives

The scale complies with the statutory EMC (Electromagnetic Compatibility) directive 2004/108/EC.

Note: Indoor use only. The displayed value may be adversely affected under extreme electromagnetic influences.



A&D Company, Limited

1-243 Asahi, Kitamoto-shi, Saitama 364-8585 JAPAN Telephone: [81] (48) 593-1111 Fax: [81] (48) 593-1119

A&D INSTRUMENTS LTD.

Unit 24/26 Blacklands Way, Abingdon Business Park, Abingdon, Oxfordshire OX14 1DY United Kingdom

Telephone: [44] (1235) 550420 Fax: [44] (1235) 550485

A&D ENGINEERING, INC.

1756 Automation Parkway, San Jose, California 95131 U.S.A. Telephone: [1] (408) 263-5333 Fax: [1] (408)263-0119

A&D Australasia Pty Ltd.

32 Dew Street, Thebarton, South Australia 5031 AUSTRALIA Telephone: [61] (8) 8301-8100 Fax: [61] (8) 8352-7409

A&D KOREA Limited 한국에이.엔.디(주)

대한민국 서울시 영등포구 여의도동 36-2 맨하탄 B/D 8층

(8th Floor, Manhattan Bldg. 36-2 Yoido-dong, Youngdeungpo-ku, Seoul, KOREA)

전화: [82] (2) 780-4101 팩스: [82] (2) 782-4280