BA Series

Liquid Crystal Display

Quick Start Guide
BA Series
Analytical Balances

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A&D Company, Ltd.

© 2021 A&D Company, Limited. All rights reserved.

No part of this publication may be reproduced, transmitted, transcribed, or translated into any language in any form by any means without the written permission of A&D Company, Limited.

The contents of this publication and the specifications of the instrument covered by this publication are subject to change for improvement without notice.

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CT

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1. Quick Start Guide

Thank you for purchasing an A&D BA series analytical balance. This quick start guide describes how to install the BA series balance, and it covers the basic functions and main operations of the balance. For further information about the BA series balance, please refer to the separate instruction manuals listed in "1.1. Detailed Manuals".

1.1. Detailed Manuals

Detailed manuals are available for download from the A&D website (<https://www.aandd.jp>):

BA Series Instruction Manual

This is the manual to help you understand the functions and operations of the BA series in detail and make full use of them.

Ionizer (static eliminator) Instruction Manual

This is the manual describing the functions of the AX-ION-25 ionizer and instructions on use.


2. Assembly and Installation

Caution


- A BA series analytical balance is a precision instrument, and it should be carefully unpacked. It is advisable to store the packing materials so that they can be used when transporting the balance.
- The contents of the package vary depending on the product. Refer to the illustration of the packing contents and make sure that everything is included.
- Do not connect the AC adapter to the balance until the balance is assembled and installed.
- Use the dedicated AC adapter specified for the balance.
- Do not connect the included AC adapter to other devices.
- If you use the wrong AC adapter, the balance and other equipment may not work properly.
- Make sure that the AC adapter is unplugged before you connect the ionizer.
- The BA series analytical balance is composed of the weighing unit and display unit. It is not possible to replace just one of the units.

2.1. Unpacking


Packing contents common to BA series

Packing contents for BA-6E / BA-6DE

Packing contents for BA-225 / BA-225D / BA-125D


2.2. Assembly

The BA-6E is used here for the example of assembly. Other models can be assembled in the same way.

Step 1 Pull out the breeze break locking handles.

Step 2 Place the breeze break unit on the weighing unit.


Step 3 Push in the breeze break locking handles to secure the units to each other.

Step 4 Assemble the weighing pan.

For BA-6E / BA-6DE, assemble the five parts of the weighing pan in the correct position by referring to the cross-sectional view.

For BA-225 / BA-225D / BA-125D, assemble the four parts of the weighing pan in the correct position by referring to the cross-sectional view.


Cross-sectional view
BA-6E / BA-6DE


Cross-sectional view
BA-225 / BA-225D
BA-125D

Step 5 Assemble the glass breeze break.

Assemble the BA-6E / BA-6DE glass breeze break in the weighing chamber. Align the side panes with the notch of the base stand.

Example of the small breeze break with the right side open.

Step 6 Insert the breeze break top door from the handle side.

Step 7 Assemble the breeze break doors on the left and right sides.


7a Insert the breeze break doors on the left and right sides.

7b Secure the breeze break doors with the latches.



Step 8 Assemble the breeze break front glass.

8a Insert the breeze break front glass.


8b Secure the breeze break front glass with the latches.

- Step 9 Using the included cable A and cable B, connect the "weighing unit" and "breeze break unit" and the "display unit" and "weighing unit" respectively, paying attention to the direction of the arrow on each cable.
- Step 10 Connect the breeze break unit and ionizer with the included cable C.
- Caution** The ionizer operates when the AC adapter is connected to the balance.
- Step 11 Connect the included AC adapter to the balance.

BA-225 / BA-225D / BA-125D


BA-6E / BA-6DE

3. Installation Considerations, Weighing Preparation and Precautions


Prepare the following installation conditions in order to bring out the full performance of the balance.

- ❑ Install the balance in an environment where the temperature and humidity are not excessive. The best operating temperature is about $20^{\circ}\text{C} \pm 2^{\circ}\text{C}$ at about 45% to 60% RH relative humidity.
- ❑ Install the balance where it is free of dust.
- ❑ The weighing table should be solid. (An anti-vibration table or stone table is ideal)
- ❑ Place the balance on a horizontal table, and make sure that it is not tilted.
- ❑ Install the balance in a stable location, avoiding vibration and shock. Corners of rooms on the first floor are best, as they are less prone to vibration.
- ❑ Install the balance where it is not affected by heaters or air conditioners.
- ❑ Install the balance where it is not exposed to direct sunlight.
- ❑ Install the balance away from equipment which produces magnetic fields.
- ❑ Level the balance with the leveling feet and bubble spirit level. Refer to "3.1. How to adjust the level of the balance".
- ❑ Be sure to warm up the balance before use for at least an hour, or at least four hours for BA-6E / BA-6DE, with the AC adapter connected to the power supply.
- ❑ Adjust the sensitivity of the balance before using it for the first time or after having moved it to another location so that accurate weighing can be performed. For details, refer to the detailed instruction manual listed in "1.1. Detailed Manuals".

Caution Do not install the balance where flammable or corrosive gas is present.

3.1. How to adjust the level of the balance


- ❑ Level the balance by adjusting the leveling feet so that the bubble of the bubble spirit level is centered in the red circle.


When the bubble is off to the left. Turn the leveling foot on the front right in the clockwise direction.	When the bubble is off to the right. Turn the leveling foot on the front left in the clockwise direction.
When the bubble is off to the backward position. Turn both leveling feet on the front in the clockwise direction at the same time.	When the bubble is off to the forward position. Turn both leveling feet on the front in the counter clockwise direction at the same time.

4. Display and Key Panel (Basic Operation)


Lit display

Blinking display

Keys

Key operation

Key operations affect how the balance functions. Normal key operation during measurement is "Press and release the key immediately" or "Press the key". Please do not "Press and hold the key (for 2 seconds)" unless required.

Press the key
(Press and release
the key immediately.)

Press and hold
the key for 2
seconds.

Key	When pressed and released	When pressed and held (for 2 seconds)
	Turns the display on and off. When the display is turned off, only the standby indicator is displayed. When the display is turned on, weighing is possible. With the password function, "Enter password at the start of weighing" will be displayed when the display is turned on. For details, refer to the separate instruction manual listed in "1.1. Detailed Manuals". The ON:OFF key is active at any time, and pressing this key during operation always turns off the display.	Switches the IR sensor on and off.
	In weighing mode, press the button to turn the digit for readability on and off. In counting or percent mode, press the button to enter the sample storing mode.	□ Enters the function table mode. Please refer to the separate instruction manual listed in "1.1. Detailed Manuals". □ Runs the repeatability check function when pressed and held for another 2 seconds after the function table menu is displayed. For details, refer to the separate instruction manual listed in "1.1. Detailed Manuals".
	Switches the weighing units stored in the function table.	Enters "Changing Weighing Speed" mode.
	Performs sensitivity adjustment of the balance using the internal mass.	Displays other items of the sensitivity adjustment menu.
	Stores or outputs the weighing value when stable according to the function table settings. (Data is output at factory setting.)	Enters mode to change the unit mass registration number in counting mode. The following settings are possible with the function table: □ Outputs "Title block" and "End block" for GLP/GMP report. □ Displays the data memory menu.
	Performs tare operation.	-
	Sets the displayed value to zero.	-
	IR sensor (touchless sensor). It reacts when you bring your hand close to it. Opening and closing of the breeze break door are assigned. For details, refer to "5. IR Sensor and Auto Door".	-

5. IR Sensors and Automatic Doors

5.1. IR sensors

BA series analytical balances are equipped with IR sensors that allow operation without directly touching the balance display. At factory setting, the IR sensors on the left and right of the display are assigned to open and close the breeze break doors. By pressing and holding the **ON:OFF** key (for about 2 seconds), you can switch the IR sensors on and off with a single touch.

5.2. Automatic doors

BA series analytical balances are equipped with automatic doors that allow you to open/close the breeze break without touching the doors.

The IR sensors on the left and right of the display are assigned settings to open/close the breeze break doors on the right and left respectively.

At factory setting, each breeze break door opens to the position it was previously opened to.

It is also possible to change which sensors the breeze break doors are assigned to and for the doors to be fully open or partially open in the function table of the balance.

For details, refer to the separate instruction manual listed in “1.1. Detailed Manuals”.

Opening the breeze break door

Step 1. If you want to open the breeze break door on the left side (or on the right side) when it is closed, hold your hand over the right IR sensor (or the left IR sensor).

Step 2. The detection buzzer sounds and the breeze break door on the left (or right) side opens.

Closing the breeze break door


Step 1. If you want to close the breeze break door on the left side (or on the right side) when it is open, hold your hand over the right IR sensor (or the left IR sensor).

Step 2. The detection buzzer sounds and the breeze break door on the left (or right) side closes.

Operation overview

The left IR sensor opens/closes the door on the right.

The right IR sensor opens/closes the door on the left.

Example of BA-6E / BA-6DE

6. Sensitivity Adjustment

6.1. Automatic sensitivity adjustment

The sensitivity of the balance is automatically adjusted using the internal mass according to ambient temperature change, set time, or interval time. For details, refer to the separate instruction manual listed in "1.1. Detailed manual".

6.2. Sensitivity adjustment using the internal mass

Sensitivity adjustment using the internal mass can be performed with a single touch.





- Step 1. With nothing on the weighing pan, plug in the balance with the AC adapter in order to warm up the balance for at least 1 hour (or 4 hours for BA-6E / BA-6DE).
- Step 2. Press the **CAL** key to display **[AL in]**.
- Step 3. The balance performs sensitivity adjustment using the internal mass. Do not allow vibration or drafts to affect the balance.
- Step 4. After sensitivity adjustment, the "sensitivity adjustment record" is output and stored as set in the function table.
For details, refer to the separate instruction manual listed in "1.1. Detailed manual".
- Step 5. When finished, the balance automatically returns to weighing mode.

6.3. Sensitivity adjustment using an external weight

Using your external weight, adjust the sensitivity of the balance.

For details, refer to the separate instruction manual listed in "1.1. Detailed manual".

7. Error Displays (Error Codes)

Display	Code	Description
	EC,E11	Stability error Weighing value is unstable and "zero display", "sensitivity adjustment", etc. cannot be executed. Check around the pan. Refer to the separate instruction manual listed in "1.1. Detailed Manuals". Improve the environment of the installation location (vibration, drafts, static electricity, etc.). To return to weighing mode, press the  key.
		Out of the setting range The value to be set exceeds the setting range. Set again within the setting range.
		Malfunction of the internal memory element of the balance If this error continues to be displayed, repair is necessary. Please contact your local dealer for repair.
	EC,E16	Internal mass error Applying the internal mass does not yield a change in the mass value as specified. Confirm that there is nothing on the pan and perform the operation from the beginning.
	EC,E17	Internal mass error The internal mass application mechanism does not function properly. Perform the operation from the beginning.
 		Abnormality in the internal memory data of the balance If this error continues to be displayed, repair is necessary. Please contact your local dealer for repair.
		Weighing unit error Data from the weighing unit has not been received. Check the connection between the weighing unit and the display unit, and turn the power off and then on again.
		Breeze break unit error Communication with the breeze break is not connected. Reconnect the cable between the breeze break unit and the weighing unit.
		Serial number error The serial numbers of the weighing unit and the display unit do not match. Reconnect with the correct combination.
	EC,E20	Sensitivity adjustment weight error (Positive value) The sensitivity adjustment weight is too heavy. Check around the pan. Check the sensitivity adjustment mass value. To return to weighing mode, press the  key.

Display	Code	Description
	EC,E21	Sensitivity adjustment weight error (Negative value) The sensitivity adjustment weight is too light. Check around the pan. Check the sensitivity adjustment mass value. To return to weighing mode, press the key.
		Overload error A sample beyond the balance weighing capacity has been placed on the pan. Remove the object from the pan.
		Weighing pan error The mass value is too light. Check that the weighing pan is installed correctly. Set the weighing pan correctly. Calibrate the balance.
		Sample mass error The sample is too light to be stored as a sample mass for the counting mode or percent mode. The sample cannot be used.
 		Unit mass error The sample unit mass for the counting mode is too light. Storing and using it for counting may cause a counting error. Add samples to reach the specified number and press the key. Pressing the key without adding samples will put the balance in counting mode. But, for accurate counting, be sure to add samples.
		Clock battery error The clock backup battery has been depleted. Press any key and set the time and date. Even if the clock backup battery is depleted, the clock and calendar function works normally as long as the balance is powered with the AC adapter. If this error appears frequently, contact your local dealer.
		Power supply voltage fault The voltage supplied from the AC adapter is abnormal. Please check if the problem is the AC adapter supplied with the balance.
 		Repeatability error The standard deviation (SD) of repeatability exceeded 50 digits. Review the installation environment of the balance. in repeatability display. in minimum weighing value (reference value) display. For details, refer to the separate instruction manual listed in "1.1. Detailed Manuals".
 Blinking (alternately)		Full memory The number of stored weighing values has reached the upper limit. In order to store new weighing values, it is necessary to delete data. For details, refer to the separate instruction manual "1.1. Detailed Manuals".
 Blinking (alternately)		Full memory The stored sensitivity adjustment / calibration test history has reached 50 results. In order to store a new result, the oldest history will be deleted. For details, refer to the separate instruction manual "1.1. Detailed Manuals".

Display	Code	Description
	EC,E00	Communications error A protocol error occurred in communications. Check the format, baud rate, and parity.
	EC,E01	Undefined command error An undefined command was found. Check the transmitted command.
	EC,E02	Not ready The received command cannot be executed: (e.g.) Q command was received when not in the weighing mode. (e.g.) Q command was received while rezeroing. Adjust the delay time to transmit a command.
	EC,E03	Timeout error With the timeout parameter set to t-UP 1 , there was a wait time of approximately 1 second or more while receiving command characters. Check the communication.
	EC,E04	Character length error The number of characters in the received command has exceeded the limit. Check the command to transmit.
	EC,E06	Format error The description of the received command is incorrect: (e.g.) The number of digits of numerical values is incorrect. (e.g.) There are alphabet characters among the numerical values. Check the transmitted command.
	EC,E07	Parameter setting error The value of the received command has exceeded the allowed value. Check the setting range of the numerical value of the command.

7.1. Other display

This is the automatic sensitivity adjustment notice (the  mark blinking).

If the balance is not used for several minutes with this mark blinking, the balance automatically performs sensitivity adjustment using the internal mass. (The blinking period depends on the operating environment.)

Advice


Although it is possible to continue using the balance even while this mark is blinking, use after sensitivity adjustment is completed is advisable in order to maintain the weighing accuracy.

8. Ionizer (AX-ION-25)

The AX-ION-25 ionizer comes standard with the BA-6E / BA-6DE. The ionizer removes static electricity by irradiating the target object with positive or negative ions, which are generated from four discharge electrodes by DC corona discharge. Ordinarily, insulators such as powders, filters, and weighing paper tend to be charged when the humidity is 45% RH or less, and an error of a few milligrams may occur during weighing. By removing static from the weighing object with the ionizer, it is possible to eliminate the error in the weighing value due to charging, and correct weighing can be performed. For details, refer to the separate instruction manual in "1.1. Detailed Manuals".


8.1. Using the ionizer

- Step 1. Make sure that the AC adapter of the balance is not connected, and then connect the balance and the ionizer (as noted in "2. Assembly and Installation").
- Step 2. Connect the AC adapter, plug it in, and then turn on the power. The power lamp of the ionizer lights up.
- Step 3. Place the target object within the effective range of static elimination.
- Step 4. Hold your hand over the IR sensor (touchless infrared proximity sensor) on the front of the ionizer to start static elimination. The ION lamp indicating that static elimination is in progress lights up.
- Step 5. Refer to the figure below for the effective range of static elimination, and perform static elimination. At factory setting, when 3 seconds have passed from the start of static elimination, the ION lamp indicating that static elimination is in progress turns off and static elimination stops.


Effective range of static elimination

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1756 Automation Parkway, San Jose, California 95131, U.S.A.
A&D AUSTRALASIA PTY LTD
A&D KOREA Limited