HR-AZ / HR-A

HR-250AZ
HR-251AZ
HR-150AZ
HR-100AZ

HR-250A
HR-251A
HR-150A
HR-100A

Quick Start Guide
Analytical Balance

English (EN)

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1WMPD4004338

1. Introduction

This quick start guide describes the basic functions and operations of the balance (HR-AZ / HR-A series). The functions and operating procedures of the balance are explained in details in the "**HR-AZ / HR-A series Analytical Balance INSTRUCTION MANUAL**". This manual can be downloaded from our website <https://www.aandd.jp>.


© 2021 A&D Company, Limited. All rights reserved.

- No part of this publication may be reproduced, transmitted, transcribed, or translated into any language in any form by any means without the written permission of A&D Company, Limited.
- The contents of this manual and the specifications of the instrument covered by this manual are subject to change for improvement without notice.
- Other trademarks and trade names are those of their respective owners.

2. Installation and Precautions

The balance is a precision instrument. Unpack it carefully. It is recommended to keep the packing materials to be used for transporting the balance in the future.

The packing contents differ depending on the balance model. See the following illustrations to confirm that everything is included.

Cautions:



- Verify that the AC adapter type is correct for your local voltage and receptacle type.
- Use the dedicated AC adapter specified for the balance.
- Do not use the AC adapter provided with the balance for other instruments.
- If you use a wrong AC adapter, the balance and other equipment may not operate properly.

2.1. Precautions for installing and weighing


⚠ Caution Do not install the balance where flammable or corrosive gas is present.

To obtain the optimum performance from the balance and acquire accurate weighing data, note the following:

- ❑ Install the balance in an environment where the temperature and humidity are not excessive. The best operating temperature is about $20^{\circ}\text{C} \pm 2^{\circ}\text{C}$ at about 45% to 60% RH relative humidity.
- ❑ Install the balance in a location that it is free of dust.
- ❑ The weighing table should be solid. (An anti-vibration table or stone table is ideal)
- ❑ Place the balance on a horizontal table, and make sure that it is not tilted.
- ❑ Install the balance in a stable place avoiding vibration and shock. Corners of rooms on the first floor are best, as they are less prone to vibration.
- ❑ Install the balance where it is not affected by heaters or air conditioners.
- ❑ Install the balance where it is not exposed to direct sunlight.
- ❑ Install the balance away from equipment that produces magnetic fields.
- ❑ Be sure to connect the AC adapter to the balance and the power supply before use for at least an hour, warm up the balance.
- ❑ Perform the sensitivity adjustment before using it for the first time or after having moved it to another location so that accurate weighing can be performed.
- ❑ Level the balance by adjusting the leveling feet so that the bubble of the bubble spirit level is centered in the red circle.


Bubble spirit level


When the bubble is off center to the left

Turn the leveling foot on the front right in the clockwise direction.


When the bubble is off center to the right

Turn the leveling foot on the front left in the clockwise direction.


When the bubble is off center to the rear

Turn both leveling feet on the front in the clockwise direction at the same time.

When the bubble is off center to the front


Turn both leveling feet on the front in the counter clockwise direction at the same time.

3. Display symbols and key operations

Display symbols

- Number of statistical data (Statistical calculation mode)
- The ratio between the weighing value and the weighing capacity, in a percentage, in the weighing mode. (Capacity indicator)
- Parameters of function.

Key Operations

The behavior of the balance differs depending on whether you "Press and release the key immediately" or "Press and hold the key (for 2 seconds)".

The normal weighing operation is to "Press and release the key immediately". Do not hold down the key unless you need to.

Press the key (Press and release the key immediately.)

Press and hold the key (for approx. 2 seconds)

Key	When pressed	When pressed and held (for 2 seconds)
	<ul style="list-style-type: none"> □ The balance is able to be turned on or off. The standby indicator is displayed when the display is turned off. The weighing mode is enabled when the display is turned on. A password is required when the balance is turned on using the password function. □ This key is always available. If this key is pressed during operation, the display will be turned off. 	
	<ul style="list-style-type: none"> □ In the weighing mode, readability changes. □ In the counting or percent mode, the balance enters the storing mode. 	The balance enters the function table mode.
	<ul style="list-style-type: none"> □ The weighing units selected on the function table are switched in order. 	The response adjustment mode is performed.

Key	When pressed	When pressed and held (for 2 seconds)
	<ul style="list-style-type: none"> <input type="checkbox"/> The HR-AZ series can perform sensitivity adjustment using the built-in weight. <input type="checkbox"/> In the function table mode, it works as cancel key. 	The sensitivity adjustment mode is performed.
	<ul style="list-style-type: none"> <input type="checkbox"/> The stable weighing data is outputted according to the function table. <input type="checkbox"/> In each operation procedure, it works as determination key. 	<ul style="list-style-type: none"> <input type="checkbox"/> No function at factory settings. <input type="checkbox"/> "Title block" and "End block" is outputted according to GLP/GMP settings on the function table.
	The current measurement value is regarded as zero and zero is displayed.	

4. Sensitivity adjustment

4.1. One-touch sensitivity adjustment using the internal mass for HR-AZ series

This function performs the sensitivity adjustment of the HR-AZ series using the built-in weight (the internal mass).

Procedure

- Step 1. Connect the AC adapter to the balance and warm up the balance for at least one hour with nothing on the weighing pan.
- Step 2. Place the accessory breeze break on the balance.
Press the **CAL** key to display **CAL , n**.
- Step 3. The sensitivity adjustment is performed using the internal mass.
Do not allow vibration or drafts to affect the balance.
- Step 4. If you choose the parameter to use the GLP output of the function table, a sensitivity adjustment report is outputted after the sensitivity adjustment.
- Step 5. The balance will automatically return to weighing mode after sensitivity adjustment.

4.2. Sensitivity adjustment using an external weight

This function performs the sensitivity adjustment using an external weight.

Procedure

- Step 1. Connect the AC adapter to the balance and warm up the it for at least one hour with nothing on the weighing pan.

Step 2. Press and hold the **CAL** key until **CALout** is displayed for approximately two seconds. Then release the key.

Step 3. The balance displays **CAL 0**.

- If you want to change the sensitivity adjustment weight, press the **SAMPLE** key and proceed to step 4. Refer to "8. CALIBRATION" in the instruction manual for sensitivity adjustment weights that can be selected.
- If you use the sensitivity adjustment weight value stored in the balance, proceed to step 5.

Step 4. Specify the weight value as follows :

SAMPLE keyAll segments blink : Weight selection mode.
The last three digits blink : Weight value adjustment.

RE-ZERO key (+) All segments blink : Select weight.

MODE key (-) The last three digits blink : Weight value adjustment.

PRINT keyWeight value is stored in memory. Data is maintained even if removing power supply.

CAL keyOperation is canceled and display returns to **CAL 0**.

Step 5. Confirm that there is nothing on the pan and press the **PRINT** key. The balance measures the **zero point**. Do not allow vibration or drafts to affect the balance. The balance displays the sensitivity adjustment weight value.


Step 6. Place the displayed **sensitivity adjustment weight** on the pan and press the **PRINT** key. The balance measures the sensitivity adjustment weight. Do not allow vibration or drafts to affect the balance.

Step 7. Remove the weight from the pan after displaying **End**.










Step 8. If the GLP output is set, "Calibration Report" is outputted.



Step 9. The balance will automatically return to the weighing mode.

Step 10. Place the sensitivity adjustment weight on the pan and confirm that the measurement value is within ± 2 digits. If it is out of that range, check the ambient conditions and repeat from steps 2.

5. Error codes

Display	Error code	Check and solutions
	EC,E11	<p>Stability error "Zero display" or "sensitivity adjustment" cannot be performed due to unstable weighing value. Check the area around the weighing pan. Refer to "2.1. Precautions for installing and weighing". Prevent vibration, drafts, static electricity and magnetic fields from influencing the balance. Press the CAL key to return to weighing mode.</p>
		<p>Out of range error The input value is out of range. Input a value within the range.</p>
	EC,E17	<p>Internal mass error (Only for the HR-AZ series) The mechanism of the internal mass did not work properly. Perform the operation again.</p>
	EC,E20	<p>Sensitivity adjustment weight error (Positive value) The sensitivity adjustment weight is too heavy. Confirm that the weighing pan is properly installed. Confirm the sensitivity adjustment weight value. Press the CAL key to return to the weighing mode.</p>
	EC,E21	<p>Sensitivity adjustment weight error (Negative value) The sensitivity adjustment weight is too light. Confirm that the weighing pan is properly installed. Confirm the sensitivity adjustment weight value. Press the CAL key to return to the weighing mode.</p>
		<p>Overload error The weighing value is beyond the weighing capacity. Remove the sample from the pan.</p>
		<p>Weighing pan error The weight value is too light. Confirm that the weighing pan and pan support are properly installed. Press the ON:OFF key two times to return to the weighing mode. If the error still persists, perform the sensitivity adjustment.</p>
		<p>Sample mass error The sample unit mass value for the counting mode or the percent mode cannot be stored because of too light. Use a sample that is heavier.</p>
  		<p>Unit mass error The sample unit mass value for the counting mode is too light. If unit mass value is saved as is, it may cause a counting error. Add samples of the specified number and press the PRINT key to acquire accurate weighing.</p>
		<p>Balance internal error If this error appears persistently, contact the local A&D dealer.</p>

Display	Error code	Check and solutions
		<p>Battery error of clock The back-up battery power for the clock has been completely depleted. Adjust the date and time after pressing one of the keys. If connecting the balance to power supply, even if the back-up battery power for the clock has been depleted the balance performance will be maintained. When this error appears frequently, repair is required.</p>
		<p>Low battery error Voltage of the power supply has dropped. Confirm whether the correct AC adapter is used. If a rechargeable battery is used, replace it.</p>
	EC,E00	<p>Communications error A protocol error occurred in communications. Confirm the format, baud rate and parity.</p>
	EC,E01	<p>Undefined command error An undefined command was received. Confirm the command.</p>
	EC,E02	<p>Not ready A received command cannot be processed. Check the timing of the transmitted command. e.g. The balance received a Q command in other than the weighing mode. e.g. The balance received a Q command while processing a RE-ZERO command.</p>
	EC,E03	<p>Timeout error When  is set, the balance did not receive the next character of a command within the time limit of one second. Confirm the communication.</p>
	EC,E04	<p>Excess characters error The balance received excessive characters in a command. Confirm the command.</p>
	EC,E06	<p>Format error A command includes incorrect data. Confirm the command. e.g. The numeric data is incorrect. e.g. The alphanumeric data is incorrect.</p>
	EC,E07	<p>Parameter setting error The received data exceeds the acceptable range. Confirm the parameter range of the command.</p>
Other errors		If the errors described above cannot be resolved or other errors are displayed, contact your local A&D dealer.

6. Specifications

	HR-100AZ HR-100A	HR-150AZ HR-150A	HR-250AZ HR-250A	HR-251AZ HR-251A	
Weighing capacity	102 g	152 g	252 g		
Maximum display	102.0084 g	152.0084 g	252.0084 g	252.008 g	
				62.0009 g *1	
Readability	0.1 mg			1 g	
				0.1 g	
Repeatability (Standard deviation)	0.1 mg	0.2 mg (200 to 250 g) 0.1 mg (0 to 200 g)		0.5 mg	
				0.1 mg	
Linearity	±0.2 mg	±0.3 mg		±1 mg	
				±0.3 mg	

*1 The balance allows weighing using a precision range, even with a heavy tare placed on the pan. (Smart range function)

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A&D AUSTRALASIA PTY LTD
A&D KOREA Limited
A&D Instruments India Private Limited