BH Series

(This document is included in the packing contents for the balance.)

Quick Start Guide
BH Series
Analytical Balances

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A&D Company, Ltd.

Warning Definition

The warnings described in this manual have the following meanings:

A potentially hazardous situation which, if not avoided, may result in minor or moderate injury or damage to the instrument.

Caution

Cautions to use the device correctly.

Note

Information or cautions to use the device correctly.

About This Manual

- (1) No part of this manual may be reprinted, copied, modified, or translated to another language without the prior written consent of A&D Company, Limited (A&D).
- (2) The contents of this manual are subject to change without notice.
- (3) Please contact A&D if you notice any uncertainty, errors, omissions, etc. in this manual.
- (4) A&D bears no liability for any loss or lost profits due to the operation of this product, and for direct, indirect, special, or consequential damages resulting from any defect in this product or this manual, even if advised of the possibility of such damage. Furthermore, A&D assumes no liability for claims of rights from third parties. Concurrently, A&D assumes no liability whatsoever for software or data losses.

© 2025 A&D Company, Limited

- ☐ Microsoft®, Windows®, Word®, and Excel® are trademarks of the Microsoft group of companies.
- ☐ The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by A&D is under license.
- ☐ iOS is the name of the operating system of Apple Inc. iOS is a trademark or registered trademark of Cisco in the U.S. and other countries and is used under license.
- ☐ Apple, the Apple logo and iPhone are trademarks of Apple Inc.
- ☐ App Store is a service mark of Apple Inc.
- ☐ Android™, Google Play and the Google Play logo are trademarks of Google LLC.
- ☐ Other product names and company names mentioned in this manual are trademarks or registered trademarks of their respective companies in Japan or other countries and regions.

1. Quick Start Guide

Thank you for purchasing an A&D BH series analytical balance. This quick start guide provides instructions on how to install the BH series balance and covers its basic functions and main operations. For details, refer to the separate instruction manual listed in "1.1. Detailed Manual".

1.1. Detailed Manual

The detailed manual is available for download from the A&D website (<https://www.aandd.jp>):

BH Series Instruction Manual

This manual provides detailed information on the functions and operations of the BH series, enabling you to fully utilize its features.

Scan the QR code to access the manual.

2. Assembly and Installation

CAUTION


- ☐ The BH series analytical balance is a precision instrument and should be carefully unpacked. It is advisable to store the packing materials for future transport of the balance.
- ☐ The contents of the package vary depending on the product. Refer to "2.1. Unpacking" to ensure all items are included.

CAUTION

- ☐ Confirm that the AC adapter type is correct for your local voltage and receptacle type.
- ☐ Use only the dedicated AC adapter specified for the balance.
- ☐ Do not connect the included AC adapter to other devices.
- ☐ Using the wrong AC adapter may cause the balance and other equipment to malfunction.

2.1. Unpacking

Packing contents and parts names for BH series:

- | | |
|--|---|
| 1 Main unit | 16 Cleaning brushes (large, small) |
| 2 Weighing pan | 17 AD-1689 tweezers (for 0.01 mg models* ¹ only) |
| 3 Low breeze break ring | 18 Spatula (for 0.01 mg models* ¹ only) |
| 4 High breeze break ring (for 0.01 mg models* ¹ only) | 19 Weigh boat (10 pcs) |
| 5 Breeze break bottom plate | 20 Pan support boss |
| 6 Display protection cover (PET plastic)* ² | 21 Backlit LCD |
| 7 Breeze break front pane | 22 Bubble spirit level |
| 8 Breeze break top pane | 23 Key |
| 9 Breeze break left side pane | 24 Leveling foot |
| 10 Breeze break right side pane | 25 Serial number |
| 11 Breeze break rear pane | 26 USB Interface (Type C)* ³ |
| 12 Quick Start Guide | 27 RS-232C interface |
| 13 AC adapter* ⁴ | 28 AC adapter input jack |
| 14 AC adapter ID labels | 29 Grounding terminal |
| 15 USB cable (approx. 2 m) | 30 External input terminal |

*¹ BH-225 / BH-225D

*² Attached to the main unit.



*³ For communication only.




*⁴ Accessories vary depending on the destination region.

2.2. Assembly

⚠ CAUTION

- ❑ Perform the following steps with the AC adapter removed from the main unit.

Step	Description	Parts diagram
1	Attach the breeze break bottom plate (5), breeze break ring (3 or 4), and weighing pan (2) to the main unit (1).	
2	Insert the breeze break rear pane (11) into the lower back groove of the main unit (1). Then, push the top of the pane into the latches until it clicks into place.	


Step	Description	Parts diagram
3	Insert the breeze break top pane (8) into the front groove of the breeze break frame from the front of the main unit (1).	 The diagram shows a perspective view of the main unit (1) with the breeze break frame. A top pane (8) is shown being inserted into the front groove of the frame. Arrows indicate the direction of insertion.
4	Insert the breeze break front pane (7) into the lower front groove of the main unit (1). Then, push the top of the pane into the latches until it clicks into place.	 The diagram shows a perspective view of the main unit (1) with the breeze break frame. A front pane (7) is shown being inserted into the lower front groove of the frame. Arrows indicate the direction of insertion and the latching mechanism.
5	For the breeze break side panes (9, 10), push the latches on the breeze break frame at the back of the main unit (1) upwards until they click. Insert the panes into the grooves of the breeze break frame from the back, ensuring the handles face outward. After inserting the breeze break side panes, return the latches to their original position.	 The diagram shows a perspective view of the main unit (1) with the breeze break frame. A side pane (9) is shown being inserted into the side groove of the frame. A side pane (10) is shown being inserted into the side groove of the frame. Arrows indicate the direction of insertion and the latching mechanism. A callout shows a latch being pushed upwards.

Step	Description	Parts diagram
6	Attach the AC adapter ID labels (14) to the AC adapters (13). CAUTION ❑ Ensure the AC adapter ID labels are attached to avoid using the wrong AC adapters.	
7	Insert the AC adapter (13) with the attached ID labels (14) into the AC adapter input jack (26) on the back of the main unit (1). Then, plug the other end into an outlet. CAUTION ❑ Be sure to warm up the balance for at least an hour before use.	

3. Installation Considerations, Preparation and Precautions

Ensure the following installation conditions to achieve optimal performance of the balance:


- ☐ The ideal operating temperature is approximately $20^{\circ}\text{C} \pm 2^{\circ}\text{C}$ with a relative humidity of 45% to 60%.
- ☐ Install the balance in a dust-free environment.
- ☐ The weighing table should be solid. An anti-vibration table or stone table is ideal.
- ☐ Place the balance on a horizontal surface, ensuring it is not tilted.
- ☐ Install the balance in a stable location, avoiding areas with vibration and shock. Corners of rooms on the ground floor are preferable as they are less prone to vibration.
- ☐ Avoid installing the balance near equipment such as heaters or air conditioners.
- ☐ Ensure the balance is not exposed to direct sunlight.
- ☐ Keep the balance away from equipment that produces magnetic fields.
- ☐ Level the balance using the leveling feet and bubble spirit level.
- ☐ Warm up the balance for at least an hour before use, with the AC adapter connected to the power supply.
- ☐ Perform a sensitivity adjustment of the balance before using it for the first time or after relocating it to ensure accurate weighing. For details, refer to the separate instruction manual listed in "1.1. Detailed Manual".

CAUTION

- ☐ Do not install the balance in areas where flammable or corrosive gases are present.


3.1. How to Adjust the Level of the Balance

Level the balance by adjusting the leveling feet so that the bubble of the bubble spirit level is centered in the black circle.


When the bubble is off to the left:

Turn the leveling foot on the front right clockwise.


When the bubble is off to the right:

Turn the leveling foot on the front left clockwise.


When the bubble is off to the back:

Turn both leveling feet on the front clockwise simultaneously.


When the bubble is off to the front:

Turn both leveling feet on the front counterclockwise simultaneously.


4. Display and Key Panel (Basic Operation)

Lit display:


Blinking display:

No.	Description
1	Processing indicator
2	Stabilization indicator
3	USB connection mark
4	When lit: Standby indicator When blinking: Automatic sensitivity adjustment notice
5	Displays weighing value and items
6	Unit display
7	Net mark
8	Gross mark
9	Preset tare mark
10	IR sensor indicator
11	When lit: Number of data instances (Data memory function) When blinking: Displays the function table setting value.
12	Response indicator (Lights up for about 30 seconds after the start of weighing)
13	Shock indicator
14	Gross zero mark
15	When lit: Interval output standby mode When blinking: Interval output mode


Keys

No.	Name
1	ON:OFF key
2	CAL key
3	MODE key
4	SAMPLE key
5	PRINT key
6	TARE key
7	RE-ZERO key
8	Left IR sensor
9	Right IR sensor

Key operations

Key operations affect how the balance functions. Normal key operation during measurement is “Press and release the key immediately” or “Press and hold the key (for 2 seconds)”. Please do not “Press and hold the key (for 2 seconds)” unless required.

Press the key
(Press and release the
key immediately.)

Press and hold the key
for 2 seconds.

Main weighing operations

Key	When pressed and released	When pressed and held (for 2 seconds)
	Turns the display on and off. When the display is turned off, only the standby indicator is displayed. When the display is turned on, weighing is possible. The ON:OFF key is active at any time, and pressing this key during operation always turns off the display.	Switches the IR sensor on and off.
	Performs sensitivity adjustment of the balance using the internal weight.	Displays the menu related to sensitivity adjustment.
	Switches the weighing units stored in the function table.	Enters “Changing Weighing Speed” mode.
	In weighing mode, press the key to turn the digit for readability on and off. In counting or percent mode, press the key to enter the sample storing mode.	☐ Enters the function table mode. For details, refer to the separate instruction manual listed in “1.1. Detailed Manual”. ☐ Runs the repeatability check function when pressed and held for another 2 seconds after the function table menu is displayed. For details, refer to the separate instruction manual listed in “1.1. Detailed Manual”.
	Stores or outputs the weighing value when stable according to the function table settings. (Data is output at factory setting.)	Enters mode to change the unit weight registration number in counting mode. The following settings are possible with the function table: ☐ Output “Title block” and “End block” for GLP/GMP report. ☐ Display the data memory menu.
	Performs tare operation.	
	Sets the displayed value to zero.	
IR SENSOR 	This is an IR sensor (touchless sensor). It reacts when you bring your hand close to it. Opening and closing of the breeze break door are assigned. For details, refer to “5. IR Sensors and Auto Doors”.	

5. IR Sensors and Auto Doors

5.1. IR Sensors

The BH series balances are equipped with IR sensors that enable operation without directly touching the display. By default, the IR sensors on the left and right sides of the display are assigned to open and close the breeze break door(s). By pressing and holding the **ON:OFF** key (for about 2 seconds), you can switch the IR sensors on and off with a single touch.

5.2. Auto Doors

The BH series balances feature a breeze break with auto doors that can be opened and closed without touching them. The IR sensors on the left and right of the display are assigned to open and close the breeze break. The door(s) connected to the joint(s) and handle(s) open and close automatically. By default, the breeze break opens to the last position it was opened to. You can also set it to fully open or half-open in the function table of the balance. Additionally, if you change the connection(s) of the joint(s), it is advisable to perform an auto door test with the function table. For details, refer to the separate instruction manual listed in “1.1. Detailed Manual”.


Opening the breeze break door(s)

Step	Description
1	When the breeze break is closed, place your hand over the right (or left) IR sensor.
2	The detection buzzer sounds, and the breeze break door(s) with the joint(s) and handle(s) connected open.

Closing the breeze break door(s)

Step	Description
1	When the breeze break is open, place your hand over the right (or left) IR sensor.
2	The detection buzzer sounds, and the breeze break door(s) with the joint(s) and handle(s) connected close.

Operation example

6. Sensitivity Adjustment

6.1. Automatic Sensitivity Adjustment

The sensitivity of the balance is automatically adjusted using the internal weight based on ambient temperature changes, set times, or interval times. For details, refer to the separate instruction manual listed in “1.1. Detailed Manual”.

6.2. Sensitivity adjustment using the internal weight

Sensitivity adjustment using the internal weight can be performed with one key press.












Step	Description
1	Be sure to warm up the balance with nothing on the weighing pan for at least an hour.
2	Press the [CAL] key. The balance displays [RL in] .
3	The internal weight is used to automatically adjust the sensitivity. Do not apply vibration and the like to the balance.
4	After sensitivity adjustment, the balance will output the “sensitivity adjustment report” and store the data as set in the function table. For details, refer to the separate instruction manual listed in “1.1. Detailed Manual”.
5	The balance returns automatically to weighing mode.













6.3. Sensitivity adjustment using an external weight

An external calibration weight is used to adjust the sensitivity.

For details, refer to the separate instruction manual listed in “1.1. Detailed Manual”.


7. Error Displays (Error Codes)

Display	Error code	Description and possible countermeasure
		Overload error The weighing value exceeds the balance's weighing capacity. Remove the object from the pan. If the error persists, repair is necessary.
		Weighing pan error The weighing value is too light. The weighing pan is not installed correctly. Set the weighing pan correctly. Perform a sensitivity adjustment. If the error persists, repair is necessary.
		Power supply voltage fault The voltage supplied from the AC adapter is abnormal. Check that the AC adapter is the one supplied with the balance.
	EC, E11	Stability error The weighing value is unstable, and "zero display", "sensitivity adjustment", etc. cannot be executed. Check around the pan. For more details, refer to the separate instruction manual listed in "1.1. Detailed Manual". Improve the environment of the installation location (vibration, drafts, static electricity, etc.). To return to weighing mode, press the CAL key.
		Entry value error The value entered exceeds the setting range. Enter a value within the setting range.
	EC, E16	Internal weight error Raising and lowering the internal weight does not yield a change in the mass value greater than that specified. Confirm that there is nothing on the pan and perform the operation from the beginning. If the error is not resolved after retrying, repair is necessary.
	EC, E17	Internal weight error The internal weight application mechanism does not function properly. Perform the operation from the beginning. If the error is not resolved after retrying, repair is necessary.
	EC, E20	Calibration weight error (Positive value) The calibration weight is too heavy. Check around the pan. Check the calibration mass value. To return to weighing mode, press the CAL key.
	EC, E21	Calibration weight error (Negative value) The calibration weight is too light. Check around the pan. Check the calibration mass value. To return to weighing mode, press the CAL key.
		Breeze break unit error The auto door is not functioning properly. Check if there are any obstructions around the door that might be hindering its operation. After confirming, either perform a door check with the function table or reconnect the power.
		Sample mass error The sample is too light to be stored as a sample mass for the counting mode or percent mode. The sample cannot be used.

Display	Error code	Description and possible countermeasure
  		Unit weight error The sample mass for the counting mode is too light. Storing and using it for counting may cause a counting error. Add samples to reach the specified number and press the PRINT key. Although pressing the PRINT key without adding samples will put the balance in counting mode, make sure to add samples for accurate counting.
 		Repeatability error The standard deviation (SD) of repeatability exceeded 50 d. ^{*1} Review the installation environment of the balance.  "SD Error" appears in repeatability display.  displayed in minimum weighing value (reference value) display. For details, refer to the separate instruction manual listed in "1.1. Detailed Manual". ^{*1} "d" represents scale division. Example: If readability is 0.1 mg, 1 d is 0.1 mg. (50 d is 5.0 mg)
 Blinking (alternately)		Full memory The number of stored weighing values has reached the upper limit. In order to store a new weighing value, it is necessary to delete data. For details, refer to the separate instruction manual listed in "1.1. Detailed Manual".
 Blinking (alternately)		Full memory The stored sensitivity adjustment / calibration test history has reached 50 results. In order to store a new result, the oldest history will be deleted. For details, refer to the separate instruction manual listed in "1.1. Detailed Manual".
		Clock battery error The clock backup battery has been depleted. Press any key and set the time and date. Even if the clock backup battery is depleted, the clock and calendar function works normally as long as the balance is powered with the AC adapter. If this error appears frequently, repair is necessary.
		Malfunction of the internal memory element of the balance If this error continues to be displayed, repair is necessary.
		Weighing sensor error If this error continues to be displayed, repair is necessary.
		Weighing sensor error Set the weighing pan correctly. If this error continues to be displayed, repair is necessary.

Display	Error code	Description and possible countermeasure
Error 8 Error 9		Abnormality in the internal memory data of the balance If this error continues to be displayed, repair is necessary.
	EC, E00	Communications error A protocol error occurred in communication. Check the format, baud rate, etc.
	EC, E01	Undefined command error An undefined command was found. Check the transmitted command.
	EC, E02	Not ready The received command cannot be executed. (e.g., Q command was received when not in weighing mode.) (e.g., Q command was received while re-zeroing.) Adjust the delay time for to transmitting a command.
	EC, E03	Timeout error With the timeout parameter set to t-UP 1 , there was a wait time of approximately 1 second or more while receiving command characters. Check the communication.
	EC, E04	Character length error The number of characters in the received command has exceeded the limit. Check the command to transmit.
	EC, E06	Format error The description of the received command is incorrect. (e.g., The number of digits of numerical values is incorrect.) (e.g., There are alphabet characters among the numerical values.) Check the transmitted command.
	EC, E07	Parameter setting error The value of the received command has exceeded the allowed value. Check the setting range of the numerical value of the command.
Other error displays		If there are any other error displays, or if any of the above errors continues to appear, repair is necessary.

7.1. Other display

The automatic sensitivity adjustment notice (the  indicator blinking) indicates that the automatic sensitivity adjustment will start.

If the balance is not used for a certain period of time with this indicator blinking, the balance automatically performs sensitivity adjustment using the internal weight. (The blinking period depends on the operating environment.)

Tips

- ☐ Although it is possible to continue using the balance even while this indicator is blinking, use after sensitivity adjustment is completed is advisable in order to maintain the weighing accuracy.

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




步骤	描述	零件图
3	从主机（1）正面将防风顶窗格（8）插入防风框架的前凹槽。	 The diagram illustrates the installation of the top windproof window frame (8) into the front slot of the main unit (1). The main unit (1) is shown from the front, and the top frame (8) is being inserted into the top slot.
4	将防风前窗格（7）插入主机（1）的前部下凹槽。然后，将窗格顶部推入锁栓，直到其卡入到位。	 The diagram illustrates the installation of the front windproof window frame (7) into the front lower slot of the main unit (1). The main unit (1) is shown from the front, and the front frame (7) is being inserted into the lower slot.
5	对于防风侧窗格（9、10），向上推主机（1）背面防风框架上的锁栓，直到其发出咔嚓声。确保把手朝外，将窗格从背面插入防风框架的凹槽中。 插入防风侧窗格后，将锁栓放回原位。	 The diagram illustrates the installation of the side windproof window frames (9 and 10) into the side slots of the main unit (1). The main unit (1) is shown from the side, and the side frames (9 and 10) are being inserted into the side slots. A detailed inset shows the locking mechanism: a handle is pushed up to lock the frame in place, and then returned to its original position after the frame is inserted.

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A&D AUSTRALASIA PTY LTD
A&D KOREA Limited
A&D Instruments India Private Limited