BH-T Series

(This document is included in the packing contents for the balance.)

Quick Start Guide
BH-T Series
Analytical Balances

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A&D Company, Ltd.

Warning Definition

The warnings described in this manual have the following meanings:

A potentially hazardous situation which, if not avoided, may result in minor or moderate injury or damage to the instrument.

CAUTION

Cautions to use the device correctly.

Note

Information or cautions to use the device correctly.

About This Manual

- (1) No part of this manual may be reprinted, copied, modified, or translated to another language without the prior written consent of A&D Company, Limited (A&D).
- (2) The contents of this manual are subject to change without notice.
- (3) Please contact A&D if you notice any uncertainty, errors, omissions, etc. in this manual.
- (4) A&D bears no liability for any loss or lost profits due to the operation of this product, and for direct, indirect, special, or consequential damages resulting from any defect in this product or this manual, even if advised of the possibility of such damage. Furthermore, A&D assumes no liability for claims of rights from third parties. Concurrently, A&D assumes no liability whatsoever for software or data losses.

© 2025 A&D Company, Limited

- ☐ Microsoft®, Windows®, Word®, and Excel® are trademarks of the Microsoft group of companies.
- ☐ The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by A&D is under license.
- ☐ iOS is the name of the operating system of Apple Inc. iOS is a trademark or registered trademark of Cisco in the U.S. and other countries and is used under license.
- ☐ Apple, the Apple logo and iPhone are trademarks of Apple Inc.
- ☐ App Store is a service mark of Apple Inc.
- ☐ Android™, Google Play and the Google Play logo are trademarks of Google LLC.
- ☐ Other product names and company names mentioned in this manual are trademarks or registered trademarks of their respective companies in Japan or other countries and regions.

1. Quick Start Guide

Thank you for purchasing an A&D BH-T series analytical balance. This quick start guide provides instructions on how to install the BH-T series balance and covers its basic functions and main operations. For more detailed information, please refer to the separate instruction manual listed in "1.1. Detailed Manuals".

1.1. Detailed Manuals

Detailed manuals are available for download from the A&D website (<https://www.aandd.jp>):

BH-T Series Instruction Manual

This manual provides detailed information on the functions and operations of the BH-T series, enabling you to fully utilize its features.

Scan the QR code to access the manual.

Ionizer (Static Eliminator) Instruction Manual

This manual describes the functions of the AD-1683A ionizer and provides instructions for its use.

2. Assembly and Installation

CAUTION


- ☐ The BH-T series analytical balance is a precision instrument and should be carefully unpacked. It is advisable to store the packing materials for future transport of the balance.
- ☐ The contents of the package vary depending on the product. Refer to "2.1. Unpacking" to ensure all items are included.
- ☐ Do not connect the AC adapter to the balance until it is fully assembled and installed.
- ☐ Once the balance is ready, connect the AC adapter. An automatic operation check will be performed.

CAUTION

- ☐ Confirm that the AC adapter type is correct for your local voltage and receptacle type.
- ☐ Use only the dedicated AC adapter specified for the balance.
- ☐ Do not connect the included AC adapter to other devices.
- ☐ Using the wrong AC adapter may cause the balance and other equipment to malfunction.

2.1. Unpacking

Packing contents and parts names for BH-T series:

- | | |
|--|---|
| 1 Main unit | 17 AD-1689 tweezers (for 0.01 mg models ^{*1} only) |
| 2 Weighing pan | 18 Spatula (for 0.01 mg models ^{*1} only) |
| 3 Low breeze break ring | 19 Weigh boat (10 pcs) |
| 4 High breeze break ring (for 0.01 mg models ^{*1} only) | 20 AD-1683A ionizer |
| 5 Breeze break bottom plate | 21 Pan support boss |
| 6 Display protection cover (PET plastic) ^{*2} | 22 Touchscreen |
| 7 Breeze break front pane | 23 Bubble spirit level |
| 8 Breeze break top pane | 24 Key |
| 9 Breeze break left side pane | 25 Leveling foot |
| 10 Breeze break right side pane | 26 Serial number |
| 11 Breeze break rear pane | 27 USB Interface (Type C) ^{*3} |
| 12 Quick Start Guide | 28 RS-232C interface |
| 13 AC adapters ^{*5} × 2 pcs | 29 AC adapter input jack |
| 14 AC adapter ID labels × 2 pcs | 30 Grounding terminal |
| 15 USB cable (approx. 2 m) | 31 External input terminal |
| 16 Cleaning brushes (large, small) | 32 Wired LAN port |
| | 33 USB interface (Type A) ^{*4} |

^{*1} BH-225TE/BH-225DTE

^{*2} Attached to the main unit.

^{*3} For communication only.



^{*4} For USB drive only.




^{*5} Accessories vary depending on the destination region.



2.2. Assembly

CAUTION

- ❑ Perform the following steps with the AC adapter removed from the main unit.

Step	Description	Parts diagram
1	Attach the breeze break bottom plate (5), breeze break ring (3 or 4), and weighing pan (2) to the main unit (1).	
2	Insert the breeze break rear pane (11) into the lower back groove of the main unit (1). Then, push the top of the pane into the latches until it clicks into place.	


Step	Description	Parts diagram
3	Insert the breeze break top pane (8) into the front groove of the breeze break frame from the front of the main unit (1).	 The diagram shows the main unit (1) with the breeze break frame. A top pane (8) is shown being inserted into the front groove of the frame. Arrows indicate the direction of insertion.
4	Insert the breeze break front pane (7) into the lower front groove of the main unit (1), then push the top of the pane into the latches until it clicks into place.	 The diagram shows the main unit (1) with the breeze break frame. A front pane (7) is shown being inserted into the lower front groove of the frame. Arrows indicate the direction of insertion and the latching mechanism.
5	For the breeze break side panes (9, 10), push the latches on the breeze break frame at the back of the main unit (1) upwards until they click. Insert the panes into the grooves of the breeze break frame from the back, ensuring the handles face outward. After inserting the breeze break side panes, return the latches to their original position.	 The diagram shows the main unit (1) with the breeze break frame. A side pane (9) is shown being inserted into the groove of the frame. A side pane (10) is shown being inserted into the groove of the frame. Arrows indicate the direction of insertion and the latching mechanism. A callout shows the latches being pushed upwards until they click.

Step	Description	Parts diagram
6	Attach the AC adapter ID labels (14) to the AC adapters (13). CAUTION ❑ Ensure the AC adapter ID labels are attached to avoid using the wrong AC adapters.	 The diagram illustrates the process of attaching identification labels (14) to AC adapters (13). It shows two examples: one for a 'BH' adapter and one for an 'AD1683A' adapter. Each example includes a label with the text '識別シール' (Identification Sticker), 'ACアダプタ、コード等に貼って下さい' (Attach to AC adapter, cord, etc.), and 'Identification Sticker. Attach Label to the AC Adapter and its Cable.' The labels are shown being attached to the adapters and their cables.
7	Insert an AC adapter (13) with the attached ID labels (14) into the AC adapter input jack (29) on the back of the main unit (1). Then, plug the other end into an outlet. CAUTION ❑ Be sure to warm up the balance for at least an hour before use.	 The diagram shows the AC adapter (13) being inserted into the AC adapter input jack (29) on the back of the main unit (1). The adapter has the ID label (14) attached to it.

3. Installation Considerations, Preparation and Precautions

Ensure the following installation conditions to achieve optimal performance of the balance:


- ☐ The ideal operating temperature is approximately $20^{\circ}\text{C} \pm 2^{\circ}\text{C}$ with a relative humidity of 45% to 60%.
- ☐ Install the balance in a dust-free environment.
- ☐ The weighing table should be solid. An anti-vibration table or stone table is ideal.
- ☐ Place the balance on a horizontal surface, ensuring it is not tilted.
- ☐ Install the balance in a stable location, avoiding areas with vibration and shock. Corners of rooms on the ground floor are preferable as they are less prone to vibration.
- ☐ Avoid installing the balance near equipment such as heaters or air conditioners.
- ☐ Ensure the balance is not exposed to direct sunlight.
- ☐ Keep the balance away from equipment that produces magnetic fields.
- ☐ Level the balance using the leveling feet and bubble spirit level. Refer to “3.1. How to Adjust the Level of the Balance”.
- ☐ Warm up the balance for at least an hour before use, with the AC adapter connected to the power supply.
- ☐ Perform a sensitivity adjustment of the balance before using it for the first time or after relocating it to ensure accurate weighing. For details, refer to the separate instruction manual listed in “1.1. Detailed Manuals”.





CAUTION

- ☐ Do not install the balance in areas where flammable or corrosive gases are present.

3.1. How to Adjust the Level of the Balance

Level the balance by adjusting the leveling feet so that the bubble of the bubble spirit level is centered in the black circle.


When the bubble is off to the left: Turn the leveling foot on the front right clockwise. 	When the bubble is off to the right: Turn the leveling foot on the front left clockwise. 
When the bubble is off to the back: Turn both leveling feet on the front clockwise simultaneously. 	When the bubble is off to the front: Turn both leveling feet on the front counterclockwise simultaneously. 

4. On-Screen Buttons and Key Panel (Basic Operation)


Standby screen

Weighing screen

Key/button operation

No.	Name	Description
①	 button	Runs a quick performance test. This test automatically checks the performance of the balance by loading and unloading the internal weight. For details, refer to the separate instruction manual listed in “1.1. Detailed Manuals”.
②	APP button	Displays the application setting screen. For details, refer to the separate instruction manual listed in “1.1. Detailed Manuals”.
③	Left IR sensor function change button	Changes the operation set for the IR sensors. By default, this button opens/closes the breeze break door(s).
④	Left IR sensor	This is one of the IR sensors (touchless sensor). When you bring your hand close to it, it reacts and the breeze break opens/closes. For details, refer to the separate instruction manual listed in “1.1. Detailed Manuals”
⑤	Stabilization indicator	Appears when the weighing value of the balance is stable.
⑥	 key	Turns the screen display on and off. The  key is active during any operation. When the displayed screen is turned off, the standby screen appears. When the screen display is turned on, the weighing screen appears.
⑦	 key	Displays the weighing screen. The  key is active during any operation.
⑧	 button	Outputs data to the device connected to the balance. For details, refer to the separate instruction manual listed in “1.1. Detailed Manuals”.
⑨	 button	Displays the sensitivity adjustment/calibration test screen. For details, refer to the separate instruction manual listed in “1.1. Detailed Manuals”.
⑩	Communication device setting button	Displays the communication device setting screen. Configures the settings of the connected communication device.
⑪	Right IR sensor function change button	Changes the operation set for the IR sensors. By default, this button opens/closes the breeze break door(s).
⑫	Right IR sensor	This is one of the IR sensors (touchless sensor). When you bring your hand close to it, it reacts and the breeze break opens/closes. For details, refer to “5.2. Auto Doors”.
⑬	Capacity, readability	Displays the capacity and readability of the balance.
⑭	Date/time display	Displays the current date and time.
⑯	 key	Displays the log-in screen. The  key is active at any time, and pressing this key during operation always displays the log-in screen. For details, refer to the separate instruction manual listed in “1.1. Detailed Manuals”.
⑰	 key	Displays the MENU screen. The  key is active during any operation.
⑱	 button	Sets the displayed value to zero.

Warning display

No.	Warning display	Name	Description	Display priority
⑮		Shock indicator	Displayed by the shock detection function.	High
		Static elimination recommended	Displayed when the humidity inside the balance is 45% or less. (Lights up for about 30 seconds after the start of weighing)	Low

5. IR Sensors and Auto Doors

5.1. IR Sensors

The BH-T series balances are equipped with IR sensors that enable operation without directly touching the display. By default, the IR sensors on the left and right sides of the display are assigned to open and close the breeze break door(s). The left and right **DOOR** buttons can be used to turn the IR sensors on/off, or to change other functions.

5.2. Auto Doors

The BH-T series balances feature a breeze break with auto doors that can be opened and closed without touching them. The IR sensors on the left and right of the display are assigned to open and close the breeze break. The door(s) connected to the joint(s) and handle(s) open and close automatically. By default, the breeze break opens to the last position it was opened to. You can also set it to fully open or half-open in the function table of the balance. Additionally, if you change the connection(s) of the joint(s), it is advisable to perform an auto door test with the function table. For details, refer to the separate instruction manual listed in “1.1. Detailed Manuals”.


Opening the breeze break door(s)

Step	Description
1	When the breeze break is closed, place your hand over the right (or left) IR sensor.
2	The detection buzzer sounds, and the breeze break door(s) with the joint(s) and handle(s) connected open.

Closing the breeze break door(s)


Step	Description
1	When the breeze break is open, place your hand over the right (or left) IR sensor.
2	The detection buzzer sounds, and the breeze break door(s) with the joint(s) and handle(s) connected close.

Operation example

6. Check Function

The check function performs both daily and regular inspections.

Additionally, the balance's performance can be easily verified by checking the repeatability with the internal weight.

The check function can be selected from the menu screen. For details, refer to the separate instruction manual listed in "1.1. Detailed Manuals".


7. Sensitivity Adjustment/Calibration Test

7.1. Automatic Sensitivity Adjustment

The sensitivity of the balance is automatically adjusted using the internal weight based on ambient temperature changes, set times, or interval times. For details, refer to the separate instruction manual listed in “1.1. Detailed Manuals”.

7.2. Sensitivity Adjustment Using the Internal Weight

The sensitivity adjustment/calibration test screen is displayed with the  CAL button on the weighing screen.

No.	Name	Description
①	External sensitivity adjustment	An external calibration weight is used to adjust the balance. Instructions are displayed on the screen for how to adjust the sensitivity.
②	Internal sensitivity adjustment	The internal weight is used to automatically adjust the balance.
③	Internal calibration test	The internal weight is used to check the accuracy of weighing, and the result is output. Note: The balance is not adjusted.
④	External calibration test	An external calibration weight is used to check the accuracy of weighing, and the result is output. Note: The balance is not adjusted.

Do not allow vibration or drafts to affect the balance during sensitivity adjustment/calibration test. For details, refer to the separate instruction manual listed in “1.1. Detailed Manuals”.

8. Typical Errors (Error Codes)

Display	Error code	Description and possible countermeasure
E		Overload error The weighing value exceeds the balance's weighing capacity. Remove the object from the pan.
-E		Weighing pan error The weighing value is too light. The weighing pan is not set correctly. Set the weighing pan correctly. Perform a sensitivity adjustment.
LoWVOLT		Power supply voltage fault The voltage supplied from the AC adapter is abnormal. Check that the AC adapter is the one supplied with the balance.
Error 1	EC, E11	Stability error The weighing value is unstable, and "zero display", "sensitivity adjustment", "calibration test", etc. cannot be executed. Check around the pan. For details, refer to the separate instruction manual listed in "1.1. Detailed Manuals". Improve the environment of the installation location (vibration, drafts, static electricity, etc.).
Error 2		Entry value error The value entered exceeds the setting range. Enter a value within the setting range.
Error 6	EC, E16	Internal weight error Raising and lowering the internal weight does not yield a change in the mass value greater than that specified. Confirm that there is nothing on the pan and perform the operation from the beginning. If this error continues to be displayed, repair is necessary.
Error 7	EC, E17	Internal weight error The internal weight application mechanism does not function properly. Perform the operation from the beginning. If this error continues to be displayed, repair is necessary.
CAL E	EC, E20	Calibration weight error (Positive value) The calibration weight is too heavy. Check around the pan. Check the calibration mass value.
-CAL E	EC, E21	Calibration weight error (Negative value) The calibration weight is too light. Check around the pan. Check the calibration mass value.
door Err		Breeze break unit error The auto doors are not functioning properly. Check if there are any obstructions around the doors that might be hindering their operation. After confirming, either perform a door check with the function table or reconnect the power.
SD Error		Repeatability error With the minimum weighing warning function, the standard deviation (SD) of repeatability exceeded 50 d. Review the installation environment of the balance. For details, refer to the separate instruction manual listed in "1.1. Detailed Manuals". ☐ SD Error "SD Error" appears in repeatability display.
rtc PF		Clock battery error The clock backup battery has been depleted. Set the time and date. Even if the clock backup battery is depleted, the clock and calendar function will work normally as long as the balance is powered with the AC adapter. If this error appears frequently, repair is necessary.
Error 3		Malfunction of the internal memory element of the balance If this error continues to be displayed, repair is necessary.

Display	Error code	Description and possible countermeasure
Error 8		Abnormality in the internal memory data of the balance If this error continues to be displayed, repair is necessary.
Error 9		Abnormality in the internal memory data of the balance If this error continues to be displayed, repair is necessary.
	EC, E00	Communications error A protocol error occurred in communication. Check the format, baud rate, etc.
	EC, E01	Undefined command error An undefined command was found. Check the transmitted command.
	EC, E02	Not ready The received command cannot be executed. (e.g., Q command was received when not in weighing mode.) (e.g., Q command was received while re-zeroing.) Adjust the delay time for transmitting a command.
	EC, E03	Timeout error With the command timeout setting, there was a wait time of approximately 1 second or more while receiving command characters. Check the communication.
	EC, E04	Character length error The number of characters in the received command has exceeded the limit. Check the command to transmit.
	EC, E06	Format error The description of the received command is incorrect. (e.g., The number of digits of numerical values is incorrect.) (e.g., There are alphabet characters among the numerical values.) Check the transmitted command.
	EC, E07	Parameter setting error The value of the received command has exceeded the allowed value. Check the setting range of the numerical value of the command.
Other error displays		If any other errors appear, or if any of the above errors continue to appear, repair is necessary.

9. Ionizer AD-1683A

The AD-1683A ionizer eliminates static electricity by emitting positive or negative ions, generated from four discharge electrodes through DC corona discharge. Typically, insulators such as powders, filters, and weighing paper become charged when the humidity is 45% RH or lower, potentially causing weighing errors of a few milligrams. The ionizer neutralizes static electricity on the object weighed, preventing errors and ensuring accurate measurements. For more information, refer to the detailed instruction manual listed in “1.1. Detailed Manuals”.

9.1. Using the Ionizer


Step	Description
1	Connect the AC adapter, plug it in, and turn on the power. The power lamp on the ionizer will light up.
2	Place the target object within the effective range of static elimination.
3	Hold your hand over the IR sensor (touchless infrared proximity sensor) on the front of the ionizer to start static elimination. The ION lamp will light up, indicating that static elimination is in progress.
4	Refer to the figure below for the effective range of static elimination, and perform static elimination. By default, static elimination will stop, and the ION lamp will turn off 3 seconds after the start of static elimination.

Ionizer AD-1683A

IR sensor (touchless infrared proximity sensor)


Power lamp (lights up when the power is turned on)

ION lamp (lights up during static elimination)

Effective range of static elimination

Top view

Side view


Figure. Static elimination range

10. Specifications

	BH-225TE	BH-225DTE
Weighing capacity	220 g	220 g 51 g
Maximum display	220.00084 g	220.0008 g 51.00009 g
Readability	0.01 mg (0.00001 g)	0.1 mg (0.0001 g) 0.01 mg (0.00001 g)
Repeatability (Standard deviation)	0.015 mg (50 g) 0.03 mg (200 g)	0.1 mg (200 g) 0.025 mg (50 g)
Linearity	±0.10 mg	±0.2 mg
Communication	RS-232C (printer, PLC, etc.), USB type A (USB drive, for data storage), USB type C (PC), LAN (TCP/IP), Stereo jack (external switch), <i>Bluetooth</i> ® (PC, etc.)	

	BH-224TE	BH-324TE
Weighing capacity	220 g	320 g
Maximum display	220.0084 g	320.0084 g
Readability	0.1 mg	
Repeatability (Standard deviation)	0.09 mg	0.1 mg
Linearity	±0.2 mg	
Communication	RS-232C (printer, PLC, etc.), USB type A (USB drive, for data storage), USB type C (PC), LAN (TCP/IP), Stereo jack (external switch), <i>Bluetooth</i> ® (PC, etc.)	

11. Compliance

11.1. Compliance with FCC Rules

This device contains transmitter module FCC ID: 2A6NFEC2832

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

(1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. (FCC = Federal Communications Commission in the U.S.A.)

FCC CAUTION: Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. The antenna used for this transmitter must not be co-located or operating in conjunction with any other antenna or transmitter. This product is certified as type of the portable device with FCC Rules. To maintain compliance with RF Exposure requirement, please use within specification of this product.

11.2. IC

IC RADIATION EXPOSURE STATEMENT FOR CANADA

This device contains transmitter module IC: 28568-EC2832

This device complies with Innovation, Science and Economic Development Canada license-exempt RSS standards. Operation is subject to the following two conditions: (1) This device may not cause interference. (2) This device must accept any interference, including interference that may cause undesired operation of the device.

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11.3. Bluetooth®

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by A&D is under license.

Other trademarks and trade names are those of their respective owners.

THIS PAGE INTENTIONALLY LEFT BLANK.

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




步骤	描述	零件图
3	从主机（1）正面将防风顶窗格（8）插入防风框架的前凹槽。	 The diagram illustrates the installation of the top windproof window frame (8) into the front slot of the main unit (1). The main unit is shown from the front, and the top frame is being inserted into the top slot.
4	将防风前窗格（7）插入主机（1）的前部下凹槽，然后将窗格顶部推入锁栓，直到其卡入到位。	 The diagram illustrates the installation of the front windproof window frame (7) into the front lower slot of the main unit (1). The main unit is shown from the front, and the front frame is being inserted into the lower slot.
5	对于防风侧窗格（9、10），向上推主机（1）背面防风框架上的锁栓，直到其发出咔嚓声。确保把手朝外，将窗格从背面插入防风框架的凹槽中。插入防风侧窗格后，将锁栓放回原位。	 The diagram illustrates the installation of the side windproof window frames (9, 10) into the side slots of the main unit (1). The main unit is shown from the side, and the side frames are being inserted into the side slots. A detailed inset shows the locking mechanism: a lock pin is pushed up to lock the frame in place, and then pushed back down to release it.

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




步驟	說明	零件圖
3	將頂部防風板 (8) 從主機 (1) 的前方插入防風架的前側凹槽。	 The diagram illustrates the installation of the top wind deflector (8). It shows the main unit (1) with a circular lens. A rectangular plate (8) is shown being inserted into a slot on the front of the unit. Arrows indicate the direction of insertion.
4	將前防風板 (7) 插入主機 (1) 前側下方的凹槽，然後將防風板頂端推入門鎖，直到固定到位。	 The diagram shows the front wind deflector (7) being inserted into the front slot of the main unit (1). The deflector is shown at an angle, with arrows indicating it should be pushed forward and then the top edge should be pushed into the lock mechanism.
5	接著安裝左右防風板 (9、10)，然後將主機 (1) 後側防風架的門鎖向上推，直到固定到位。將防風板從後側插入防風架的凹槽中，手把務必朝外。插入左右防風板後，將門鎖推回原位。	 The diagram illustrates the installation of the side wind deflectors (9 and 10). It shows the main unit (1) with a side slot. A detailed inset shows a close-up of the lock mechanism being pushed up. The side deflectors (9 and 10) are shown being inserted into the side slot. Arrows indicate the direction of insertion and the movement of the lock mechanism.

步驟	說明	零件圖
6	將交流電轉接器 ID 標籤 (14) 貼到交流電轉接器 (13)。 警告 ❑ 務必貼上交流電轉接器 ID 標籤，避免用錯交流電轉接器。	 Diagram illustrating the identification labels (14) for the AC adapters (13). The labels are shown for two models: BH and AD1683A. The labels are attached to the AC adapters and their cables. The labels are shown both separately and being applied to the adapters.
7	將貼有 ID 標籤 (14) 的交流電轉接器 (13) 插入主機 (1) 後側的交流電轉接器輸入插孔 (29)。然後將另一端插入插座。 警告 ❑ 使用天平前，務必先預熱至少一小時。	 Diagram illustrating the connection of the AC adapter (13) to the main unit (1). The AC adapter (13) is shown being inserted into the AC adapter input port (29) of the main unit (1). The AC adapter (13) has the ID label (14) attached to it.

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A&D AUSTRALASIA PTY LTD
A&D KOREA Limited
A&D Instruments India Private Limited