# AD-8121

MULTI-FUNCTION PRINTER

# INSTRUCTION MANUAL

Instruction-AD-8121-v.1.b 92.05.12.OGA

MULTI-FUNCTION PRINTER



# Table of Contents

FCC Rules	Page iii
Section A	
Features	Page A•2
Unpacking your AD-8121	Page A•3
Items Supplied	Page A•3
Handling Precautions	Page A•4
External View	Page A•5
Installing the Battery	Page A•6
installing the Ink Ribbon Cassette and Printer Paper	Page A•7
Installing the Ink Ribbon Cassette	Page A•7
Installing the Printer Paper	Page A•9
Section B	
Control Functions	Page B•2
Keyboard	Page B•2
DIP Switches	Page B•4
Print Format	Page B•5
Data Format When Printing In Mode 1 or Mode 2	Page B•6
Section C	
Setup	Page C•2
Power Supply	Page C•3
Battery Operation	Page C•3
Section D	
Calendar/Clock Function	Page D•2
Date/Time Settings	Page D•3
Example of Date/Time Settings	Page D•5
Interval Printing and Chart Printing	Page D•6
Using theInterval Printing in the Normal Mode	•
Using the Interval Printing in the Chart Mode	-
Section E	
Statistical Calculation Mode	Page F•2

Selection of Standard Mode or	
Statistical Calculation Mode	Page E•2
Samples and Contents of	
Statistical Calculation Mode Printing	Page E•3
Section F	
Dump Print Mode (MODE 3)	Page F•2
Printing in MODE 3	Page F•2
Online and Offline	Page F•2
Other Keys	Page F•2
Specifications	Page F•3
General Specifications	Page F•3
Input Specifications	Page F•3
Statistical Calculation Specifications	Page F•4

# FCC Rules


Please note that this equipment generates, uses and can radiate radio frequency energy. This equipment has been tested and has been found to comply with the limits of a Class A computing device pursuant to Subject J of Part 15 of FCC rules. These rules are designed to provide reasonable protection against interference when equipment is operated in a commercial environment. If this unit is operated in a residential area it might cause some interference and under these circumstances the user would be required to take, at his own expense, whatever measures are necessary to eliminate the interference.

(FCC = Federal Communications Commission in the U.S.A.)

### **Section A**

This section describes features of the AD-8121, the items supplied, handling precautions, external view, and how to install the batteries, printer paper and ink ribbon cassette.

Features	A•2
Unpacking Your AD-8121	
Items Supplied	
Handling Precautions	
External View	
Installing the Battery	
Installing the Ink Ribbon Cassette and Printer Paper	
Installing the Ink Ribbon Cassette	
Installing the Printer Paper	



### Thank you for your AD purchase!

This is an Instruction Manual for the AD-8121 Multi-Function printer. The AD-8121 is a product of years of design, development, and in-field testing. It is designed to withstand harsh environmental conditions - and each printer is subjected to several levels of quality control before it leaves the factory. Every care has been taken during the manufacturing process of this printer to ensure that it will perform accurately and reliably for many years.

The AD-8121 Multi-Function Printer is mainly for use with A&D's electronic balances, counting balances, platform scales and weighing indicators. A highly reliable printer mechanism is perfect for industrial use. You can record and statistically analyze the measuring data obtained from the weighing instruments. It is easy to operate and measurements can be processed on the spot. You can operate the AD-8121 using either AC adaptor or alkaline batteries. The AD-8121 Multi-Function Printer provides fast dependable printing.

# Features


- Impact dot matrix printer for superior readability
- A full range of statistical functions: Weight data, Total weight data, Counting data Total counting data, Number of operations, Standard deviation, Chart, Year, Month, Date, Hour, Minute and Second
- Interval printing with built-in Timer/Clock function: 5s, 10s, 30s, 1m, 5m, 10m, 30m
- Charting selectable to show change in weight, Interval time, Start time and Stop time
- Battery operation permits use anywhere
- Clock, calendar and pre-set functions are stored in non-volatile memory
- Compact and light weight
- Accepts RS-232C and Current loop input from A&D's Electronic Balances, Scales and Indicators
- Quiet and speedy printing




# **Unpacking Your AD-8121**




### Items Supplied





### **Accessories (Sold seperately) & Option**



\* Option: OP-01 Adaptor Cable for current loop (50 cm /20 inches)



Adaptor cable for current loop




# Handling precautions


In order to obtain the best possible performance from the AD-8121 instruction, read this operation manual thoroughly before operation.

Be careful when handling this unit because it is a high-precision electronic device. Note the following:

- As foreign matter in the printer section may cause the printer to jam. Use this
  unit where it will not be subjected to dirt, dust, water, or oil, and where it will not
  be subjected to direct sunlight or high temperature and high humidity.
- In order to protect the high-precision mechanism and electrical circuits of the AD-8121 from being damaged, prevent the unit from being struck or having an excessive force applied to it during its operation. Do not attempt to disassemble it.
- Be sure to connect/disconnect the plug connectors of the cables only when the power is OFF.
- Only use the specified printer paper and AC adaptor.
- Only use alkaline cells.
- Remove the batteries when not using the printer for a long time to prevent potential damage from battery leakage.
- Wait at least two seconds to switch on again after turning off the power supply switch.




- If a malfunction occurs, refer to this manual to check the power source, switch operations, and cable connections.
  - When the exterior of this unit has stains, lightly wipe them off using a dry cloth or a cloth dampened in neutral detergent. Do not use an organic solvent for cleaning.






# Installing the Battery


• Step 1. Open the battery box.



• Step 2. Install the batteries.



• Step 3. Close the battery box.







# Installing the Ink Ribbon and Printer Paper

### Installing the Ink Ribbon Cassette


• Step 1. Slide the printer paper cover in the direction of the arrow ①, then remove it in the upward direction.




Step 2. Remove the printer cover.



• Step 3. Installing the ink ribbon cassette supplied with your printer.







Replace the printer cover and printer paper cover in the reverse order of their removal.

### Installing the Printer Paper


• Step 1. Slide the printer paper cover in the direction of the arrow ①, then remove it in the upward direction.



• Step 2. Remove the printer cover.




• Step 3. Cut the leading edge of the printer paper at a right angle to the feeding direction and remove wrinkles and folds from the paper if there are any.



Installing the Ink Ribbon and Printer Paper


● Step 4. Turn on the power. Press the FEED key while lightly feeding the printer paper edge into the paper inlet until it comes out of the outlet an appropriate length. Insert the supplied shaft into the printer paper core, then set the shaft onto the bearings of the printer.



# **Section B**

This section describes the control function, print procedure, print format and data format when printing in Mode 1 or Mode 2.

Control Functions	B•2
Keyboard	B•2
DIP Switches	
Print Format	B•5
Data Format When Printing In Mode 1 or Mode 2	







# **Control Functions**



### Keyboard





### **DIP Switches**

#### MODE



DIP Switch No.	ltem	OFF	ON
NO. 1	Mode selection	MODE 1	MODE 2
³NO. 2	Mode selection	MODE1/MODE2	MODE 3
NO. 3	Unstable data	No printing	Printing
NO. 4	Data input specifications	RS-232C	Current loop

□ DIP switch No. 1 and No. 2 (Mode selection)
The printer has three types of print modes (MODE 1, MODE 2 and MODE 3) as described below.

Table 1

	MODE 1	MODE 2	MODE 3
Printing	Prints data when data is received.	Prints data in conjunction with the DATA key or interval setting.	Prints data received as it is.
Data output mode for connected weighing instrument Key mode	Key mode Auto printing	Stream mode	Key mode Auto printing
Interval printing	No	Yes	No
Statistical calculation	Yes	Yes	No
Others	Data on such A & D pr weighing instruments can when using the standard of	only be printed correctly	See Section F•2 "Dump Print Mode," for details.

 Use the DIP switches on the right side of the printer to select the mode.

	Switch No.	4	3	2	1	\$153 ¢
	MODE 1	-	-	OFF	OFF	
	MODE 2		-	OFF	ON	
	MODE 3	-	-	ON	-	
L						

- ☐ DIP switch No. 3 (Printing unstable data)

  If "No printing" is selected (with DIP switch No. 3 set OFF), unstable data will be printed as asterisks (\*). (If weight data is "OVER", asterisks (\*) will be printed regardless of the switch No. 3 setting.)
- ☐ DIP switch No. 4 (Data input specifications)

  This switch is used to select the RS-232C interface or current loop.

  See Section F•3, "Input Specifications," for details.

### **門** Print Format

In MODE 1 and MODE 2, the AD-8121 only prints data (on such A&D products as balances and weighing instruments) set in the standard output format. When the standard format is used, the AD-8121 prints data as shown below. (To print output data as it is, set DIP switch No. 2 ON to select MODE 3. See Section F•2,"Dump Print Mode," for details.)

#### Printing Example in the MODE1 or MODE2

The printed data consists of alphanumeric characters and special symbols.

"/" is used for the calendar, as in Year/Month/Date.

":" is used for the time, as in Hour: Minute: Second. AM is used to indicate morning hours and PM (or P) is used to indicate afternoon and evening hours.

Communication error message (\*C ERROR\*) will be printed in case of an incorrect baud rate setting, improper cable connection, or no reception when the DATA key is pressed in MODE 2.

#### - Printing Example -

1:47:49 PM
Weight
NET weight
GROSS weight
TARE weight
Unit weight
Quantity
Accumulated quantity
Unstable data
Communication
Error

01:47:4	19 PM	
WT	123.456	g
NT	6432.15	kg
GS	9876.50	t
TR	200.55	g
UW	15.283	kg
QT	20000	PC
AQ	6543210	PC
	•	
*C ERF	ROR*	

#### Error messages

If an error occurs, one of the messages listed in Table-2 will be printed.


Table-2


Error message	Probable cause
*	<ul><li>① Unstable or "over" data was input.</li><li>② Data which is different from the specified format was input.</li></ul>
*C ERROR*	<ol> <li>Inproper cable connection.</li> <li>Baud rate do not match. (all should be set to 2400 bps.)</li> <li>No power for a balance, etc.</li> </ol>
*F ERROR*	① Printer malfunction due to noise.
*NO DATA*	① No data received from a balance, etc.



### Data Format When Printing In Mode 1 or Mode 2

To perform printing using Mode 1 or Mode 2, a defined data format must be used. When you connect the AD-8121 with instruments other than A&D products, you must match your instrument's output to A&D format. This format is shown below.





#### Header 1, Header 2

Decide which two characters will indicate the data type.



Select these from Table 1 on the next page.

Header 1, header 2, and numeric data are respectively separated by ", ". You can omit header 2.

#### Numeric data

Numeric data is a character string which contains a decimal point (•) and requires polarity sign (+ or -). The positive sign (+) will not be printed. Up to 7 digits can be used for the numeric part.

#### <u>Unit</u>

Up to three alpha characters or spaces can be used. They will be printed right - justified as in the example. The specified character string will be converted automatically as shown in Table 2 on the next page.

When there is a terminator after the numeric data, the defined character string will be printed automatically by header 1 and header 2, see Table 1.

#### **Terminator**

The printer starts printing when it receives a terminator. The terminator is fixed as CR•LF (0DH, 0AH).

Table 1

Abbreviation for the type of weight data	Header 1	Header 2	Character string printed automatically when there is no data indicating unit.
AQ	AQ	*1	цРС
GS	ST	GS	Ц <b>П</b> Т
NT	ST	NT	⊔Kg
QT	QT	*1	∐PC
TR	ST	TR	ULIT
ÜW	UW	*1	UU9
WT	ST	*1	ППâ

<sup>&</sup>quot;\*1" indicates that you can omit header 2.

Data except data format shown above is printed as " \* ".

Table 2


Received data	Printing
*2G	<u>ппа</u>
KG, Kg, kG	Цkg
LB, Lb, IB	Шlb
OZ, Oz	□oz

<sup>&</sup>quot; \*2 " indicates any character.

# Section C

This section describes the setup and power supply.

Setup	C•2
Power Supply	
Battery Operation	






### Setup




Verify that power of both the AD-8121 and the weighing instrument are turned off before connecting the units together.

1) Connection to weighing Instruments.

For data supply from an RS-232C interface, plug the 25 pin connector into the output connector of the Weighing Instrument and plug the smaller connector into




- Connecting the specified AC Adaptor
   For AC operation, connect the specified AC Adaptor to the adaptor jack on the printer.
- 3) Connecting the specified Foot switch (SW:128) For data acquisition by pressing the foot switch. Connect the foot switch to the EXT. P jack on the printer.



#### **IMPORTANT**

Please order the proper AC Adaptor you need according to the AC Adaptors list.



#### **AC ADAPTORS**

Model #	Voltage	Plug
AX-TB105	100	Α
AX-TB106	120	Α
AX-TB107	220	-
AX-TB107-BF	220	BF
AX-TB107-S	220	S
AX-TB121	220	С
AX-TB122	240	-
AX-TB122-BF	240	BF
AX-TB122-C	240	С
AX-TB169	240	S
		L



The power supplied for this unit is either from the AC Adaptor or the alkaline batteries. When the AC adaptor is used while the alkaline batteries are installed, the power will be supplied from the AC adaptor.



When the power to this unit is turned off, all the statistical calculation data stored in memory will be erased. Also, if the AC adaptor is connected/disconnected while the unit is operating (disregarding whether the unit is operated by the AC adaptor or the batteries), the data accumulated may be erased.

### **Battery Operation**

- Pay close attention to the battery polarity. (See Section. A•6, "Installing the Battery" for details.)
- Only use alkaline cells. Manganese cells offer inadequate power to operate the printer.
- A new battery will last about 3500 data operations under the following conditions.

Temperature conditions	25°C/77°F
Print conditions	Prints the time. Prints WT,+8888.888g every five seconds in the non-statistical calculation modes (MODE 1). Inputs data via the RS-232C interface.




When not using the printer for a long period, be sure to remove the batteries to prevent potential damage from battery leakage.

# Section D

This section describes the calendar/Clock function (date/time setting) and chart printing function.

Calendar/Clock Function	D•2
Date/Time Settings	
Example of Date/Time Settings	
Interval Printing and Chart Printing	
Using Interval Printing in the Normal Mode	
Using Interval Printing in the Chart Mode	






### Calendar/Clock Function

The AD-8121 has a built-in Calendar/Clock backed up by a lithium battery. The Calendar/Clock continues to operate with the power turned off.


- Three kinds of calendar formats are available: Year/Month/Day, Month/Day/ Year, Day/Month/Year
  - The time is set in 12 hours mode. (01:00:00 to 12:59:59)
- The printing or non-printing of the date and time can be selected. Otherwise, the date and time will be printed on each line independently. In the statistical calculation mode, the date is not printed.
- In MODE 3, the date and time can be printed by using control codes sent from the data source.
- Use the keys to set the date and time. The setting can be checked by printing it out.
- Leap years are not considered. Therefore, February 29th is always March 1st.

### Date/Time Settings



Press the SET key while holding down the
 RSLT key to enter the "Date/time" setting

In this mode, the keys function as follows.




SET ! (Set key) Enters the specified data and shifts to the next setting item. CL (Clear key) Selects the date or time for printing or no printing. CE (Cancel key) Selects the date format or AM or PM. STAT. Each time this key is pressed, the target digit of the (Stat key) setting value is incremented by one. RSLT. (Result key) Each time this key is pressed, the target digit is shifted to the right. DATA (Data key) Prints the data set for checking. FEED (Feed key) Each time this key is pressed, the target digit of the setting value is decremented by one.

#### Calendar/Clock Function


 In the date/time setting mode, the date is printed using the stored date format.

The date format will be rotated through the three types as follows by pressing the CE key.



Pressing the CL key while setting the date or time lets you select printing or no printing.

During"Printing" status, the green LED on the panel will light.



• When the SET key is pressed, the following information regarding printing and no printing is printed out.

"Printing" \*DATE PRINT ON\* \*TIME PRINT ON\*

"No printing" \*DATE PRINT OFF\* \*TIME PRINT OFF\*

• In the date/time setting mode, the following is printed when the SET key is pressed after the Date/time has been set.

\*The following is printed after \*DATE PRINT ON\* or \*DATE PRINT OFF\*.

\*DATE PRINT ON\*

PM 01:46:46

or

\*DATE PRINT OFF\*

AM 10:16:13

\*The following is printed after \*TIME PRINT ON\* or \*TIME PRINT OFF\*

\*TIME PRINT ON\*

91/11/05 01:46P

or

**\*TIME PRINT OFF\*** 

91/11/05 10:16A



The "second" is cut off when the setting time is changed. However, the "second" remains as it is when no setting time is changed ( RSLT. or FEED is not pressed).

For example,

SET AM 09:12:26

STAT.

(Left 2 minutes)

SET 91/07/12 09:14A

The date/time setting procedure is described on the next page.

# **Date/Time Settings Procedure**

Suppose that the date/time is set using the following example:

Date: Month/Date/Year, 11/03/91

Time: 10:05 AM Date: No printing Time: Printing

The data enclosed in parentheses "( )" is not printed.

SET + RSLT. Start Y/M/D91/11/01 LED ON (Previous setting) CE Date format changing (M/D/Y11/01/91) G:WT/R:QT LED on M/D/Y11/01/91 DATA Print checking (Green) ON LINE RSLT. Shifts ^ 1 digit to the right. (11/01/91)RSLT. Shifts ^ 1 digit to the right. (11/01/91)Shifts ^ 1 digit to the right. RSLT. (11/01/91) G:WT/R:QT LED off STAT. +1 (11/02/91)ON LINE **LEDOFF** STAT. +1 (11/03/91)CL No printing the date SET Enter the date \*DATE PRINT OFF\* PM01;23:45 LEDON G:WT/R:QT LED on CE AM/PM selection (Green) (AM 01:23:45) ON LINE STAT. +1 (11:23:45)RSLT. Shifts ^ 1 digit to the right. (11:23:45) FEED -1 (10:23:45)RSLT. Shifts ^ 1 digit to the right. (10:23:45) FEED -1 (10:13:45) FEED -1 (10:03:45)DATA Print checking AM 10:03:45 RSLT. Shifts ^ 1 digit to the right. (10:03:45) STAT. +1 (10:04:45) LEDOFF (Non-statistical +1 STAT. (10:05:45)calculation) SET Date/time setting G:WT/R:QT completed **\*TIME PRINT ON\*** - LED off ON LINE

11/03/9110:05A



### **Interval Printing and Chart Printing**


The AD-8121 can print input data at constant intervals using the built-in Calendar/Clock function. Characters (Normal mode) and charts (Chart mode) can also be printed. An interval of 5, 10 or 30 seconds, and 1, 5, 10 or 30 minutes can be set. Interval printing can be done in MODE 2. Set output data print timing of the connected weighing instruments to "Stream mode".

### ===

### Using the Interval Printing in the Normal Mode

• Press the CE key.

Now, "No chart printing" is selected. See Table-4.



#### Table-4

Target digit	Key	Print
No chart printing	CE	*CHART OFF*
10 ¹ - 10 °	CL	6543210
10 <sup>2</sup> - 10 <sup>1</sup>	SET	6543210
10 <sup>3</sup> - 10 <sup>2</sup>	RSLT.	6543210
10 <sup>4</sup> - 10 <sup>3</sup>	STAT.	6543210
10 <sup>5</sup> - 10 <sup>4</sup>	DATA	6543210
10 <sup>6</sup> - 10 <sup>5</sup>	FEED	6543210

(Table-4 also shows that the CE key selects "No chart printing", and that the CL, SET, RSLT., STAT., DATA and FEED keys set the target digit in the chart printing mode.)

#### Interval setting

This function is used for automatically logging data from the weighing instruments at regulat intervals. Hold down the SET key and press the FEED key to enter the interval setting mode.

\*INTERVAL TIME\* is then printed.

The length of the interval can be determined by

SET

STAT.


CE

RSLT.

pressing any of the keys (see the table below for the correspondence between times and keys) to be pressed subsequently.

Table-5

Interval	Key	
5 seconds	FEED	
10 seconds	DATA	
30 seconds	STAT.	
1 minute	RSLT.	
5 minutes	SET	
10 minutes	CL	
30 minutes	CE	



#### Print start command after interval setting

• When the interval is set, press the DATA key. Interval printing then begins.



If no data is received, the printer waits until data is received without an error message.

#### Exiting the interval printing mode

• Press the SET key to exit the interval printing.



When the interval printing is in hold status, press the SET key to exit this mode.

# When interval printing is terminated in the normal mode (character printing)

When interval printing is terminated in the normal mode, \*INTERVAL OFF\* is printed.



To restart interval printing, hold down the SET key and press the FEED key, then set a new interval again.



### Using the Interval Printing in the Chart Mode

Hold down the SET key and press the STAT. key.
 \*CHART MODE\* is then printed. Now the printer is in the chart mode.
 See Table-6 for the key functions used to set the target digits of chart generating.

Table-6,

Target digit	Key	Print
10¹ -10°	CL	6543210
10 <sup>2</sup> -10 <sup>1</sup>	SET	6543210
10 <sup>3</sup> -10 <sup>2</sup>	RSLT.	6543210
10 <sup>4</sup> -10 <sup>3</sup>	STAT.	6543210
10 <sup>5</sup> -10 <sup>4</sup>	DATA	6543210
10 <sup>6</sup> -10 <sup>5</sup>	FEED	6543210

(A target digit setting in chart printing is assigned to each key.

The target digits are indicated by cursor "^" in the "Print" column in Table-6.

Interval setting


Hold down the SET key and press the FEED key to enter the interval setting mode.


\*INTERVAL TIME\* is then printed.

The length of the interval can be determined by pressing any of the keys (see the table below for the correspondence between times and keys) to be pressed subsequently.

Table-7

Interval	Key
5 seconds	FEED
10 seconds	DATA
30 seconds	STAT.
1 minute	RSLT.
5 minutes	SET
10 minutes	CL
30 minutes	CE





### Print start command after interval setting


• When the interval is set, press the DATA key. Interval printing then begins.



If no data is received, the printer waits until data is received without an error message.

 In the chart mode, the starting time, first sample of data and target digits (i.e., WT83353.11g) are printed when data is received.
 See the printing example below.

### Printing example in the chart printing mode




### Pause in chart printing

◆ To pause chart printing, press the SET key. Once printing is paused (STOP), the stop time, final sample data, and chart generated target digits are marked with "^". See the printing example above.


The key assignment in the chart printing mode after pressing the SET key is described on the next page.

### Key assignment in the chart printing mode after pressing the SET key

- Press the DATA key to start printing again in the chart mode.
- Press the CL or CE key to terminate interval printing and print \*CHART OFF\*.
   For print-out in the chart printing mode again, set the interval again.
- Press the STAT. key to exit the interval printing and print \*CHART OFF\*.
   After this, you will be in the statistical calculation mode.



- Press the FEED key to feed one line of paper in the printer.
- The SET key can be used with the FEED, STAT., or RSLT. key to set the interval, target digits of chart printing, and date/timer.




- The printer does not function when the RSLT. key is pressed.
- The printer outputs \*NO DATA\* when data is not received after chart printing is paused.

# Section E

This section describes the statistical calculation function.

Statistical Calculation Mode	E•2
Selection of Standard Mode or	
Statistical Calculation Mode	E•2
Statistical Calculation Printing Mode	F•3





## Statistical Calculation Mode



No statistical calculations will be available in MODE 3.

The AD-8121 provides the result of statistical calculation as well as printing normal weight, percentage or counting data.


Weight or quantity/counting can be selected for calculation.

Additilnally, time can be printed for each sample of data.



### Selection of Standard Mode or Statistical Calculation Mode

 To select the standard mode or statistical calculation mode (for weight or quantity/counting), press the STAT. key.
 The selected LED on the printer goes on.



In the date print mode, the date will be printed when the printer is set from the standard mode to the weighing statistical calculation mode.



If previous data remains in memory, the mode is held even if the STAT. key is pressed. In this case, press the CL key to delete the remaining statistical data. When the data is deleted, \*CLEAR\* is printed. (See the printing example on the next page.)

## Statistical Calculation Printing Mode

- If previous data remains in memory, the mode is held even if the STAT. key is pressed. In this case, press the CL key to delete the remaining statistical data. When the data is deleted, the printer will print \*CLEAR\*..... (1)
- When data is input, a sample number is added automatically and printed. For the "Time print" setting, the time will be printed after the sample number. (Note that the date will not be printed for each item of data.) ···· (2)
- If an error occurs, the entry may be deleted by pressing the CE key. \*CANCEL\* is then printed. The last data is deleted from the statistical calculation. .... (3)
- If the position of the decimal point changes or if there is a change in the weight unit (e.g. ib, oz, ozt, g, kg, t, dwt, ct, mm, TL, GN, %, PC etc.) such data will be rejected from statistical calculations. .... (4)
- Press the RSLT. key to obtain the calculation results. The number of data items (N) and total (TOTAL) are then printed. .... (5) For the "Time print" setting, the time will be printed on the line preceding that of the number of data items (N).
- Press the RSLT. key again to print the maximum data (MAX), minimum value (MIN), average data (X), difference between maximum and minimum data (R), and standard deviation ( $\sigma$ ). .... (6)
- To continue calculation, input the next data....(7)
- The maximum number of data entries possible is 999. If 1,000 data entries are made the printer will automatically make a full print-out of all the statistical information for the preceding 999 entries, clear its memory and then start again by printing the 1,000th entry as entry number one.

Printing Example -

(1)	*CLEA	.R*		
(2)	NO. WT	1	10:09:52 178.632	
	10:10:0	741		y
(4)	WT	,,,,,	22.481	g
( ' '	NO.	2	10:10:30	
	WT		178.668	
	NO.	3	10:10:31	Ă
	WT		178.654	g
	10:10:5			
(4)	WT		178.6537	
	NO.	4	10:11:09	
(0)	WT		178.596	g
(3)	*CANC			۸ .
	WT	4	10:11:44 178.640	Α
		5	10:11:59	<b>9</b> Δ
	WT	•	178.599	
(5)	DATE	91/	11/03	
	N TOTAL		5	
	101112		893.193	g
(7)	NO.	6	10:12:30	Α
	WT		178.623	g
	NO.	7		
	WT		178.647	g
	DATE	91/	11/03	
	N TOTAL		7	
	TOTAL			
	TOTAL		1250.463	α
(6)			1250.463 178.663	g q
(6)	MAX MIN		1250.463 178.663 178.599	g
(6)	MAX		178.663	g g


σ R 0.0224 g

0.069 g

# Section F

This section describes dump printing and the AD-8121 specifications.

Dump Print Mode (MODE 3)	F•2
Printing in MODE 3	
On-line and Off-line	
Other Keys	F•2
Specifications	
General Specifications	
Input Specifications	
Statistical Calculation Specifications	F•4





# **Dump Print Mode (MODE 3)**



### Printing in MODE 3



Set DIP switch No. 2 ON.

 In MODE 3, the received data is printed as it is, no statistical calculations will be available in MODE 3.

The data consists of ASCII code 20 (H) to 7F (H) and characters shown to the right can be printed. (Comma (,) and 2C (H) cannot be printed-out.)

The specified code is processed as follows.

1B (H) 44 (H) Prints the date.

1B (H) 54 (H) Prints the time.

Up to 16 characters can be printed per line. For 17 or more characters, a line is generated after 16 characters are printed.

Set an interval of 1.1 seconds or more between the first sample of data (consisting of one line) and the next sample of data.

\*\* MULTI \*\*

\*\* FUNCTION \*\*

\*\* PRINTER \*\*

\*\* AD-8121 \*\*

\*\* CHARACTER \*\*

!"#\$%%'()\*+ --/

9123456789:;<=>?


@ABCDEFGHIJKLMNO
PQRSTUVWXYZ[\]^\_

'abcdef9hiJklmno
Parstuvwxyz(!)~"



### On-line and Off-line

In the dump print mode (MODE 3), the received data is printed as it is. The printer can also be set for no printing even when data is received.



Printing and no printing are alternately selected each time the DATA key is pressed in MODE 3. "Online" denotes printable status when data is received, and "offline" denotes non-printable status when data is received.

Table-8

MODE 3 (DIP switch No. 2 ON)

Online LED on, green data.
 Offline Does not print received data.

The green LED on the printer goes on or off, accordingly.

### **四** Other Keys

● In MODE 3, paper can be fed by pressing the FEED or DATA key. All keys except the FEED and the DATA keys are disabled.

(The EXT.P key functions the same as the DATA key.)

G:WT/R:QT

ON LINE

G:WT/R:QT

ON LINE

LED on

(Green)

LED off



## **General Specifications**

Model .	AD-8121
Printer Type	Mechanical dot impact printer
Character Size	5 x 7 dots 2.5 (H) x 1.8 (W) mm 0.1(H) x 0.7 (W) inches
Printing Speed	1 line/sec.
Paper Feed Speed	1 line/sec.
Number of Character	16 characters/line
Power Supply	AC/DC Adaptor (9V 700mA)/ Alkaline batteries (Approx. 7.000 lines)
Lithium Battery Life	Approx. 5 years
Printer Head Life	Approx. 500,000 lines
Operating Temperature Range	0°C to 40°C / 14°F to 104°F
Operating Humidity Range	80%RH or less (No condensation)
Storage TemperatureRange	-10°C to 50°C / 14°F to 122°F
Data Input Section	RS-232C or current loop
External Control	Printing only
Dimensions	180 (W) x 160 (D) x 80.5 (H) mm 7.1 (W) x 6.3 (D) x 3.2 (H) inches
Weight	Approx. 420g / 0.92lb (Without paper or batteries)


Recording Paper	PP143
Paper Quality	Fine
Paper Thickness	0.07 mm
Paper Width	44.5 mm
Length	Approx. 50 m
Maximum RolledDiameter	ø 65mm

## ==

## **Input Specifications**

Communication Method	RS-232C or current loop
Baud Rate	2400 bps
Data Bit	7 bits
Parity Bit	1 bit (EVEN)
Stop Bit	1 bit
Code	ASCII
Input Connector	XG4C-T0003 (Omron)

Pin No	Signal
3	RXD (Received data)
4	To pin 5
5	To pin 4
7	SG (Signal ground)
9	Current loop
10	Current loop



AD-8121 INPUT CONNECTOR



## Statistical Calculation Specifications

N: Number of data (999 samples Max.)

TOTAL: Cumulative total

MAX: Maximum data

MIN : Minimum data

X: Average data

 $\left( \sigma = \sqrt{\frac{N \cdot \Sigma(x \, i)^2 - (\Sigma x \, i)^2}{N \, (N - 1)}} \right)$ 

 $\sigma$ : Standard deviation

R: Range of data (Difference between

maximum and minimum data)

• Specifications are subject to change for improvement without notice.



A&D Company, Limited
3-23-14 Higashi-Ikebukuro, Toshima-ku, Tokyo 170 Japan
Telephone: [81] (03) 5391-6132 Fax: [81] (03) 5391-6148 Telex: 2422816 AANDD J

A&D ENGINEERING, INC. 1555 McCandless Drive, Milpitas, CA. 95035 U.S.A. Telephone: [1] (408) 263-5333 Fax: [1] (408) 263-0119

A&D INSTRUMENTS LTD. Abingdon Science Park, Abingdon, Oxford OX14 3YS England Telephone: [44] (0235) 550420 Fax: [44] (0235) 550485

A&D MERCURY PTY. LTD. 32 Dew Street, Thebarton, South Australia 5031 Australia Telephone: [61] (08) 352-3033 Fax: [61] (08) 352-7409

A&D KOREA Limited
3rd Floor Hanam Bldg 44-27 Yoldo-dong Youngdeungpo-ku Seoul, Korea
Telephone: (82) (02) 784-4264 Fax: [82] (02) 784-6557