

Additional Information (FZ-*i*/*i*WP/CT/GD series)

Contents

This manual describes additional functions for the FZ-*i* series, FZ-*i*WP series, FZ-CT series and FZ-GD series.

- Date and time output The Date and time can be added to the weighing data format.
- The Function switch..... The description of the switches to protect parameters.
- Correcting the internal mass value using an external calibration weight value. The correction method can be selected from either the correction 1 or correction 2.
 - The internal mass Calibrating the balance with internal mass, measuring an external calibration weight, calculating the correction value, and storing it in the balance.
 - The internal mass Calibrating the balance with an external calibration weight, correcting and storing the internal mass value with automatic calibration using the internal mass.

Caution

- The FX-*i* series, FX-*i*WP series, FX-CT series and FX-GD series do not have an internal mass and the date and time function.
- In the factory settings, the internal mass correction function can not be changed to protect against unintentional operations. When using the function, change the function switch to "Permit". Refer to "The Function Switch" on the next page or "Permit Or Inhibit" of the instruction manual.

Date and Time Output

With the parameters of "Time/Date output (*S-td*)" in "Data output (*dout*)", **Date** and **time** can be added to the data format of weighing data.


The Function Table For Output Of Date And Time

<i>dout</i> Data out	<i>S-td</i> Time/Date output	0	No output (Factory settings)
		1	Time only
		2	Date only
		3	Time and date

The Data Format Of Date And Time


Date *dout S-td 2 or 3*

This format can be changed in "Clock (Refer to *CL Adj*)". The year is output in a four-digit format.



Time *dout S-td 1 or 3*

This format outputs time in 24-hour format.



The Function Switch

The balance has function switches to protect the parameters against unintentional operations. When using the internal mass correction function, change each function switch to "Permit".

Procedure



Step 1 Press the **ON:OFF** key to turn off the display.

Step 2 While pressing and holding the **PRINT** key and the **SAMPLE** key, press the **ON:OFF** key to display **P5**.

Step 3 Press the **PRINT** key. Then the balance displays the function switches.

Step 4 Set the switches using the following keys.

- SAMPLE** key To select a switch to change the parameter. The selected switch blinks.
- RE-ZERO** key To change the parameter of the switch selected.
 - 0** To inhibit the change.
 - 1** To permit the change.
- PRINT** key To store the new parameter and return to the weighing mode.
- CAL** key To cancel the operation and return to the weighing mode.



The Function Table

- 0** To inhibit the change.
- 1** To permit the change. (Factory settings)

Internal Mass Correction

- 0** To inhibit the correction. (Factory settings)
- 1** To permit the correction.


The Internal Mass Correction 1

- The correction 1 corrects the internal mass using an external calibration weight. Calibrate the balance with the internal mass. Measure an external calibration weight. Calculate the correction value. And store it in the balance.
- The correction range of the internal mass value is as follows:

Model	External calibration weight	Allowance range
FZ-300CT	50g	±0.050g
FZ-600CT FZ-120i FZ-120iWP FZ-120GD	100g	
FZ-200i FZ-200iWP FZ-200GD FZ-300i FZ-300iWP FZ-300GD	200g	
FZ-1200i FZ-1200iWP FZ-1200GD	1kg	±0.50g
FZ-2000i FZ-2000iWP FZ-2000GD FZ-3000i FZ-3000iWP FZ-3000GD	2kg	


Procedure

- Example : Correcting the weight value by +0.003g in 100g using the FZ-120iWP. If correcting the weight value by +0.003g in 50g, the weight is changed into 100g, and the correction value is +0.006g.



- Step 1 Perform one-touch calibration. Measure an external calibration weight. Calculate the correction value.
- Step 2 Change the function switch of "Internal Mass Correction" and "The Function Table" to "Permit".
- Step 3 Press and hold the **[SAMPLE]** key to enter the function table and release the key when **bR5FnC** is displayed.
- Step 4 Press the **[SAMPLE]** key several times until **[5 in 1]** is displayed.
- Step 5 Press the **[PRINT]** key to correct the internal mass value using the following keys.

[RE-ZERO] key To increase the value.	Value is displayed -50 digits after +50 digits.
[MODE] key To decrease the value.	
[PRINT] key To store the new parameter and return to the weighing mode.	
[CAL] key To cancel the operation and return to the weighing mode.	
- Step 6 Press the **[CAL]** key to return the weighing mode.
- Step 7 Press the **[CAL]** key to perform the one-touch calibration.
- Step 8 Measure an external calibration weight and confirm that the weighing error is within the specified value in the "Specifications" of the instruction manual. If not, perform the internal mass correction and calibration.



The Internal Mass Correction 2

- The correction 2 corrects the internal mass using an external weight.
Calibrate the balance with an external calibration weight. Correct and store the internal mass value with automatic calibration using the internal mass.
- The corrected mass value is maintained in memory even if the AC adapter is removed.
The target range and allowance of an external weight is as follows:

Model	External calibration weight	Allowance range
FZ-300CT	50g* 20g	-0.050g ~ +0.050g
FZ-600CT	100g* 50g	
FZ-120i FZ-120iWP FZ-120GD	200g* 100g 50g	
FZ-300i FZ-300iWP FZ-300GD	300g 200g* 100g 50g	
FZ-1200i FZ-1200iWP FZ-1200GD	1000g* 500g	-0.50g ~ +0.50g
FZ-2000i FZ-2000iWP FZ-2000GD	2000g* 1000g 500g	
FZ-3000i FZ-3000iWP FZ-3000GD	3000g 2000g* 1000g 500g	

* Factory settings

Procedure

- Step 1 Change the function switch of "Internal Mass Correction" and "The Function Table" to "Permit".
- Step 2 Press and hold the **[SAMPLE]** key in the weighing mode to enter the function table and release the key when **bR5FnC** is displayed.
- Step 3 Press the **[SAMPLE]** key several times until **[5 in 2]** is displayed. If **[5 in 2]** is not displayed, proceed to step 1.
- Step 4 Press the **[PRINT]** key to display **[RL 0]**.
Perform the calibration using an external calibration weight. Refer to "Calibration Using An External Weight" of the instruction manual.
- Step 5 When displaying **rERdy ?** after the calibration, Remove the weight on the weighing pan and attach the breeze break. Press the **[PRINT]** key.
- Step 6 When displaying **[RL SET]**, the balance corrects and stores the new internal mass value.
- Step 7 When displaying **[RL in]**, the balance performs calibration using the internal mass automatically.
- Step 8 The balance displays **End** and returns to the weighing mode automatically.
- Step 9 Measure an external calibration weight and confirm that the weighing error is within the specified value in the "Specifications" of the instruction manual. If not, perform the internal mass correction and calibration.

