

### INSTRUCTION MANUAL

Instruction-HF/HF-G/HR-03/05-v.1.a-94.10.03

### Serial Interface

HF-03

HF-05

HR-03

HR-05



C-Param	eter Keys and Displays	-2
(	C-Parameters Settings	-3
i	Error Codes	-6
Serial Int	erface (OP-03)	-7
!	Installation	-7
;	Specifications	-8
(	Connections to other Equipment	-9
(	Connection to an AD-8121	-9
I	Data Output	10
I	Data Format	11
I	Data Format Examples	12
	Command list	14
,	Commands to request Weighing Data	15
	Commands to control the Balance	15
1	Commands to set Data	16
1	Commands to request Data	17
ı	Command Examples	19
Current l	Loop Interface (OP-05)	21
	Installation	21
	Circuit Diagram and Connector Drawing	21


### **C-Parameter Keys and Displays**

0	)			

This mark appears when the memorized parameter is shown in the display.



The MODE key is used to select the group of C-parameters.



The SAMPLE key or RANGE key is used to select the item from the group selected by the mode key.



The RE-ZERO key is used to select a parameter for the item selected by the MODE and SAMPLE key (or RANGE key).




The PRINT key is used to save the new C-Parameter settings and to exit to the weighing mode.



The ON:OFF key cancels the new C-parameter settings and turns the display off.

- Turn the display off.
- Press and hold the RE-ZERO key and press the ON:OFF key. Release both keys.
- Press the MODE key. The balance enters the function setting mode and [ [] will be displayed.
- Set the C-parameter using the keys described on the preceding page.




### **■** C-Parameters Settings

### [ - 4 Data out Settings for options 03 and 05

	parameter	It is necessary to sele If you use auto print m band".	ct the method and condition when data is transmitted. node, please set both "Auto print polarity" and "Auto print			
, ,	<i>⁺□</i>	PRINT key mode	PRINT key command is accepted only if the display is stable. The display will blink when data is transmitted.			
Prink Data out mode	1	Auto print mode A	Data is transmitted when the display is stable and meets the conditions of "Auto print polarity" and "Auto print band". This standard value is the zero-point. Refer to the figure at the bottom of the page.			
	2	Auto print mode B	Data is transmitted when display is stable and meets the condition of "Auto print polarity" and "Auto print band". This standard value is the last weighing data. Refer to the figure at the bottom of the page.			
	3	Streem mode	Data is transmitted continuously.			
   8P-P	parameter	This parameter sets weighing data is able	This parameter sets the polarity condition from the standard value that weighing data is able to transmit.			
Auto print polarity	<b>.</b> 0	Data is able to be transmitted when the weighing data is more than the standard value. (It is called "Plus polarity")				
	1	Data is able to be transmitted when the weighing data is less than standard value. (It is called "Minus polarity")				
	2	Data is able to be transmitted when the weighing data is less than or more than the standard value. (It is called "Absolute polarity")				
	parameter	This parameter sets the band width condition from the standard value that weighing data is able to transmit.				
RP-Ь	*0	Data is able to be tran value more than 10 dig	Data is able to be transmitted when the weighing data deviates from the standard value more than 10 digits.			
Auto print band	1	Data is able to be tran value more than 100 d	smitted when the weighing data deviates from the standard ligits.			
	2	Data is able to be transmitted when the weighing data deviates from the standard value more than 1000 digits.				

<sup>\*</sup> Factory setting.

#### **AUTO PRINT MODE:** O Stable data and print



### [-4 Continued

PAUSE		Selects the use of a pause. If your printer requires a pause, set this to 1.
Data pause	*0	No pause.
		Using a pause.
AL-F		Selects whether or not to use Auto feed after printing. (AD-8121A)
Auto feed	*D	Not using auto feed
	1	Using auto feed
Ar-d		Selects whether to use auto re-zero after transmitting.
Automatically	*0	Not using re-zero
re-zero after data out	- 1	Using re-zero
inFo	:	Selects whether you verify the calibration using printer, a computer or no verification after calibration.
Verifying the calibration	•0	No verification of the calibration
		Verifying the calibration using the AD-8121 printer.
	2	Verifying the calibration using a computer.

<sup>\*</sup> Factory setting.

### [ - 5 Serial Interface Settings for options 03 and 05.

		Parameter definition and use.
. 00		600 baud
<i>⊾PS</i>	1	1200 baud
Baud rate	* 2	2400 baud
	3	4800 baud
	4	9600 baud
		Parameter definition and use.
bE-Pr	* 0	7 bits, Even parity check
Parity bit	1	7 bits, Odd parity check
	2	8 bits, no parity check
[r-LF		Parameter definition and use. (Common setting to both transmitting and receiving.)
Terminator	* 0	C <sub>R</sub> , L <sub>F</sub>
	1	C <sub>R</sub>
		Selects a weighing data format. Refer to the book for the option used.
	* 0	A&D standard
ESPE	1	Dump print format
Data format	2	KF format
	3	MT format
		Selects the maximum waiting time between receiving command characters.
E-UP Receive time	0	NO limit
	* 1	One second.  If the time is over, the balance canceles receiving the command and outputs the error code.
E-Lad Error code	•0	Error codes or <ak> signal are not output.</ak>
Euroi code	1	Error codes or <ak> signal are output.</ak>
CE5		Selects the use of the control lines CTS and RTS.
CTS control	*0	Not using CTS or RTS.
	1	Using CTS and RTS. Keep the RTS line (active) high while the computer receive data. The balance will set CTS low if it is busy.

<sup>\*</sup> Factory setting.

## **Error Codes**


Error code	Description of the error			
E00	Communications error There is a protocol error in communications. Check the format, baud rate and parity.			
E01	Undefined error The command received was not in the list of comm this balance.	ands recognized by		
E02	Balance not ready The command can not be received due to a timing is busy. Commands that the balance can not resp when the balance is in the counting mode.	error or the balance ond to such as "Q",		
E03	Time over error The balance did not receive next character of a complimit of one second.	mand within the time		
E04	Excess characters error  The command has more characters than is required or the range of the data is beyond what the balance will accept. Example; when the calibration weight entered is greater than the range of the balance.			
E05	Terminator error A command is followed by other than CR or CR, LF (carriage return and a line feed).			
E06	Format error A command that should include numerical data has none or the data is in the wrong place in the command.			
E07	Out of range error The data entered exceeds the range that the balan	nce will accept.		
E10	Internal operation error The balance is in an abnormal operating condition.	Error D		
E11	Stability error The balance can not stabilize due to vibration or other environmental problem.	Error 1		
E20	Calibration error The calibration weight is to heavy.	[RL E		
E21	Calibration error The calibration weight is to light.	-[RL E		
E22	Zero out of range error The balance can not zero the display as the zero point offset exceeds the range.			

# Serial Interface (OP-03)


Option 03, is a serial interface that consists of two sections, a bi-directional EIA RS-232C interface and a 20mA passive current loop.

### **Installation**


Turn off the balance and remove the AC adapter.




Remove two screws securing the blank cover on the rear of the balance.



Remove the cable attached to this panel and connect it to the socket on the option board as shown.



Install the option board using the screws removed in step 1.



### **Specifications**

Transmission system: EIA RS-232C, 20mA current loop (passive)

Transmission form

: Asynchronous, bi-directional, half duplex

Data format

: Baud rate : 600, 1200, 2400, 4800, 9600 bps

Data

7 or 8 bits

Parity

Even, Odd (7 bit)

None (8 bit)

Stop bit

: 1 or 2 bits : ASCII

Code

DATA	RS-232C levels	Current loop	
1	- 5V to - 15V	20mA	
0	+ 5V to + 15V	0mA	

#### Pin connections

#### RS-232C

Pin No.	Signal name	Direction	Description
1	FG	-	Frame ground
2	RxD	Input	Receive data
3	TxD	Output	Transmit data
4	RTS	Input	Ready to send
5	CTS	Output	Clear to send
6	DSR	Output	Data set ready
7	GND	-	Signal ground
8 - 25	n/a	-	-

#### **Current loop**


Pin No.	Signal name	
. 1	N.C.	
2	N.C.	
3	Loop	
4	N.C.	
5	Loop	
CASE	Frame GND	



13 12 11 10 9 8 7 6 5 4 3 2 1 \0.00000000000000

25 24 23 22 21 20 19 18 17 16 15 14

Current loop specifications Maximum voltage-----25V



### **Connection to other Equipment**

- The current loop is of the passive type. It requires an external source of 20mA DC.
- The RS-232C is of the DCE type (Data Communications Equipment) and can use standard DCE cables.
- When connecting to another piece of equipment, consult the manual for that equipment for proper settings and connections.

### Connection to an AD-8121

The following balance functions must be set to use the AD-8121 printer

"C" function	Settings
C-4 Print 0, 1, 2, 3 C-4 AP-P 0, 1, 2 C-4 AP-B 0, 1, 2 C-5 BPS 2 C-5 BE-Pr 0 C-5 Cr-LF 0 C-5 E4PE 0, 1 C-5 E-Cod 0 C-5 CE5 0	Select a print mode Select the polarity for the auto-print mode Set the auto-print band Select "2400bps" Select "7 bits, Even parity check" Select "CR, LF" Select "A&D Standard format" or "Dump print format" Set the receive timing to 1 second Select "Error codes are not output" Select "Not using CTS and RTS"

### **■** Data Output

There are four modes to control the transmission of the weighing data.

#### **Key Mode**

When you press the PRINT key, the balance transmits the weighing data when the display is stable (the stability indicator is on). When the data is transmitted the display will blink one time.

#### Auto-print Mode A

The balance transmits the weighing data when the display is stable (the stability indicator is on), meets the conditions of "Auto-print polarity" and "Auto-print band". The reference for the auto-print band is the zero point. When the data is transmitted the display will blink one time.

#### **Auto-print Mode B**

The balance transmits the weighing data when the display is stable (the stability indicator is on), meets the conditions of "Auto-print polarity" and "Auto-print band". The reference for the auto-print band is the last weighing data printed. When the data is transmitted the display will blink one time.

C-4	Pr int	2		Auto-print mode B
C-4	AP-P	X	***************************************	Auto-print polarity $X = 0$ , 1, 2
C-4	RР-Ь	Χ	• • • • • • • • • • • • • • • • • • • •	Auto-print band X = [] / 2

#### Stream Mode

The balance transmits the weighing data continuously.

C-1 
$$5PEEdX$$
 ...... Display update rate.  $X = 0$ ,  $I$ ,  $Z$  C-5  $BPS$   $X$  ..... Baud rate.  $X = 0$ ,  $I$ ,  $Z$ ,  $Z$ ,  $Y$ 

#### NOTE:


When the baud rate is set to 600 or 1200bps and the refresh rate of the display is set to high speed, the balance is unable to transmit the data completely.

### Data Format

There are four formats for transmission of the weighing data. The setting of *[-5 LYPE* selects the data format.

#### 


This format is used when the peripheral equipment is capable of receiving A&D format. If an AD-8121 is to be used and you are sending a 15 character data string (no excluding the terminator), set the printer to mode 1 or 2.



- · A two character header indicates the status of the stability.
- The weighing data (with leading zeros) plus sign and decimal point, followed by a three character "unit of weight" make up the body of the data.
- A terminator consisting of  $C_R$ ,  $L_F$  to indicate to the peripheral equipment that all of the data has been sent.
- Header: Stable header is ST, Stable header for counting mode is QT
   Unstable header is US
   Overload header is OL

#### Dump Print Format EYPE |

This format is used when the peripheral equipment is not capable of receiving A&D format. If an AD-8121 is to be used and you are sending a 16 character data string (no excluding the terminator), set the printer to mode 3.




- A two character header indicates the status of the stability if not overloaded or the display is zero.
- The weighing data (with leading zeros replaced by spaces) plus sign and decimal point, followed by a three character "unit of weight" make up the body of the data.
- A terminator consisting of  $C_R$ ,  $L_F$  to indicate to the peripheral equipment that all of the data has been sent.
- Header: Stable header is WT, Stable header for counting mode is QT
   Unstable header is US

#### KF Format

#### F 7PE

This is the Karl-Fischer moisture meter format and is used when the peripheral equipment can not communicate using A&D format.




- The data consists of 14 characters (no excluding the terminator).
- The sign of the weighing data is first if the balance is not in overload. The sign is omitted if the balance is at zero.
- · The sign is followed by the weighing data (with leading zeros replaced by spaces) and decimal point. The weight data is followed by the unit if the balance is stable.
- A terminator consisting of  $C_{\rm R}, L_{\rm F}$  to indicate to the peripheral equipment that all of the data has been sent.
- · The unit presents at stable. The unit does not present at unstable.

#### MT Format

#### FABE

· The data length will be changed by the unit or overload.




- · The weighing data is proceeded by a header of two characters. If stable, one character and a space are transmitted.
- · The minus sign will be next if the weighing data is negative. The sign is omitted if the weighing data is positive or at zero. Leading zeros are replaced by spaces.
- If the balance is in overload, the weighing data is omitted.
- Header: Stable header is S .......

Unstable header is SDD

### **Format**

Data format examples use HF-30(

STABLE A&D | S 0.000g D.P. 0 0 0 KF MT **UN-STABLE** A&D|U -832 IOs D.P. I KF MT



UNITS (HF -G series can not change the unit and HR series can use without animal weighing mode)

Unit and display	sign /	A&D	D.P.	KF	MT
g	g	uu g	u g	uguu	<u>_</u> g
mg mode	mg	∟ m g	∟mg	_mg_	_ m g
Counting mode	pcs	⊔ P C	□ P C	□ p c s	∟PCS
Precent mode	%	니니%	ㅁㅁ%	_%	%
Ounce (Avoir)	<b>UZ</b>	0 z	0 Z	⊔ O Z ⊔	_ 0 Z
Pound	lb	1 b	1 b	_ 1 b _	_ 1 b
Pound Ounce	_ DZ	_ 0 Z	O Z	_ O Z _	_ O Z
Troy Ounce	ozt	o z t	ozt	u ozt	u ozt
Metric Carat	ct	_ ct	ct	_ ct_	_ct
Momme	mo <b>m</b>	m o m	m o m	⊔ m o m	∟ m o
Pennyweight	dwt	d w t	dwt	u d w t	∟d w t
Grain	GN	□GN	□GN	_gr_	∟GN
Tael (HK general, Sing.)	ti	_ t 1	t 1	_ t 1 s	_ t 1
Tael (HK, jewelry)	tı	t 1	t 1	_ t 1 h	_ t 1
Tael (China)	ti	t 1	t 1	_ t 1 t	_ t 1
Tael (Taiwan)	ti	t 1	t 1	L t 1 C	_ t 1
Tola (India)	t		பட t	_ t o 1	_ t
Messghal	·M	m e s	mes	<u></u> МS	∟ m
Animal mode	â	L L g	u u g	шдшш	<u>ு</u> g
•					

<sup>□</sup> Space, ASCII 20H

<sup>&</sup>lt;sup>С</sup>в Carriage Return, ASCII 0DH

L<sub>F</sub> Line Feed, ASCII 0AH

### Command list

Command examples use HF-300.

#### **Command list**

Commands to	request weighing data Description of the command				
Q	Request for weighing data immediately				
SI	Request for weighing data immediately				
S	Request for weighing data when stable				
SIR	Request for weighing data continuously				
С	Clear the SIR command				
	control the balance				
Р	Display ON/OFF (same as ON:OFF key)				
ON	Display ON				
OFF	Display OFF				
R	Display RE-ZERO (same as RE-ZERO key)				
CAL	Perform calibration				
U <sup>∙</sup> ²	Change unit of weigh (same as MODE key)				
PRT	Print (same as PRINT key)				
RNG "	Change Range (same as SAMPLE key or RANGE key)				
SMP *2	Enter the sample weight in counting mode and percent mod				
(same as SAMPLE) key or RANGE key)					
Commands to	set values				
CW:	Set the calibration weight				
	Set the digital tare weight				
PT:	Set the digital tare weight				
ID:	Set the digital tare weight Set the ID number				
ID: UW: *²					
ID:	Set the ID number				
ID: UW: *²	Set the ID number Set the unit weight for counting mode				
ID: UW: " <sup>2</sup> PW: " <sup>2</sup> FC:	Set the ID number  Set the unit weight for counting mode  Set the 100% weight for percent mode  Set a function parameter				
ID: UW: "2 PW: "2 FC: Commands to	Set the ID number Set the unit weight for counting mode Set the 100% weight for percent mode Set a function parameter  request Data				
ID: UW: '2 PW: '2 FC: Commands to	Set the ID number  Set the unit weight for counting mode  Set the 100% weight for percent mode  Set a function parameter  request Data  Request for the calibration weight				
ID: UW: "2 PW: "2 FC: Commands to ?CW ?PT	Set the ID number  Set the unit weight for counting mode  Set the 100% weight for percent mode  Set a function parameter  Prequest Data  Request for the calibration weight  Request for the digital tare weight				
ID: UW: "2 PW: "2 FC:  Commands to ?CW ?PT ?ID	Set the ID number  Set the unit weight for counting mode  Set the 100% weight for percent mode  Set a function parameter  request Data  Request for the calibration weight  Request for the digital tare weight  Request for the ID number				
ID: UW: "2 PW: "2 FC:  Commands to ?CW ?PT ?ID ?UW "2	Set the ID number  Set the unit weight for counting mode  Set the 100% weight for percent mode  Set a function parameter  request Data  Request for the calibration weight  Request for the digital tare weight  Request for the ID number  Request for the unit weight for counting mode				
ID: UW: "2 PW: "2 FC:  Commands to ?CW ?PT ?ID ?UW "2 ?PW "2	Set the ID number  Set the unit weight for counting mode  Set the 100% weight for percent mode  Set a function parameter  Prequest Data  Request for the calibration weight  Request for the digital tare weight  Request for the ID number  Request for the unit weight for counting mode  Request for the 100% weight for percent mode				
ID: UW: "2 PW: "2 FC:  Commands to ?CW ?PT ?ID ?UW "2 ?PW "2 ?FC	Set the ID number  Set the unit weight for counting mode  Set the 100% weight for percent mode  Set a function parameter  Prequest Data  Request for the calibration weight  Request for the digital tare weight  Request for the ID number  Request for the unit weight for counting mode  Request for the 100% weight for percent mode  Request for a function parameter				
ID: UW: "2 PW: "2 FC:  Commands to ?CW ?PT ?ID ?UW "2 ?PW "2 ?FC ?UT "2	Set the ID number  Set the unit weight for counting mode  Set the 100% weight for percent mode  Set a function parameter  Prequest Data  Request for the calibration weight  Request for the digital tare weight  Request for the ID number  Request for the unit weight for counting mode  Request for the 100% weight for percent mode  Request for a function parameter  Request the current unit of weight				
ID: UW: "2 PW: "2 FC:  Commands to ?CW ?PT ?ID ?UW "2 ?PW "2 ?FC	Set the ID number  Set the unit weight for counting mode  Set the 100% weight for percent mode  Set a function parameter  Prequest Data  Request for the calibration weight  Request for the digital tare weight  Request for the ID number  Request for the unit weight for counting mode  Request for the 100% weight for percent mode  Request for a function parameter				

HF-G and HR series use RANGE key and HF series uses SAMPLE key.

<sup>&</sup>lt;sup>12</sup> HF-G series can not use.

-		
1	Com	mands to request Weighing Data
	Q	Query for weighing data (the balance will respond with the weighing data immediately)  Command QCRLF  Reply ST, +0012.783gCRLF
	SI	Send the weighing data immediately (same as Q) (the balance will respond with the weighing data immediately)  Command  SICRLF  Reply  US, +0012.783gC_RLF
	S	Send the weighing data when it is stable (the balance display will blink when the data is transmitted)  Command SCRLF  Reply ST, +0027.835 LgCRLF
	SIR	Send the weighing data continuously (the balance sends the data in stream mode)  Command  SIRCRLF  Reply  US, +0027.835 L g CRLF
	C	$\begin{array}{c ccccccccccccccccccccccccccccccccccc$
3	Com	mands to control the Balance
_	P	Display ON/OFF (same as pressing the ON/OFF key) (if the balance is on it will turn off, if it is off it will turn on )  Command PCRLF
	ON	Display ON (If the balance is OFF, it will turn ON)
	OFF	Command ONCRLF  Display OFF (If the balance is ON, it will turn OFF) (if the balance is already off, nothing will happen)  Command OFFCRLF
	R	RE-ZERO the balance (same as pressing the RE-ZERO key) (the balance display will zero)  Command RCRLF
	CAL	Preform calibration (the balance will enter the calibration mode)

Command

Ū	Change the mode (same as pressing the MODE key) (the balance will shift to the next selected unit of weigh, to the counting or percent mode) Command  UCRLF  (HF and HR series can use this command)
PRT	Print (same as pressing the PRINT key) (if the balance will send data dependent on the function parameters)  Command PRTC <sub>R</sub> L <sub>F</sub>
RNG	Change range (same as pressing the RANGE key or SAMPLE key)  Command RNGCRLF
SMP	Enter the sample weight in counting mode and percent mode (same as pressing the RANGE key or SAMPLE key)  Command SMPCBLE

### Commands to set Data

CW: Command CW: 0300.034 L g CR LF

When in the calibration mode, this command is used to set the actual value of the calibration weight. The calibration weight is reset to the standard value upon entry into the calibration mode.

Caution: only grams can be used as the unit of weight when setting the value. Please set the calibration weight after entering the calibration mode.

PT: Command PT: 0045.670 \_\_\_\_ g C\_RL\_F

This command is used to set the digital tare weight.

Caution: pressing the RE-ZERO key or commanding RE-ZERO will clear this value.

ID: Command ID: 123-ABCCRLF

When using GLP (Good Laboratory Practice), this command is used to set the ID number to be printed out of the verification of calibration. The ID can consist of seven characters, A through F, numbers, 0 through 9, spaces or the negative sign (-).

UW: Command UW: 001.23400 g CRLF

When using counting mode, this command is used to set the unit weight. The command that delete the unit weight is as follows:

Command UW: C GREF (HF and HR series can use this command)

PW: Command PW: 034.560 U g CRLF

When using percent mode, this command is used to set the 100% weight. The command that delete the 100% weight is as follows:

Command PW: C GRLF (HF and HR series can use this command)

FC: Setting a "C" parameter

Group number

Item number

Parameter

Command

FC: 0 1; 3 CRLF

C-C Cond 3 Sets display stability to strong

FC: 5 0; 3 CRLF

FC: 5 3; 1 CRLF

Sets the baud rate to 4800bps

FC: 5 3; 1 CRLF

Sets the print format to D.P.

	-	· · · · · · · · · · · · · · · · · · ·	·		-			
Grou	ıp Number		Item a	and Item N	umber			
	Group	0	1	2	3	4	5	6
0	[-[] Environment	SEB-b Stability band	Lond Resp. / Environ.	Erc Zero tracking				·
1	[- ] Display	SPEEd Refresh rate	Par nE Decimal point	P-on Auto start function				
2	C-2 Auto re-zero	Arto re-zero an/off	Ar-b Auto re-zero band	Pr-L Detection time				
3	[-]	ERL Calibration inhibit						
4	C-4 Data out	Prink Data out mode	RP-P Auto print polarity	RP-b Auto print band	PRUSE Data pause	RL-F Auto feed	Rr-d Zero after data out	, nFo
5	E-5 Serial interface	LPS Baud rate	be-Pr	Er-LF Terminator	LYPE Data format	L-UP	E-Cod	CES control
6	C-6	Response/environment is common data with the condition of resonnse accessable using the keyboard. If a value is set in the C parameters, it will be changed if new conditions of response are set.						
7	[-7							
8	[-8	will be changed if new conditions of response are set.						
9	[-9 Parameter control	Pn ID protect	PF Parameter protect					

### Commands to request Data

? C W Request for the calibration weight.

'்"= a space

Command ? C W C<sub>R</sub> L<sub>F</sub>

Reply CW, +0300.034 \_\_\_\_g c<sub>R</sub>L<sub>F</sub>

Caution, if the calibration mode has been entered since the calibration weight was last entered, the value returned will be the standard calibration weight. The unit of weight will always be grams.


?PT Request for the digital tare weight.

Command ? PTC<sub>R</sub>L<sub>F</sub>

Reply PT, +0045.670 \_\_\_\_gC\_RL\_F


Caution, if the RE-ZERO key has been pressed since the last digital tare was registered, the value returned will be a positive value equivalent to the weight of any item on the pan prior to pressing the RE-ZERO key. The unit of weight will be the unit of weight currently in use.

The fact that the value returned after RE-ZERO has been pressed is the actual tare weight, means this can be used to read the tare weight into a program.



### **Command Examples**

### ON Command Command examples use HF-300. Computer HF series balance (HF-G and HR series) Time OFF state ON command Confirmation of AK CR LF command receipt Display test Command compleation Display zero





Please note that there needs to be a delay time of one second between the balance acknowledgment <AK> and the next command transmission to the balance. Line 124 sets the delay time. To change the delay, set "100" to another value. Example: 124 FOR I = 1 TO 200: NEXT I.

#### Example using a BASIC language program

1...

1...

LINE INPUT #1, AK\$ 123

124 FOR I = 1 TO 100: NEXT I

125 PRINT #1, "Q"


Receive <AK>

Delay time (the '100' sets the delay)

TX: 'Q' command

HF/HF-G/HR series Options v.1.a


#### OP-03 | CAL | Command (Calibration)



# Current Loop Interface (OP-05)

Option 05 is a 20mA current loop interface. As this option is similar to option 03, please refer to that section for characteristics and use.

This option is passive and requires an external current source to operate. It is a transmit only device and as such is best used with external displays or printers. Many such devices supply the current required by this interface.



### Installation


The installation of this option is the same as option 03, please refer to that section for the installation.

### Circuit Diagram and Connector Drawing

Pin	Discription			
1	No connection			
2	No connection			
3	Current loop (positive or negative)			
4	No connection			
5	Current loop (positive or negative)			
6	No connection			
7	No connection			
Shell	Frame ground			



7 Pin DIN Connector



Specifications
Maximum voltage-----25V