MS-74AT / MX-53AT

(This document is included in the packing contents of the product.)

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Please see the separate
instruction manual [here](#)

SAFETY NOTES AND OTHER NOTES

Word and symbol (for harm to persons / for damage to property):

 WARNING	If this notation is ignored and mishandled, it indicates that a person may be killed or seriously injured.
 CAUTION	If this notation is ignored and mishandled, it indicates that a situation where a person is injured or property damage is expected to occur.

Graphic and symbol:

	 This symbol indicates mandatory actions that must be followed. Specific mandatory actions are indicated by text or illustrations within or near  .
---	--

General message:

Caution	This is a description of important points for correct usage.
Note	'High possibility of inappropriate handling' or 'general advice on using the product'.
	Indicates high temperature caution.

About This Manual

- (1) No part of this manual may be reprinted, copied, modified, or translated to another language without the prior written consent of A&D Company, Limited (A&D).
- (2) The contents of this manual are subject to change without notice.
- (3) Please contact A&D if you notice any uncertainty, errors, omissions, etc. in this manual.
- (4) A&D bears no liability for any loss or lost profits due to the operation of this product, and for direct, indirect, special, or consequential damages resulting from any defect in this product or this manual, even if advised of the possibility of such damage. Furthermore, A&D assumes no liability for claims of rights from third parties. Concurrently, A&D assumes no liability whatsoever for software or data losses.

© 2024 A&D Company, Limited

- ☐ The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by A&D is under license.
- ☐ Other product names and company names mentioned in this manual are trademarks or registered trademarks of their respective companies in Japan or other countries and regions.

About this instruction manual

Thank you for purchasing an A&D MS-74AT or MX-53AT moisture analyzers. This manual describes the safety precautions, installation instructions, basic functions and main operations of the moisture analyzer. For more detailed contents, please refer to the separate instruction manual from the A&D website.

Safety Information!

This product is a moisture analyzer that heats and dries samples using a halogen lamp based on the principle of thermogravimetric analysis, evaporating moisture to measure moisture content, solid content, etc. from the mass change. Do not use it for any other purposes.

This manual contains basic information for handling this product. Please read and understand this manual thoroughly before use.

This product is designed to be used by skilled or experienced operators.

Modifying or disassembling this product, or using it in a manner not specified in this manual or the separate instruction manual, may compromise safety. A&D assumes no responsibility for any issues arising from such actions.

1. Handling precautions


1.1. Installation precautions

⚠ WARNING

--Fatal accidents caused by fires and explosions--

- **Measure in an appropriate environment.**
 - ❑ Do not use the product in a hazardous environment such as in a flammable gas atmosphere.
 - ❑ Strictly adhere the following installation environment conditions:
Temperature: 5 to 40°C, Humidity: 85%RH or less (without condensation)
 - ❑ Do not use the product in a poorly ventilated installation environment. Using the product in enclosed spaces where heat dissipation is obstructed can lead to unexpected temperature increases or poisoning symptoms.
- **⚠ Do not place flammable items around the product.**
 - ❑ The parts of the moisture analyzer become hot during and immediately after measurement. There is a risk of ignition.
 - ❑ Never place anything on top of the heater cover.
 - ❑ Do not operate the device while cables or other objects are hanging on it.


Heater cover

⚠ WARNING

--Risk of death or malfunction due to electric shock--

- **Check the power supply voltage.**
The appropriate power supply voltage is indicated on the back of the moisture analyzer.
Make sure that it matches the power supply voltage you are using.
- **Use the standard power cable provided when connecting the moisture analyzer to the power supply.**
- **Ground the moisture analyzer.**
- **Be careful of water exposure. The moisture analyzer is not waterproof. If water gets inside the product, it may cause electric shock or malfunction.**

Caution

--Precautions for accurate measurement--

- **The mass sensor uses a high-precision electromagnetic balance system. Please pay attention to the following points.**
 - ❑ Install on a sturdy table.
 - ❑ Avoid vibrations.
 - ❑ Install in a location that is not exposed to wind from air conditioners or similar sources.


1.2. Precautions for use

⚠ WARNING

-- Risk of death due to fire, explosion, or other toxic hazards --

- **Do not measure hazardous specimens.**
 - ❑ Do not heat samples that may explode or ignite, or heat samples that generate harmful substances, as doing so is extremely dangerous. Samples with unknown properties are likewise dangerous.
 - ❑ If the sample surface dries first and the internal pressure increases, it may burst. This is dangerous, so do not perform such measurement.
 - ❑ If the sample ignites, immediately unplug the product and take appropriate actions.
Note that the case of the moisture analyzer is made of flame-retardant material (UL94 V-0).
 - ❑ Do not place flammable items nearby.
 - ❑ The parts of the moisture analyzer become hot during and immediately after measurement.
Do not place flammable items nearby as doing so may cause a fire.
 - ❑ Never place anything on top of the heater cover.

Heater cover

⚠ WARNING

--Blindness and other accidents--


- Do not look directly at the halogen lamp during heating. Doing so may cause eye pain and vision problems.
- Wear appropriate protective equipment such as protective clothing, safety glasses, and protective gloves while measuring.
- Do not modify or disassemble the moisture analyzer. Doing so can cause malfunction, electric shock, or fire. If you suspect a malfunction, please consult the store where you purchased the product.
- Do not expose the moisture analyzer, power cable, or included accessories to extreme temperatures, strong chemical vapors, moisture, shock, vibration, or strong electromagnetic fields. Follow the operating conditions specified in “[7. Specifications](#)”

⚠ Caution


--⚠ Burns caused by hot surfaces--

- Follow the correct operating methods and procedures.


- ❑ Position the breeze break, tray holder, sample tray, and sample tray handle properly.
- ❑ Be sure to open and close the heater cover by holding the heater cover handle.

- ❑ During the measurement, the shaded area in the figure becomes hot. Only touch the heater cover handle and the sample tray handle during operation.

- ❑ Do not touch the inside of the heater cover or hot parts such as the sample tray directly with your hands immediately after heating. Doing so may result in burns. During and immediately after measurement, these parts are hot. When handling the moisture analyzer, use the designated handles and provided tools, and avoid touching the glass or metal parts near the halogen lamp (heater) in the glass housing.
- ❑ In particular, the glass housing part becomes hot and may cause burns if touched.

- ❑ Immediately after measurement, the sample tray and the part around the sample tray handle are hot. Cool it in an appropriate location. When handling the sample tray, use tweezers or similar tools.

CAUTION

--Injuries caused by improper handling--

- Dispose of used halogen lamp unaltered.
Breaking a halogen lamp can cause glass fragments to scatter and lead to injury.
- Be careful not to drop, bump, or scratch glass components, including halogen lamp.
They may break, leading to injury.
- The edges of the sample tray are sharp. Be careful not to cut yourself.
- Use the handle when opening and closing the heater cover. Be careful not to pinch your fingers.
- When moving the moisture analyzer, ensure that it is cool and lift it without tilting.
Moving it while tilted may cause the heater cover to open unexpectedly, leading to injury.

CAUTION

--Equipment failure due to improper handling--

- It is recommended to replace halogen lamp after their rated lifespan (approximately 5000 hours) as failure to do so may cause damage.
- Be careful not drop, bump, or scratch glass parts, including halogen lamp. They may break, leading to injury.
- Make sure that dust or water does not enter the inside of the moisture analyzer.
- Use only A&D accessories and attachments.
- If the moisture analyzer behaves unexpectedly, unplug the power and plug it in again.
If it does not work properly, please request a repair.

Caution

- **Precautions for high temperature heating**
 - ☐ Do not continuously heat at a temperature setting of 200°C for more than 30 minutes. The safety mechanism may activate and cut off the power supply to the halogen lamp.
 - ☐ For safety reasons, if heating at 160°C or higher continues for 1 hour, the maximum temperature will automatically be limited to 160°C.
 - ☐ During measurement, the STOP button on the touchscreen is always active. If any abnormalities or dangers are detected, immediately press the STOP button to halt the measurement.

1.3. Precautions after use, precautions for storage, etc.

CAUTION

-- Burns caused from hot surfaces--

- Make sure that each part of the moisture analyzer has cooled down sufficiently before performing any operations.
 - Especially when replacing the lamp, make sure the temperature around the glass housing has decreased before proceeding. For detailed lamp replacement instructions, refer to the separate 'Instruction Manual'.

--Equipment failure--

- Make sure that dust or water does not enter the inside of the moisture analyzer.
- When cleaning, use a cloth dampened with diluted neutral detergent. Do not use organic solvents or chemical wipes as they may cause damage.
- Always unplug the power cable from the outlet before performing maintenance.
- Use the dedicated packing box for transport.


Caution

- Clean the glass housing if it becomes dirty, as it may not heat properly if dirty. Refer to “[6.1. Cleaning the heater](#)” for handling instructions.
- Clean the halogen lamp if it has fingerprints or other dirt, as those can shorten its lifespan. Refer to the separate instruction manual for handling instructions.
- Do not apply shock or drop the moisture analyzer as doing so may damage the halogen lamp or mass sensor.
- If the product will not be used for a long period of time, unplug the power cable from the outlet.

2. Packing contents and part names

2.1. Part names


- Make sure all the accessories are included.
- Keep the packaging material for use during transport or repairs.

No.	Name
1	Heater cover
2	Heater cover handle
3	Display
4	Leveler
5	Keys
6	Sample tray handle
7	Foot adjuster
8	Fuse (T6.3A 250V)
9	Power inlet
10	RS-232C interface

No.	Name
11	USB interface
12	LAN interface
13	Dedicated USB Type A connector for memory stick
14	Disposable aluminum tray
15	Sample tray
16	Tray holder
17	Breeze break
18	Halogen lamp
19	Glass housing
20	Serial number

2.2. Accessories

No.	name	MS-74AT	MX-53AT	Accessory Number
21	Sample tray	20 pcs	20 pcs	AX-MXA-31
22	Tray holder	○	○	
23	Sample tray handle	2 pcs	2 pcs	AX-MXA-35
24	Breeze break	○	○	
25	Body cover	○	○	AX-MXA-39
26	Quick Start Guide	○	○	
27	Warranty card	○	○	
28	Display protective cover	○	○	AX-MXAT-38
29	Power cable	○	○	
30	USB cable 2 m ^{*2} (Type-A - Type-C)	○	○	AX-KO7919-200
31	Test sample ^{*1}	○	○	AX-MX-33
32	Fiberglass sheet	○	○	AX-MXA-32-2
33	Disposable aluminum tray	100 pcs	100 pcs	AX-MXA-30
34	Spoon	○	○	AX-MX-37
35	Tweezers	○	○	AX-MX-36
36	Cleaning brush (large/small)	○	○	AX-CLEANING-SET

○: Included -: Not included

^{*1} Sodium tartrate dihydrate 30 g

^{*2} UL certified products do not include a USB cable.

3. Screens and Controls


3.1. Standby screen

- At the standby screen, press the **ON:OFF** key, press the **HOME** key, or touch the screen to go to the home screen.


1	ON: OFF key	If at the standby screen, goes to the home screen. If at other than the standby screen, goes to the standby screen.
2	HOME key	Goes to the home screen.

3.2. Home screen

1	Program measurement button	If at the standby screen, goes to the home screen. If not at the standby screen, goes to the standby screen.
2	Log-in "User name" display	Displays the username of the logged-in user.
3	Simple measurement button	Starts the measurement by setting only the tray top temperature and the measurement accuracy.
4	Shortcut buttons	Starts the measurement by the program registered in the "Shortcut".
5	Menu button	Displays the menu screen.
6	Tray top temperature display	Displays the tray top temperature when it is high.
7	Communication device display	Displays the connected communication device.
8	Date & time display	Displays the current date and time.
9	Weighing button	Displays the weighing screen.
10	Log-in button	Displays the "Log-in" screen.

3.3. Measurement screen






1	Message display	Displays the measurement procedure.
2	<u>Program Info</u> button	Displays the program information screen.
3	<u>Cancel</u> button	Ends the measurement and returns to the home screen.
4	<u>Data output</u> button	Outputs (prints out) the weighing value and moisture content.
5	Main view	Displays the weighing value and moisture content.
6	Unit display	Displays the unit for the measurement value and moisture content.
7	Diagram display	Displays the diagrams and graphs showing the measurement procedure.
8	<u>Zero</u> button	Zeros the mass display.

4. Measurement Preparation





Step	Description
1.	Select the installation location referring to “1.1. Installation precautions”.
2.	(Adjust Level) Turn the left and right foot adjusters to center the bubble of the level within the black circle. ○ Black circle ● Bubbles Bubble spirit level level Tilted Up Foot adjuster Down Down Foot adjuster Up When the bubble is off to the left: Turn the right foot adjuster clockwise. When the bubble is off to the right: Turn the left foot adjuster clockwise. When the bubble is off to the rear position: Turn both foot adjusters clockwise at the same time. When the bubble is off to the forward position: Turn both foot adjusters counterclockwise at the same time.




Step	Description	Parts Diagram
3.	Align the slotted hole of the breeze brake with the protrusion on the main body and place it.	Breeze brake Slotted holes in the breeze brake Protrusion on the main body




Step	Description	Parts Diagram
4.	Place the tray holder so that the $\triangle$ mark on the tray holder aligns with the $\triangle$ mark on the main body.	 $\triangle$ symbol on the tray holder Tray holder $\triangle$ symbol on the protrusion of the body
5.	Place the sample tray on the sample tray holder, and then align the sample tray holder with the notch in the breeze brake. Caution ● If using a disposable aluminum tray, make sure to place it on the sample tray.	 Sample tray handle Sample tray Notch  Place on the notch
6.	Connect the power cable and ground the moisture analyzer. Note ● It is recommended to apply power for more than 30 minutes before measuring if the moisture content is below 1%.	 Power cable Grounding Methods When grounding with an outlet  Grounding terminal Power cable Outlet with grounding terminal

5. Measurement Procedure

5.1. Measurement procedure for moisture content (Simple measurement: When automatic measurement start is on)


Step	Description	Procedure
1.	Turn on the main unit, and when the standby screen is displayed, press the ON:OFF key to display the home screen.	
2.	Touch the Simple measurement button to display the simple measurement screen. Set the set temperature and measurement accuracy, and then touch the Next button.	 
3.	Place the sample tray on the sample tray handle, and then place it on the tray holder. Caution The sample tray handle should be placed in the notch part of the breeze break.	
4.	Close the heater cover. When closed, the mass display is automatically reset to zero. (Avoid any disturbances such as vibration during measurement.) If the mass display deviates from zero, press the Zero button again to display zero.	 

Step	Description	Procedure
5.	Open the heater cover and place an appropriate amount of sample using the level meter as a guide. Caution • A minimum of 0.1 g of sample is required. • Make the sample as flat as possible.	
Remarks	Press the Data output button to output (print out) the current mass display.	
6.	Wait until the mass indication stabilizes, and then close the heater cover. When closed, it will automatically start heating. As the moisture in the sample evaporates, the measurement value changes.	

Step	Description	Work
7.	When the end conditions (moisture content change is below the specified level, or heating for the specified time) are met, a buzzer will sound, and the measurement will end.	
Remarks	At the results display, the following operations can be done: Switch display button Switch to the graph view. Data output button Output (print out) the data. File output button Output PDF report to USB memory Next meas. button Start the measurement with the same program. End button Go back to the home screen.	
9.	Open the heater cover and take out the sample by holding the sample tray. Note The sample tray can be washed and reused.	




6. Maintenance

- Before performing maintenance, the power plug must be unplugged from the outlet.
- Make sure that each part of the moisture analyzer has cooled down sufficiently before performing maintenance.
- The sample tray, tray holder, and breeze break can be removed.
- Clean off dirt with a cloth soaked in water or water with a small amount of neutral detergent and wrung out tightly.
- Do not use organic solvents or chemical cloths.
- To transport, please use the special packaging box.

6.1. Cleaning the heater

- If the glass housing becomes dirty, it may not heat properly, so please keep it clean. The glass housing can be easily removed by unscrewing the two screws.
- Clean any fingerprints or dirt from the halogen lamp, as they can shorten its lifespan.
- Do not touch the reflective metal plate on the back side of the halogen lamp. This can cause a discrepancy between the set temperature of the sample tray and the actual heating temperature.
- Do not touch the temperature sensor adjacent to the halogen lamp. Doing so may cause a deviation between the set temperature of the sample tray and the actual heating temperature.


6.2. Troubleshooting

If the measurement results are suspected to be incorrect


- Please refer to the separate instruction manual and try self-inspection.
- When displaying mass, place and remove weights to check the repeatability of the measured values. Taller weights may come into contact with the heater, so please use shorter weights whenever possible. If you must use a 50 g weight, measure with the heater cover open. In this case, avoid disturbances such as wind. The maximum height of weights that can be placed on the sample tray (the height from the sample tray to the glass housing) is approximately 26 mm.
- During the function test, make sure the moisture content of the included test sample is measured correctly. For detailed instructions on the function test, refer to the separate 'Instruction Manual'.
- Check to see if any wind or vibration (due to air conditioning) around the moisture analyzer. Place the product on a stable table and keep it free from wind and vibration.
- Make sure the samples properly processed.
In particular, for cutting large particle samples and using glass fiber sheets, please refer to the separate instruction manual.
- Check that the product measurement method is performed correctly.
In particular, for preheating before measurement and for changing trays during continuous measurement, please refer to the separate instruction manual.

Weight height
Up to 26 mm

If the lamp does not light up when starting the measurement, or if it takes a long time to reach the set temperature

- If the heater cover is open, power is not supplied to the halogen lamp.
- The appropriate voltage is indicated on the connector of the halogen lamp. Make sure that the rated voltage matches the power supply voltage.
- Check the fuse condition.
After unplugging the power plug from the power outlet, check the fuse on the power inlet on the backside of the main unit. For the confirmation method, refer to the separate instruction manual.
- Make sure you are not measuring at a low temperature immediately after measuring at a high temperature.
If the temperature of the tray is hotter than the set temperature, the lamp will not light up.
- In cases other than the descriptions above, it is likely that the halogen lamp has burned out. Check the condition of the lamp filament before replacing the lamp. Please refer to the separate instruction manual for the replacement method.

7. Specifications

			MS-74AT	MX-53AT
Heating method			400 W halogen lamp	
Range of sample tray temperature settings			30°C to 200°C	
Measurable sample mass			0.1 g to 71 g	0.1 g to 51 g
Measure ment accuracy repeatabi lity, standard deviation	Moisture content ^{*1}	Sample weight: 5 g or greater	0.01%	0.02%
		Sample weight 1 g or greater	0.05%	0.1%
	Mass		0.0005 g	0.001 g
Minimum display	Moisture content		0.001%, 0.01%, 0.1%	0.01%, 0.1%
	Mass		0.0001 g	0.001 g
Display			5 inches WVGA TFT color LCD display	
Interface		RS-232C	D-Sub 9-pin (male) EIA RS-232C	
		USB (Device)	Type-C (female) USB 2.0 HID CDC	
		USB (Host)	Type-A (female) USB memory only	
		LAN	RJ45 TCP/IP port number: 10001 (fixed)	
		Bluetooth®	Bluetooth 5.0 (Bluetooth Low Energy)	
Operating temperature and humidity range			5°C to 40°C 85%RH or less (non-condensing)	
Range of use			Use indoors	
Altitude			Up to 2000 m	
Power			120 V AC Version 240 V AC Version Voltage fluctuations Power load	100 V to 120 V, 50/60 Hz, 3A ^{*3} 200 V to 240 V, 50/60 Hz, 1.5A –15%, +10% Approximately 500 W ^{*2}
Overvoltage category			II	
Pollution grade			2	
Maximum power consumption			500 W	
External dimensions			215 (W) ×380 (D) ×176 (H)	
Body weight			Approx. 6 kg (excluding accessories)	

^{*1} After preheating, measure the accessory test sample (approximately 5 g of sodium tartrate) at 160°C, using standard heating, standard mode, and measurement accuracy set to **MIDDLE**). After each measurement, leave the heater cover open and let it cool at room temperature for 10 minutes.

^{*2} Please confirm that this analyzer is compatible with your local voltage, receptacle type, and power cable.

^{*3} UL certified products are 120 V Version only.

8. Disposal

In accordance with the European Directive on Waste Electrical and Electronic Equipment (WEEE) 2012/19/EU, this product cannot be disposed of as general waste. Please dispose of it in accordance with the applicable laws of each country.

Dispose of this product at a recycling collection point for electrical and electronic equipment as specified by local regulations. If you have any questions, please contact the relevant administrative department. If you transfer this product to another person, please inform them of these disposal regulations.

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MEMO

[illegible]

THIS PAGE INTENTIONALLY LEFT BLANK.

A&D Company, Limited

3-23-14 Higashi-Ikebukuro, Toshima-ku, Tokyo 170-0013, JAPAN
Telephone: [81] (3) 5391-6132 Fax: [81] (3) 5391-1566

A&D ENGINEERING, INC.

Headquarters Office: 4622 Runway Boulevard, Ann Arbor, Michigan 48108, U.S.A.
Sales Office: 47747 Warm Springs Boulevard, Fremont, California 94539, U.S.A.
Tel: [1] (800) 726-3364 Weighing Support: [1] (888) 726-5931 Inspection Support: [1] (855) 332-8815

A&D INSTRUMENTS LIMITED

Unit 24/26 Blacklands Way, Abingdon Business Park, Abingdon, Oxfordshire OX14 1DY United Kingdom
Telephone: [44] (1235) 550420 Fax: [44] (1235) 550485

A&D AUSTRALASIA PTY LTD

32 Dew Street, Thebarton, South Australia 5031, AUSTRALIA
Telephone: [61] (8) 8301-8100 Fax: [61] (8) 8352-7409

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A&D Instruments India Private Limited

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A&D KOREA Limited