# MS-74A / MX-53A / MF-53A / ML-53A

(This document is included in the packing contents of the product.)

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MS-74A / MX-53A / MF-53A / ML-53A  
水分分析仪

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MS-74A / MX-53A / MF-53A / ML-53A  
水分分析儀

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Please see the separate  
instruction manual [here](#)

# SAFETY NOTES AND OTHER NOTES

## Word and symbol (for harm to persons / for damage to property):

|                                                                                                  |                                                                                                                                              |
|--------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------|
|  <b>WARNING</b> | If this notation is ignored and mishandled, it indicates that a person may be killed or seriously injured.                                   |
|  <b>CAUTION</b> | If this notation is ignored and mishandled, it indicates that a situation where a person is injured or property damage is expected to occur. |

## Graphic and symbol:

|                                                                                   |                                                                                                                                                                                                                                                                                                                            |
|-----------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  |  This symbol indicates mandatory actions that must be followed.<br>Specific mandatory actions are indicated by text or illustrations within or near  . |
|-----------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

## General message:

|                                                                                   |                                                                                        |
|-----------------------------------------------------------------------------------|----------------------------------------------------------------------------------------|
| <b>Caution</b>                                                                    | This is a description of important points for correct usage.                           |
| <b>Note</b>                                                                       | 'High possibility of inappropriate handling' or 'general advice on using the product'. |
|  | Indicates high temperature caution.                                                    |

## About This Manual

- (1) No part of this manual may be reprinted, copied, modified, or translated to another language without the prior written consent of A&D Company, Limited (A&D).
- (2) The contents of this manual are subject to change without notice.
- (3) Please contact A&D if you notice any uncertainty, errors, omissions, etc. in this manual.
- (4) A&D bears no liability for any loss or lost profits due to the operation of this product, and for direct, indirect, special, or consequential damages resulting from any defect in this product or this manual, even if advised of the possibility of such damage. Furthermore, A&D assumes no liability for claims of rights from third parties. Concurrently, A&D assumes no liability whatsoever for software or data losses.

© 2024 A&D Company, Limited

- ☐ Other product names and company names mentioned in this manual are trademarks or registered trademarks of their respective companies in Japan or other countries and regions.

## About this instruction manual

Thank you for purchasing an A&D MS-74A, MX-53A, MF-53A, or ML-53A moisture analyzer. This manual describes the safety precautions, installation instructions, basic functions and main operations of the moisture analyzer. For more detailed contents, please refer to the separate instruction manual from the A&D website.

## Safety Information

This product is a moisture analyzer that heats and dries samples using a halogen lamp based on the principle of thermogravimetric analysis, evaporating moisture to measure moisture content, solid content, etc. from the mass change. Do not use it for any other purposes.

This manual contains basic information for handling this product. Please read and understand this manual thoroughly before use.

This product is designed to be used by skilled or experienced operators.

Modifying or disassembling this product, or using it in a manner not specified in this manual or the separate instruction manual, may compromise safety. A&D assumes no responsibility for any issues arising from such actions.

# 1. Handling precautions


## 1.1. Installation precautions

### ⚠ WARNING

#### --Fatal accidents caused by fires and explosions--

- **Measure in an appropriate environment.**
  - ❑ Do not use the product in a hazardous environment such as in a flammable gas atmosphere.
  - ❑ Strictly adhere the following installation environment conditions:  
Temperature: 5 to 40°C, Humidity: 85%RH or less (without condensation)
  - ❑ Do not use the product in a poorly ventilated installation environment. Using the product in enclosed spaces where heat dissipation is obstructed can lead to unexpected temperature increases or poisoning symptoms.
- **⚠ Do not place flammable items around the product.**
  - ❑ The parts of the moisture analyzer become hot during and immediately after measurement. There is a risk of ignition.
  - ❑ Never place anything on top of the heater cover.
  - ❑ Do not operate the device while cables or other objects are hanging on it.


Heater cover

### ⚠ WARNING

#### --Risk of death or malfunction due to electric shock--

- **Check the power supply voltage.**  
The appropriate power supply voltage is indicated on the back of the moisture analyzer.  
Make sure that it matches the power supply voltage you are using.
- **Use the standard power cable provided when connecting the moisture analyzer to the power supply.**
- **Ground the moisture analyzer.**
- **Be careful of water exposure. The moisture analyzer is not waterproof. If water gets inside the product, it may cause electric shock or malfunction.**

## Caution

#### --Precautions for accurate measurement--

- **The mass sensor uses a high-precision electromagnetic balance system. Please pay attention to the following points.**
  - ❑ Install on a sturdy table.
  - ❑ Avoid vibrations.
  - ❑ Install in a location that is not exposed to wind from air conditioners or similar sources.


## 1.2. Precautions for use

### **⚠ WARNING**

#### **-- Risk of death due to fire, explosion, or other toxic hazards --**

- **Do not measure hazardous specimens.**
  - ❑ Do not heat samples that may explode or ignite, or heat samples that generate harmful substances, as doing so is extremely dangerous. Samples with unknown properties are likewise dangerous.
  - ❑ If the sample surface dries first and the internal pressure increases, it may burst. This is dangerous, so do not perform such measurement.
  - ❑ If the sample ignites, immediately unplug the product and take appropriate actions.  
Note that the case of the moisture analyzer is made of flame-retardant material (UL94 V-0).
  - ❑ Do not place flammable items nearby.
  - ❑ The parts of the moisture analyzer become hot during and immediately after measurement.  
Do not place flammable items nearby as doing so may cause a fire.
  - ❑ Never place anything on top of the heater cover.

Heater cover

### **⚠ WARNING**

#### **--Blindness and other accidents--**


- **Do not look directly at the halogen lamp during heating. Doing so may cause eye pain and vision problems.**
- **Wear appropriate protective equipment such as protective clothing, safety glasses, and protective gloves while measuring.**
- **Do not modify or disassemble the moisture analyzer. Doing so can cause malfunction, electric shock, or fire. If you suspect a malfunction, please consult the store where you purchased the product.**
- **Do not expose the moisture analyzer, power cable, or included accessories to extreme temperatures, strong chemical vapors, moisture, shock, vibration, or strong electromagnetic fields. Follow the operating conditions specified in “6. Specifications”**

## ⚠ Caution


### --⚠ Burns caused by hot surfaces--

- Follow the correct operating methods and procedures.


- ❑ Position the breeze break, tray holder, sample tray, and sample tray handle properly.
- ❑ Be sure to open and close the heater cover by holding the heater cover handle.

- ❑ During the measurement, the shaded area in the figure becomes hot. Only touch the heater cover handle and the sample tray handle during operation.

- ❑ Do not touch the inside of the heater cover or hot parts such as the sample tray directly with your hands immediately after heating. Doing so may result in burns. During and immediately after measurement, these parts are hot. When handling the moisture analyzer, use the designated handles and provided tools, and avoid touching the glass or metal parts near the halogen lamp (heater) in the glass housing.
- ❑ In particular, the glass housing part becomes hot and may cause burns if touched.

- ❑ Immediately after measurement, the sample tray and the part around the sample tray handle are hot. Cool it in an appropriate location. When handling the sample tray, use tweezers or similar tools.

**⚠ CAUTION****--Injuries caused by improper handling--**

- Dispose of used halogen lamp unaltered.  
Breaking a halogen lamp can cause glass fragments to scatter and lead to injury.
- Be careful not to drop, bump, or scratch glass components, including halogen lamp.  
They may break, leading to injury.
- The edges of the sample tray are sharp. Be careful not to cut yourself.
- Use the handle when opening and closing the heater cover. Be careful not to pinch your fingers.
- When moving the moisture analyzer, ensure that it is cool and lift it without tilting.  
Moving it while tilted may cause the heater cover to open unexpectedly, leading to injury.

**⚠ CAUTION****--Equipment failure due to improper handling--**

- It is recommended to replace halogen lamp after their rated lifespan (approximately 5000 hours) as failure to do so may cause damage.
- Be careful not drop, bump, or scratch glass parts, including halogen lamp. They may break, leading to injury.
- Make sure that dust or water does not enter the inside of the moisture analyzer.
- Use only A&D accessories and attachments.
- If the moisture analyzer behaves unexpectedly, unplug the power and plug it in again.  
If it does not work properly, please request a repair.

**Caution**

- **Precautions for high temperature heating**
  - ☐ Do not continuously heat at a temperature setting of 200°C for more than 30 minutes. The safety mechanism may activate and cut off the power supply to the halogen lamp.
  - ☐ For safety reasons, if heating at 160°C or higher continues for 1 hour, the maximum temperature will automatically be limited to 160°C.
  - ☐ During the measurement, the STOP key is always active. If you notice any abnormalities or potential danger, immediately press the STOP key to stop the measurement.

## 1.3. Precautions after use, precautions for storage, etc.

### CAUTION

#### -- Burns caused from hot surfaces--

- Make sure that each part of the moisture analyzer has cooled down sufficiently before performing any operations.
  - Especially when replacing the lamp, make sure the temperature around the glass housing has decreased before proceeding. For detailed lamp replacement instructions, refer to the separate 'Instruction Manual'.

#### --Equipment failure--

- Make sure that dust or water does not enter the inside of the moisture analyzer.
- When cleaning, use a cloth dampened with diluted neutral detergent. Do not use organic solvents or chemical wipes as they may cause damage.
- Always unplug the power cable from the outlet before performing maintenance.
- Use the dedicated packing box for transport.


### Caution

- Clean the glass housing if it becomes dirty, as it may not heat properly if dirty. Refer to “[5.1. Cleaning the heater](#)” for handling instructions.
- Clean the halogen lamp if it has fingerprints or other dirt, as those can shorten its lifespan. Refer to the separate instruction manual for handling instructions.
- Do not apply shock or drop the moisture analyzer as doing so may damage the halogen lamp or mass sensor.
- If the product will not be used for a long period of time, unplug the power cable from the outlet.

## 2. Packing contents and part names

### 2.1. Part names


- Make sure all the accessories are included.
- Keep the packaging material for use during transport or repairs.

| No. | Name                |
|-----|---------------------|
| 1   | Heater cover        |
| 2   | Heater cover handle |
| 3   | Display             |
| 4   | Leveler             |
| 5   | Keys                |
| 6   | Sample tray handle  |
| 7   | Foot adjuster       |
| 8   | Fuse (T6.3A 250V)   |
| 9   | Power inlet         |

| No. | Name                     |
|-----|--------------------------|
| 10  | RS-232C interface        |
| 11  | USB interface            |
| 12  | Disposable aluminum tray |
| 13  | Sample tray              |
| 14  | Tray holder              |
| 15  | Breeze break             |
| 16  | Halogen lamp             |
| 17  | Glass housing            |
| 18  | Serial number            |

## 2.2. Accessories


| No. | Name                                            | MS-74A  | MX-53A  | MF-53A  | ML-53A  | Accessory Number |
|-----|-------------------------------------------------|---------|---------|---------|---------|------------------|
| 18  | Sample tray                                     | 20 pcs  | 20 pcs  | 10 pcs  | 10 pcs  | AX-MXA-31        |
| 19  | Tray holder                                     | ○       | ○       | ○       | ○       |                  |
| 20  | Sample tray handle                              | 2 pcs   | 2 pcs   | 1 pc    | 1 pc    | AX-MXA-35        |
| 21  | Breeze break                                    | ○       | ○       | ○       | ○       |                  |
| 22  | Body cover                                      | ○       | ○       | —       | —       | AX-MXA-39        |
| 23  | Quick Start Guide                               | ○       | ○       | ○       | ○       |                  |
| 24  | Warranty card                                   | ○       | ○       | ○       | ○       |                  |
| 25  | Display protective cover                        | ○       | ○       | ○       | ○       | AX-MXA-38        |
| 26  | Power cable                                     | ○       | ○       | ○       | ○       |                  |
| 27  | USB cable 2m <sup>*2</sup><br>(Type-A — Type-C) | ○       | ○       | —       | —       | AX-KO7919-200    |
| 28  | Test sample <sup>*1</sup>                       | ○       | ○       | —       | —       | AX-MX-33         |
| 29  | Fiberglass sheet                                | ○       | ○       | —       | —       | AX-MXA-32-2      |
| 30  | Disposable aluminum tray                        | 100 pcs | 100 pcs | 100 pcs | 100 pcs | AX-MXA-30        |
| 31  | Spoon                                           | ○       | ○       | —       | —       | AX-MX-37         |
| 32  | Tweezers                                        | ○       | ○       | —       | —       | AX-MX-36         |
| 33  | Cleaning brush<br>(large/small)                 | ○       | ○       | —       | —       | AX-CLEANING-SET  |

○: Included      -: Not included

<sup>\*1</sup> Sodium tartrate dihydrate 30 g

<sup>\*2</sup> UL certified products do not include a USB cable.


## 2.3. Key Switches and Displays


| No. | Name                                                                   |
|-----|------------------------------------------------------------------------|
| 1   | Program Number<br>Data Number (when using data memory)                 |
| 2   | Measurement Mode                                                       |
| 3   | Heating Pattern                                                        |
| 4   | Sample Tray Temperature Display<br>(Set value / Actual measured value) |
| 5   | Measurement Time<br>(Set value / Elapsed time)                         |
| 6   | Measurement Accuracy                                                   |
| 7   | Target Sample Mass<br>Moisture Change Rate Display                     |

| No. | Name                                                                         |
|-----|------------------------------------------------------------------------------|
| 8   | Level Meter                                                                  |
| 9   | Operation Display                                                            |
| 10  | USB Connection Mark                                                          |
| 11  | Processing Mark                                                              |
| 12  | Stability Mark for Measured Value                                            |
| 13  | Main Displays<br>(Measured Value / Moisture Content /<br>Program Name, etc.) |
| 14  | Unit Display                                                                 |



### Example of Normal Mode Display

### Display Example of Operation Display

## Functions and Operations for Key Switches

| Key Switches                                                                      |                                                                                   | Functions and Operations                                                                                       |
|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------|
|  |  | During mass display, selects the program number.                                                               |
|  |  | Selects the measurement conditions.                                                                            |
|  |  | Changes the measurement conditions.                                                                            |
|  |  | Stores the measurement conditions or outputs the measured values.                                              |
|  |  | Starts the measurement.<br>However, the measurement will not start unless there is a sample of at least 0.1 g. |
|  |  | Stops the measurement.                                                                                         |
|  |  | Zeros the mass display.<br>Hold down to turn off the display.                                                  |






Refer the instruction manual for more detailed operation






### 3. Measurement Preparation






| Step | Description                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1.   | Select the installation location referring to “1.1. Installation precautions”.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| 2.   | <p>(Adjust Level) Turn the left and right foot adjusters to center the bubble of the level within the black circle.</p> <div style="display: flex; align-items: center; justify-content: space-around;"> <div style="text-align: center;"> <p>Bubble spirit level</p> <p>level      Tilted</p> </div> <div style="text-align: center;"> <p>Up<br/>Foot<br/>adjuster<br/>Down</p> </div> <div style="text-align: center;"> <p>Down<br/>Foot<br/>adjuster<br/>Up</p> </div> </div> <div style="display: flex; justify-content: space-between; margin-top: 10px;"> <div style="width: 48%;"> <p><b>When the bubble is off to the left:</b><br/>Turn the right foot adjuster clockwise.</p> </div> <div style="width: 48%;"> <p><b>When the bubble is off to the right:</b><br/>Turn the left foot adjuster clockwise.</p> </div> </div> <div style="display: flex; justify-content: space-between; margin-top: 10px;"> <div style="width: 48%;"> <p><b>When the bubble is off to the rear position:</b><br/>Turn both foot adjusters clockwise at the same time.</p> </div> <div style="width: 48%;"> <p><b>When the bubble is off to the forward position:</b><br/>Turn both foot adjusters counterclockwise at the same time.</p> </div> </div> |





| Step | Description                                                                                   | Parts Diagram                                                                                   |
|------|-----------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------|
| 3.   | Align the slotted hole of the breeze brake with the protrusion on the main body and place it. | <p>Breeze brake</p> <p>Slotted holes in the breeze brake</p> <p>Protrusion on the main body</p> |

| Step | Description                                                                                                                                                                                                                                                                                | Parts Diagram                                                                                                                                                                                                                                                                                                                                                                                   |
|------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 4.   | Place the tray holder so that the $\triangle$ mark on the tray holder aligns with the $\triangle$ mark on the main body.                                                                                                                                                                   |  <p><math>\triangle</math> symbol on the tray holder</p> <p>Tray holder</p> <p><math>\triangle</math> symbol on the protrusion of the body</p>                                                                                                                                                                |
| 5.   | <p>Place the sample tray on the sample tray holder, and then align the sample tray holder with the notch in the breeze brake.</p> <p><b>Caution</b></p> <ul style="list-style-type: none"> <li>● If using a disposable aluminum tray, make sure to place it on the sample tray.</li> </ul> |  <p>Sample tray handle</p> <p>Sample tray</p> <p>Notch</p>  <p>Place on the notch</p>                                                                                                                                      |
| 6.   | <p>Connect the power cable and ground the moisture analyzer.</p> <p><b>Note</b></p> <ul style="list-style-type: none"> <li>● It is recommended to apply power for more than 30 minutes before measuring if the moisture content is below 1%.</li> </ul>                                    |  <p>Power Plug</p> <div data-bbox="901 1680 1380 2049"> <p><b>Grounding Methods</b></p> <p>When grounding with an outlet</p>  <p>Grounding terminal</p> <p>Power cable</p> <p>Outlet with grounding terminal</p> </div> |

## 4. Measurement Procedure


| Step | Description                                                                                                                                                                                                                                                                                               | Procedure                                                                            |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|
| 1.   | Turn on the main unit and show the mass display.                                                                                                                                                                                                                                                          |    |
| 2.   | Place the sample tray on the sample tray handle, and then place that on the tray holder.<br><br><b>Caution</b> <ul style="list-style-type: none"> <li>The sample tray handle should be placed in the notch part of the breeze break.</li> </ul>                                                           |    |
| 3.   | Close the heater cover.                                                                                                                                                                                                                                                                                   |  |
| 4.   | When the 'Stable Mark' for the mass value lights up and the displayed value stabilizes, press the <b>RESET</b> key to zero the mass display. (Avoid disturbances such as vibration during measurement.)<br>If the mass display deviates from zero, press the <b>RESET</b> key again to return it to zero. |  |
| 5.   | Open the heater cover and use the level meter as a guide to place an appropriate amount of sample. Close the heater cover.<br><br><b>Caution</b> <ul style="list-style-type: none"> <li>A minimum of 0.1 g of sample is required.</li> <li>Make the sample as flat as possible.</li> </ul>                |  |

| Step    | Description                                                                                                                                                                                               | Procedure                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|---------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Remarks | While pressing the  key in the mass display, you can check the unit and minimum display of the moisture content results. |  <p>While pressing  </p>                                                                                                                                                                                           |
| Remarks | When you press the <b>ENTER</b> key, the current measuring value is output (printed out).                                                                                                                 |  <p> <b>PRINT / ENTER</b></p>                                                                                                                                                                                                                                                                                                                                                            |
| 6.      | <p>Close the heater cover, wait until the measuring value stabilizes, and press the <b>START</b> key.</p> <p>As the moisture in the sample evaporates, the measurement value changes.</p>                 |  <p></p>   <p>Measurement time</p> <p>Change in moisture content</p> <p>Sample tray temperature</p>  <p>Measured value</p> |

| Step        | Description                                                                                                                                                                                                                                                                                                                                                                                                                                                                         | Procedure                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |         |  |            |                              |            |                              |             |               |             |       |   |       |
|-------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------|--|------------|------------------------------|------------|------------------------------|-------------|---------------|-------------|-------|---|-------|
| Remarks     | <p>If you press <b>SELECT</b> key during the measurement, you can temporarily select a different measurement unit.</p> <h3>Caution</h3> <ul style="list-style-type: none"><li>Changing the measurement unit will also change the data output. Do not switch units while outputting data during measurement.</li></ul>                                                                                                                                                               | <div></div> <div><table><tr><th>Display</th><th></th></tr><tr><td>% MOIST /W</td><td>Moisture content (Wet basis)</td></tr><tr><td>% MOIST /D</td><td>Moisture content (Dry basis)</td></tr><tr><td>% RATIO D/W</td><td>Solid content</td></tr><tr><td>% RATIO W/D</td><td>Ratio</td></tr><tr><td>g</td><td>Grams</td></tr></table></div> | Display |  | % MOIST /W | Moisture content (Wet basis) | % MOIST /D | Moisture content (Dry basis) | % RATIO D/W | Solid content | % RATIO W/D | Ratio | g | Grams |
| Display     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |         |  |            |                              |            |                              |             |               |             |       |   |       |
| % MOIST /W  | Moisture content (Wet basis)                                                                                                                                                                                                                                                                                                                                                                                                                                                        |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |         |  |            |                              |            |                              |             |               |             |       |   |       |
| % MOIST /D  | Moisture content (Dry basis)                                                                                                                                                                                                                                                                                                                                                                                                                                                        |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |         |  |            |                              |            |                              |             |               |             |       |   |       |
| % RATIO D/W | Solid content                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |         |  |            |                              |            |                              |             |               |             |       |   |       |
| % RATIO W/D | Ratio                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |         |  |            |                              |            |                              |             |               |             |       |   |       |
| g           | Grams                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |         |  |            |                              |            |                              |             |               |             |       |   |       |
| 7.          | When the end conditions (moisture content change below a specified level, or heating for a specified amount of time) are met, a buzzer will sound, and the measurement will be end.                                                                                                                                                                                                                                                                                                 | <div></div>                                                                                                                                                                                                                                                                                                                                                                                                                |         |  |            |                              |            |                              |             |               |             |       |   |       |
| 8.          | <p>At the results display, the following operations can be done:</p> <p>The selected 'text' will flash:</p> <p><b>ENTER</b> Key<br/>Output (print out) data.</p> <p> ,  Key<br/>Switch the display unit of results between % or g.</p> <p><b>SELECT</b> Key, <b>RESET</b> Key<br/>g display is activated.</p> |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |         |  |            |                              |            |                              |             |               |             |       |   |       |
| 9.          | <p>Open the heater cover. Use the sample tray handle to lift the sample plate and remove the sample.</p> <h3>Note</h3> <ul style="list-style-type: none"><li>The sample tray can be washed and reused.</li></ul>                                                                                                                                                                                                                                                                    | <div><p>Removing the sample</p></div>                                                                                                                                                                                                                                                                                                                                                                                     |         |  |            |                              |            |                              |             |               |             |       |   |       |




## 5. Maintenance

- Before performing maintenance, the power plug must be unplugged from the outlet.
- Make sure that each part of the moisture analyzer has cooled down sufficiently before performing maintenance.
- The sample tray, tray holder, and breeze break can be removed.
- Clean off dirt with a cloth soaked in water or water with a small amount of neutral detergent and wrung out tightly.
- Do not use organic solvents or chemical cloths.
- To transport, please use the special packaging box.

### 5.1. Cleaning the heater

- If the glass housing becomes dirty, it may not heat properly, so please keep it clean. The glass housing can be easily removed by unscrewing the two screws.
- Clean any fingerprints or dirt from the halogen lamp, as they can shorten its lifespan.
- Do not touch the reflective metal plate on the back side of the halogen lamp. This can cause a discrepancy between the set temperature of the sample tray and the actual heating temperature.
- Do not touch the temperature sensor adjacent to the halogen lamp. Doing so may cause a deviation between the set temperature of the sample tray and the actual heating temperature.


## 5.2. Troubleshooting

### If the measurement results are suspected to be incorrect


- Please refer to the separate instruction manual and try self-inspection.
- When displaying mass, place and remove weights to check the repeatability of the measured values. Taller weights may come into contact with the heater, so please use shorter weights whenever possible. If you must use a 50 g weight, measure with the heater cover open. In this case, avoid disturbances such as wind. The maximum height of weights that can be placed on the sample tray (the height from the sample tray to the glass housing) is approximately 26 mm.
- During the function test, make sure the moisture content of the included test sample is measured correctly. For detailed instructions on the function test, refer to the separate 'Instruction Manual'.
- Check to see if any wind or vibration (due to air conditioning) around the moisture analyzer. Place the product on a stable table and keep it free from wind and vibration.
- Make sure the samples properly processed.  
In particular, for cutting large particle samples and using glass fiber sheets, please refer to the separate instruction manual.
- Check that the product measurement method is performed correctly.  
In particular, for preheating before measurement and for changing trays during continuous measurement, please refer to the separate instruction manual.

Weight height  
Up to 26 mm











### If the lamp does not light up even after the measurement starts, or if the set temperature is reached late

- It usually takes about 6 seconds for the light to turn on after pressing the **START** key.
- When the heater cover is open, power is not supplied to the halogen lamp.
- The appropriate voltage is indicated on the connector of the halogen lamp. Make sure that the rated voltage matches the power supply voltage.
- Check the fuse condition.  
After unplugging the power plug from the power outlet, check the fuse on the power inlet on the backside of the main unit. For the confirmation method, refer to the separate instruction manual.
- Make sure you are not measuring at a low temperature immediately after measuring at a high temperature.  
If the temperature of the tray is hotter than the set temperature, the lamp will not light up.
- In cases other than the descriptions above, it is likely that the halogen lamp has burned out. Check the condition of the lamp filament before replacing the lamp. Please refer to the separate instruction manual for the replacement method.

### 5.3. Error Display

| Error Display                                                                       | Contents and Countermeasures                                                                                                                                                                                                                                                                                                                                                                                                                                     |
|-------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | <b>Overload</b><br>The sample mass has exceeded the allowable range. Reduce the sample.<br>If this error occurs with only the sample tray, request a repair.                                                                                                                                                                                                                                                                                                     |
|    | <b>Underload</b><br>The sample is too light. (The output from the mass sensor is too small)<br>After setting the tray and sample tray correctly, press the <b>RESET</b> key.<br>Adjust the sensitivity of the mass sensor.<br>If it still doesn't work, request a repair.                                                                                                                                                                                        |
|    | <b>Sample Quantity Error at Heating Start</b><br>Adjust the sample quantity (increase or decrease) and proceed with heating.<br>You can set the sample quantity as desired. Refer to the separate instruction manual.                                                                                                                                                                                                                                            |
|    | <b>Self-Inspection Errors</b><br>An error occurred during self-inspection. Request a repair.                                                                                                                                                                                                                                                                                                                                                                     |
|    | <b>Heater cover error</b><br>If this error persists, repairs are needed.                                                                                                                                                                                                                                                                                                                                                                                         |
|   | <b>Sensitivity adjustment weight defect (Positive)</b><br>The weight for adjusting the sensitivity of the mass sensor is too heavy. Check the area around the trays. Make sure that the sensitivity adjustment weight of the mass sensor is not in contact with the glass housing of the heater cover, and that the mass of the sensitivity adjustment weight of the mass sensor is appropriate. Press any key or wait 15 seconds to return to the mass display. |
|  | <b>Sensitivity adjustment weight defect (Negative)</b><br>The weight for adjusting the sensitivity of the mass sensor is too light. Check the area around the trays. Make sure that the mass of the weight for adjusting the sensitivity of the mass sensor is appropriate. Press any key or wait 15 seconds to return to the mass display.                                                                                                                      |
|  | <b>Data entry time up when adjusting the heating temperature</b><br>This indicates that no temperature input was received during the heating temperature adjustment for a certain period (5 minutes). Press any key to clear the error. If you need to adjust the heating temperature, start over again from the beginning.                                                                                                                                      |
|  | <b>Full memory</b><br>The number of stored measurement results in the data memory function has reached its limit. To store new data, delete the existing memory. Refer to the separate instruction manual.                                                                                                                                                                                                                                                       |
|  | <b>Built-in Clock Battery Error</b><br>Press any key and enter the date and time. Refer to the separate instruction manual.<br>If the error occurs again, request a repair.                                                                                                                                                                                                                                                                                      |
|  | <b>Built-in Clock Malfunction</b><br>Please request a repair.                                                                                                                                                                                                                                                                                                                                                                                                    |
|  | <b>Internal error</b><br>Turn off the power once, then turn it on again, and check.<br>Make sure that the power supply (voltage, frequency) is appropriate.<br>If the error recurs, request a repair.                                                                                                                                                                                                                                                            |

| Error Display                                         | Contents and Countermeasures                                                                                                                            |
|-------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------|
| <div>Error3</div> <div>Error8</div> <div>Error9</div> | <b>IC errors</b><br>Request a repair.                                                                                                                   |
| <div>Ht Err</div>                                     | <b>Temperature control error</b><br>Turn off the power, then leave it for 30 minutes or longer, and check again. If the error recurs, request a repair. |
| <div>LOWVOLT</div>                                    | <b>AC Power Voltage Error</b><br>Check the power voltage.<br>If using a power strip with other devices, and the power supply is fine, request a repair. |
| <div>FREQErr</div>                                    | <b>AC Power Frequency Error</b><br>Check if the power supply is appropriate.<br>If the power supply is fine, request a repair.                          |

## 6. Specifications

|                                                                             |                                   |                                        | MS-74A                                                                     | MX-53A                                                                                                                          | MF-53A                | ML-53A      |
|-----------------------------------------------------------------------------|-----------------------------------|----------------------------------------|----------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------|-----------------------|-------------|
| Heating method                                                              |                                   |                                        | 400 W halogen lamp                                                         |                                                                                                                                 |                       |             |
| Range of sample tray temperature settings                                   |                                   |                                        | 30°C to 200°C                                                              |                                                                                                                                 |                       |             |
| Measurable sample mass                                                      |                                   |                                        | 0.1 g to 71 g                                                              | 0.1 g to 51 g                                                                                                                   |                       |             |
| Measure<br>ment<br>accuracy,<br>repeatabi<br>lity,<br>standard<br>deviation | Moisture<br>content <sup>*1</sup> | Sample<br>weight:<br>5 g or<br>greater | 0.01%                                                                      | 0.02%                                                                                                                           | 0.05%                 | 0.1%        |
|                                                                             |                                   | Sample<br>weight<br>1 g or<br>greater  | 0.05%                                                                      | 0.1%                                                                                                                            | 0.2%                  | 0.5%        |
|                                                                             | Mass                              |                                        | 0.0005 g                                                                   | 0.001 g                                                                                                                         | 0.002 g               | 0.005 g     |
| Minimum<br>display                                                          | Moisture content                  |                                        | 0.001%,<br>0.01%,<br>0.1%                                                  | 0.01%,<br>0.1%                                                                                                                  | 0.05%,<br>0.1%,<br>1% | 0.1%,<br>1% |
|                                                                             | Mass                              |                                        | 0.0001 g                                                                   | 0.001 g                                                                                                                         | 0.002 g               | 0.005 g     |
| Interface                                                                   |                                   | RS-232C                                | D-Sub 9-pin (male) EIA RS-232C                                             |                                                                                                                                 |                       |             |
|                                                                             |                                   | USB                                    | Type-C (female) USB 2.0 HID CDC                                            |                                                                                                                                 |                       |             |
| Operating temperature and humidity range                                    |                                   |                                        | 5°C to 40°C 85%RH or less (non-condensing)                                 |                                                                                                                                 |                       |             |
| Range of use                                                                |                                   |                                        | Use indoors                                                                |                                                                                                                                 |                       |             |
| Altitude                                                                    |                                   |                                        | Up to 2000 m                                                               |                                                                                                                                 |                       |             |
| Power                                                                       |                                   |                                        | 120 V AC Version<br>240 V AC Version<br>Voltage fluctuations<br>Power load | 100 V to 120 V, 50/60 Hz, 3A <sup>*3</sup><br>200 V to 240 V, 50/60 Hz, 1.5A<br>−15%, +10%<br>Approximately 500 W <sup>*2</sup> |                       |             |
| Overvoltage category                                                        |                                   |                                        | II                                                                         |                                                                                                                                 |                       |             |
| Pollution grade                                                             |                                   |                                        | 2                                                                          |                                                                                                                                 |                       |             |
| Maximum power consumption                                                   |                                   |                                        | 500W                                                                       |                                                                                                                                 |                       |             |
| External dimensions                                                         |                                   |                                        | 215 (W) × 380 (D) × 176 (H)                                                |                                                                                                                                 |                       |             |
| Body weight                                                                 |                                   |                                        | Approx. 6 kg (excluding accessories)                                       |                                                                                                                                 |                       |             |

<sup>\*1</sup> After preheating, measure the accessory test sample (approximately 5 g of sodium tartrate) at 160°C, using standard heating, standard mode, and measurement accuracy set to (MID.). After each measurement, leave the heater cover open and let it cool at room temperature for 15 minutes.

<sup>\*2</sup> Please confirm that this analyzer is correct for your local voltage and receptacle type and the power cable.

<sup>\*3</sup> UL certified products are 120 V Version only.

## 7. Disposal

In accordance with the European Directive on Waste Electrical and Electronic Equipment (WEEE) 2012/19/EU, this product cannot be disposed of as general waste. Please dispose of it in accordance with the applicable laws of each country.

Dispose of this product at a recycling collection point for electrical and electronic equipment as specified by local regulations. If you have any questions, please contact the relevant administrative department. If you transfer this product to another person, please inform them of these disposal regulations.

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# MEMO

[illegible]

THIS PAGE INTENTIONALLY LEFT BLANK.

### **A&D Company, Limited**

3-23-14 Higashi-Ikebukuro, Toshima-ku, Tokyo 170-0013, JAPAN  
Telephone: [81] (3) 5391-6132 Fax: [81] (3) 5391-1566

### **A&D ENGINEERING, INC.**

Headquarters Office: 4622 Runway Boulevard, Ann Arbor, Michigan 48108, U.S.A.  
Sales Office: 47747 Warm Springs Boulevard, Fremont, California 94539, U.S.A.  
Tel: [1] (800) 726-3364 Weighing Support: [1] (888) 726-5931 Inspection Support: [1] (855) 332-8815

### **A&D INSTRUMENTS LIMITED**

Unit 24/26 Blacklands Way, Abingdon Business Park, Abingdon, Oxfordshire OX14 1DY United Kingdom  
Telephone: [44] (1235) 550420 Fax: [44] (1235) 550485

### **A&D AUSTRALASIA PTY LTD**

32 Dew Street, Thebarton, South Australia 5031, AUSTRALIA  
Telephone: [61] (8) 8301-8100 Fax: [61] (8) 8352-7409

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### **A&D Instruments India Private Limited**

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### **A&D KOREA Limited**