

Ambulatory Blood Pressure Monitor TM2440-EX, Instruction Manual, Corrections and Additions

Thank you for purchasing the Ambulatory Blood Pressure Monitor TM2440-EX. The Instruction Manual that accompanies the product has been corrected and had content added. Please read it in conjunction with the details shown below.

Instruction manual Page 22 6.1. Automatic Blood Pressure Measurement (A-BPM)

The following changes have been made.

Wrong	Correct
The initial pressurization value is set to 180 mmHg at the factory.	The initial pressurization value for automatic blood pressure measurement is AUTO. In the case of AUTO, the appropriate pressurization value is set automatically.

Instruction manual Page 37 8.2.3. Initial Pressurization Value

The following changes have been made.

Wrong	Correct
The initial pressurization value is set to 180 mmHg at the factory.	The initial pressurization value setting is AUTO (appropriate pressurization value set automatically) when the ⌚ mark is displayed and 180 mmHg when the ⌚ mark is hidden.

Instruction Manual Page 37 8.3. A-BPM Preset Programs

The following is added after "A-BPM can use only while the automated measurement can be performed. "

Addition		
<table border="1"> <thead> <tr> <th style="text-align: center;">Note</th> </tr> </thead> <tbody> <tr> <td>The Start time and Operation time are initially set calculating from the time at which the ⌚ mark was displayed. When the next automatic measurement is made, this needs to be set again.</td> </tr> </tbody> </table>	Note	The Start time and Operation time are initially set calculating from the time at which the ⌚ mark was displayed. When the next automatic measurement is made, this needs to be set again.
Note		
The Start time and Operation time are initially set calculating from the time at which the ⌚ mark was displayed. When the next automatic measurement is made, this needs to be set again.		

The following changes have been made.

Wrong	Correct		
<p>#3: Example for automated measurement. Start time : Stores a time. (0 to 23 hour) Operation time : Set to "OFF" Response : A-BPM starts the blood pressure measurement at the preset Start time and continues until A-BPM is suspended. If the EVENT switch is pressing and holding until the  mark is shown, A-BPM starts at the preset Start time.</p>	<p>#3: If the Start time for automatic measurement is specified and the Operation time is set to "OFF", long pressing the EVENT switch will start automatic blood pressure measurement from the Start time set and measurement will continue until the EVENT switch is long pressed. If the EVENT switch is long pressed again, automatic measurement will continue.</p> <table border="1" data-bbox="815 645 1441 880"> <tr> <td data-bbox="815 645 1441 696" style="text-align: center;">Note</td> </tr> <tr> <td data-bbox="815 696 1441 880">If the Operation time is set, automatic blood pressure measurement will be performed from when the EVENT switch is initially long pressed until the Operation time has passed even if the EVENT switch is long pressed midway through automatic measurement.</td> </tr> </table>	Note	If the Operation time is set, automatic blood pressure measurement will be performed from when the EVENT switch is initially long pressed until the Operation time has passed even if the EVENT switch is long pressed midway through automatic measurement.
Note			
If the Operation time is set, automatic blood pressure measurement will be performed from when the EVENT switch is initially long pressed until the Operation time has passed even if the EVENT switch is long pressed midway through automatic measurement.			
<p>#4: Example for automated measurement. Start time : Set to "OFF" Operation time : Stores time to be continued. (1 to 27 hours) Response : A-BPM starts blood pressure measurement and stops after the Operation time. If the  mark is hidden by pressing and holding the EVENT switch during the Operation time, A-BPM stops. If the  mark is shown by pressing and holding the EVENT switch again, A-BPM is performed for the Operation time.</p>	<p>#4: If the Start time for automatic measurement is set to "OFF" and the Operation time is set, long pressing the EVENT switch will implement automatic blood pressure measurement after initial automatic measurement until the set Operation time has passed. If the EVENT switch is long pressed midway through automatic blood pressure measurement, automatic measurement will stop. If the EVENT switch is long pressed again, automatic measurement will be performed until the set Operation time has passed.</p> <table border="1" data-bbox="815 1254 1441 1451"> <tr> <td data-bbox="815 1254 1441 1305" style="text-align: center;">Note</td> </tr> <tr> <td data-bbox="815 1305 1441 1451">If the Start time is set, long pressing the EVENT switch midway through automatic blood pressure measurement will stop automatic measurement and long pressing it again will resume measurement.</td> </tr> </table>	Note	If the Start time is set, long pressing the EVENT switch midway through automatic blood pressure measurement will stop automatic measurement and long pressing it again will resume measurement.
Note			
If the Start time is set, long pressing the EVENT switch midway through automatic blood pressure measurement will stop automatic measurement and long pressing it again will resume measurement.			

The following notes were added to "Example 2 Automatic measurement"

Addition
NB. Automatic blood pressure measurement will continue even if turned "OFF" and "ON" again during automatic measurement.

The following notes were added to "Example 3 Automatic measurement"

Addition
NB. Automatic blood pressure measurement will not be performed once the Operation time has been exceeded even if turned "OFF" and "ON" again during automatic measurement.