# Digital Blood Pressure Monitor

**UA-779**

***Instruction Manual***

***Manuel d'instructions***

***Manual de Instrucciones***

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## Preliminary remarks

The device conforms to the following requirements: European Directive 93/42 EEC for Medical Products Act; Medical Products Act; European Standards for Electrical Medical Equipment EN 60601-1 (General Safety Provisions), EN 60601-1-2 and EN 55011 (Electromagnetic Compatibility).

The above is evidenced by the CE mark of conformity accompanied by the reference number of a designated authority.


This device is designed for adults only.

### Environment for use




The device is for use in a hospital and / or patient's home.

## Precaution

1. Precision components were used in the construction of this device. Extremes in temperature, humidity, direct sunlight, shock or dust should be avoided.
2. Clean the device with a dry, soft cloth. Never use thinner, alcohol, benzene, or wet dusters.
3. Avoid tightly folding the cuff or storing the hose tightly twisted for long periods, as such treatment may shorten the life of the components.
4. The device, cuff and tubing are not water resistant. Prevent rain, sweat and water from wetting the device, cuff and tubing.

# Parts Identification

# Symbols

Symbols	Function / Meaning	Treatment
	Turns on or off the device	_____
	Direction guide to install battery	_____
	Direct current	_____
<b>SN</b>	Serial number	_____
	Date of manufacture	_____
	Device, cuff and tubing are designed to provide special protection against electric shocks.	_____
	Appears when measurement is in progress and flashes when the pulse is detected.	Measurement is in progress - remain as still as possible.
	Appears when the battery voltage is too low for the device to work properly.	Replace all batteries with new ones.
	Appears if the pressure value is unstable due to movement during the measurement.	Try the measurement again. Remain very still during the measurement.
	Appears if the cuff is not attached properly or the cuff is not connected to the device firmly.	<ul style="list-style-type: none"> <li>? Attach the cuff properly and try the measurement again.</li> <li>? Insert the air hose into the air socket firmly.</li> </ul>
	Appears if pulses are not detected correctly.	Attach the cuff properly and try the measurement again.


**M**

Previous measurement stored in  
MEMORY

# Preparation for Measurements

## 1. Installing / Changing the Batteries

Remove the battery cover and insert new batteries into the battery compartment as shown, taking care to observe the proper direction for each battery. Use R6P, LR6, AA, or similar type batteries.


If the symbol  appears on the display panel, this means that the batteries are depleted. Always change all batteries at the same time. If you do not intend to use the device for longer periods, remove the batteries. Only use high quality, leakproof batteries with the specified ratings. Never mix new and old batteries or batteries from different manufactures. Used batteries do not fall within the category of household waste, and must be disposed of according to the applicable regulations.

## 2. Connecting the Air Hose

Insert the air connector plug into the air socket firmly.

## 3. Attaching the Arm Cuff

Wrap the cuff around the upper arm about 2-3 cm above the elbow as shown. Place the cuff directly against the skin, as clothing may cause a faint heart beat, and result in an error. Constriction of the upper arm caused by rolling up a shirt sleeve may prevent accurate readings.

## 4. Posture While Measuring

Remain seated in a comfortable position, placing the center of the upper arm at the same height as your heart.

## 5. Measurement


During the measurement, the cuff will become very tight. Do not be alarmed as this is normal.

## 6. After the Measurement

After the measurement, remove the cuff and record your data.


# Measurements

1. Place the cuff on the arm that will be used for the measurement (preferably the left arm).
2. Press the START/STOP switch.
  - When the START/STOP switch is pressed, all of the display symbols will appear for about one second.


- When the "0" is displayed, the device is calibrating its zero point. If there is still some air in the arm cuff, it will be exhausted at this time.
- The micro pump will pressurize the cuff.  
**Note: The cuff will become very tight. Do not be alarmed, as this is normal.**
- If the START/STOP switch is pressed again during pressurization or during the measurement, the pump will stop and the air will be released from the cuff.

  3. When pressurization has completed, the automatic exhaust mechanism will gradually reduce the pressure in the cuff and the Measurement in Progress symbol will appear along with the current pressure reading. This symbol will flash in unison with the pulse of the person being monitored.

# Measurements

- Remain still while the device is in operation.
  - When a measurement is made with insufficient pressure, the device will pump again to a higher pressure.
  - If necessary, the starting pressure can be set manually. Press the START/STOP switch and hold it down until the display shows a pressure 30-40 mmHg above your normal systolic blood pressure, then release the switch. The device will start the measurement from that pressure.
4. When measurement has completed, the systolic pressure, diastolic pressure readings, and the pulse rate are displayed. The air is automatically released from the cuff.

5. Turning the Power OFF  
Press the START/STOP switch to turn the power off.
6. Subsequent Measurements  
If a subsequent measurement is required, turn off the power and turn it on again. When the "0" is displayed, the device is ready for measurement again.  
**Note: Please wait for about 10 minutes before repeating the measurement.**
7. Automatic Power OFF Function  
If this device is left on after a measurement, it will turn itself off after about 1 minute. It may be turned off at any time by pressing the START/STOP switch.


# Memory Operation

## Recalling Memory Data



Note)

This device stores up to seven of the latest measurements in memory automatically. The data is retained as long as the batteries are in the device. Once you remove the batteries, the data will be cleared.

1. With nothing appeared on the display, press the MEMORY button.
2. The newest data appears. First, the memory No., then its data, two times.
3. When the MEMORY button is pressed while data is displayed, the previous data appears. Each time the MEMORY button is pressed, the memorized data appears one by one in the order of decreasing memory number. After last data appears, the display goes blank.

With nothing displayed, press MEMORY button.

Each time MEMORY button is pressed, the previous data appears.

## Notes for Proper measurements


1. Sit in a comfortable position where you can place the arm to be used for the measurement on a table or other support that will place the center of the upper arm at about the same height as your heart.
2. Relax for about 5 or 10 minutes before starting the measurement.
3. A normal person's blood pressure varies constantly depending on what you are doing and what you have eaten. What you drink can have a very strong and rapid effect on your blood pressure.
4. As this device bases its measurements on the heart beat, if you have a very weak or irregular heart beat, the device may have difficulty determining your proper blood pressure.
5. Should the device detect a condition that is abnormal, it will stop the measurement and display an error. See the section on description of display symbols for details.
6. This blood pressure device is intended for use by adults only. Consult with your physician before using this device on a child. A child should not use this device unattended.

# How to Make Proper Measurements

## How is a Proper Measurement Made?

Blood pressure varies depending upon the conditions prevailing at the time of the measurement, consider the following so that the most accurate measurements may be made.

1. In preparation for blood pressure measurement, the subject should urinate and should remain relatively still for 5 to 10 minutes before a measurement.
2. Shirts or other garments which fit tightly on the upper arm should be removed before attaching the arm cuff.
3. Exercise, eating, drinking, and smoking before a measurement can affect the results.
4. Blood pressure varies constantly throughout the day. Measurement should be made regularly at the same time each day.
5. Do not be too impressed by the results of one measurement. Keep a record of blood pressure measurements. Many readings tell a story.
6. When making repeated measurements, the arm becomes heavy with blood, resulting in wide variations in measurement values. For this reason, repeated measurements should be made after a rest period of approximately 5 to 10 minutes or after raising the arm until the engorged condition is relieved.
7. In general, the blood pressure is low in summer (when it is hot) and high in winter (when it is cold).
8. Emotional stresses may cause an increase in blood pressure.

# About Blood Pressure

## What is Blood Pressure?

The blood pressure is at its highest levels in the large arteries near the heart and drops off towards the peripheral areas of the circulatory system.

The blood pressure varies with the beating of the heart. When the heart contracts, forcing the blood out, the pressure inside the arteries is said to be systolic. When the heart expands, the pressure of the blood inside the arteries is said to be diastolic.

---

## Hypertension

Hypertension, an abnormally high arterial blood pressure, which is most common among older adults, if left unattended, can cause many health problems including stroke and heart attack. It is wise to control the blood pressure to prevent it from becoming high, by reducing salt intake, and by controlling diet and activities. People who were born with high blood pressure can prevent the progress of heart disease by means of adequate control.

---

## Why is it a Good Thing to Measure Blood Pressure at Home

Having your blood pressure measured in a hospital or a clinic may cause you to become nervous, thus artificially raising the blood pressure reading. Also blood pressure varies with a variety of conditions. So an accurate judgment of the heart condition can not be made with a single measurement.

Having the capability to measure blood pressure at home provides a good way to determine your fundamental blood pressure. This is defined as the first blood pressure reading upon waking, before taking any food or drink.

It is possible for you to carry out your own blood pressure control at home. Take blood pressure readings with your own blood pressure device on a regular basis. Record the measurements. If the accumulated results seem to be above normal, it is recommended that a doctor review the record.

# Specifications

## Performance Specifications

Measurement method Oscillometric measurement

Measurement range Pressure : 20 - 280 mmHg  $\pm 4$  mmHg

Pulse : 40 - 200 beats / minute  $\pm 5$  %

Battery 4 x 1.5V (R6P, LR6, or AA)

Cuff Arm Circumference

Classification Type BF 

Device, cuff and tubing are designed to provide special protection against electric shocks.

Clinical test According to ANSI / AAMI SP-10 1987

	Standard deviation	Mean error
Systolic blood pressure	5.43 mmHg	2.24 mmHg
Diastolic blood pressure	5.18 mmHg	1.54 mmHg

## Environment Specifications

Operating environment +10 °C - +40 °C, less than 85 %RH

Storage environment -20 °C - +70 °C, less than 85 %RH

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**A&D Company, Limited**
**A&D INSTRUMENTS LTD.**
**A&D MERCURY PTY. LTD.**