Wrist Digital Blood Pressure Monitor

UB-401

Instruction Manual

Manuel d'instructions

Manual de Instrucciones

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Preliminary Remarks

The device conforms to European Directive 93/42 EEC for Medical Products. This is evidenced by the CE mark of conformity accompanied by the reference number of a designated authority.


The device is designed for all except neonates and infants.

Environment for use


The device is for use in a patient's home.

Precautions



1. Precision components were used in the construction of this device. Extremes in temperature, humidity, direct sunlight, shock or dust should be avoided.
2. Clean the device with a dry, soft cloth. Never use thinner, alcohol, benzine, or wet dusters.
3. Avoid tightly folding the cuff or storing the hose tightly twisted for long periods, as such treatment may shorten the life of the components.
4. The device and cuff are not water resistant. Prevent rain, sweat and water from wetting the device and cuff.
5. Measurements may be impaired if the device is used close to televisions, microwave ovens, cellular telephones, X-ray or other devices with strong electrical fields.
6. Used equipment, parts and batteries are not treated as ordinary household waste, and must be disposed of according to the applicable regulations.

Parts Identification

(An example data display)


Symbols

Symbols	Function / Meaning	Treatment
	Turns on and off the device	_____
	Direction guide to install battery	_____
	Direct current	_____
SN	Serial number	_____
2000	Date of manufacture	_____
	Type BF: Device, cuff and tubing are designed to provide special protection against electrical shocks.	_____
	Appears when measurement is in progress and flashes when the pulse is detected.	Measurement is in progress -- remain as still as possible.
	Previous measurement stored in MEMORY	_____
	Appears when the battery voltage is too low for the device to work properly.	Replace all batteries with new ones.
	Appears if the systolic and diastolic measurements are within 10 mmHg of each other.	Fasten the cuff correctly, and try the measurement again.
	Appears if the pressure value is unstable due to movement during the measurement.	Try the measurement again. Remain very still during the measurement.
	Appears if the pressure value did not increase during pump up.	Fasten the cuff correctly, and try the measurement again.
	Appears when the cuff is not fastened correctly.	Fasten the cuff correctly, and try the measurement again.
	PUL.DISPLAY ERROR	Fasten the cuff correctly, and try the measurement again.
SYS.	Systolic blood pressure in mmHg	_____
DIA.	Diastolic blood pressure in mmHg.	_____
PUL./min	Pulse beats per minute	_____
	The label of the medical device by the EC directive	_____
	The certification trade mark registered to the Australian Communication Authority by the Trademarks office.	_____

Using the Monitor

1. Installing the batteries

1. Slide the battery cover as shown to remove it from the device.
2. Insert two Micro batteries in the box, with positive (+) and negative (-) terminals matching those indicated on the cuff.
3. Slide the battery cover back on the device.

CAUTION

- Insert the batteries as shown on the cuff. If not, the device will not work.
- When  (LOW BATTERY mark) appears in the display, replace all batteries with new ones at the same time. Do not mix old and new batteries. It may cause the battery life to be shortened, or the device to malfunction.
-  (LOW BATTERY mark) does not appear when the batteries are drained.
- Battery life varies with the ambient temperature and may be shorter at low temperatures.
- Remove batteries if the device is not to be used for a long period of time. The batteries may leak and cause a malfunction.
- Use the specified batteries only. The batteries provided with the device are for testing the monitor performance and may have a shorter life.


Using the Monitor

2. Fastening the cuff

1. Wrap the cuff around your wrist about 1 cm above your hand as shown in the figure at the right.
2. Fasten the cuff tightly using the Velcro strip

(Note)

For proper measurements, fasten the cuff tightly and measure on a bare wrist.


3. How to take proper measurements

For best accuracy in blood pressure measurement:

- Sit comfortably at a table. Rest your arm on the table.
- Relax for about 5 to 10 minutes before measurement.
- Raise your hand so that the cuff is at the same level as your heart.
- Remain still and keep quiet during measurement.
- Do not measure right after physical exercise or a bath.

Rest twenty or thirty minutes before taking measurement.

- Measure your blood pressure at about the same time every day.

4. Measurement

During measurement, it is normal for the cuff to feel very tight. (Do not be alarmed.)

5. After measurement

After measurement, press the START button to turn off the power.
Remove the cuff and record your data.

(Note)

The device is provided with an automatic power shut-off function which turns off the power automatically one minute after measurement.

Allow at least ten minutes between measurements on the same person.


Measurements

Note)

The UB-401, once used, will provide an inflation appropriate to the user, based on the previous measurement.

If your systolic pressure is expected to exceed 150 mmHg, read "For those whose systolic pressure is over 150 mmHg" on page 7.

1. Sit comfortably with the cuff at the same level as your heart and relax.

2. Press the START button.

All display symbols appear briefly, then the display changes as indicated in the figure at the right as measurement begins.


The cuff starts to inflate. It is normal for the cuff to feel very tight.

Note)

If you wish to stop inflation at any time, press the START button again.

3. When inflation is complete, deflation starts automatically and the (heart mark) appears to indicate that measurement is in progress.

Once the pulse is detected, the mark flashes with each pulse beat.


Note)

If an appropriate pressure is not obtained, the device starts to inflate again automatically.

4. When measurement is complete, the systolic and diastolic pressure readings and pulse rate are displayed.

The cuff exhausts the remaining air and deflates completely.

5. Press the START button again to turn off the power.

Remove the cuff.

Note)

The device is provided with an automatic power shut-off function.

Allow at least ten minutes between measurements on the same person.


Measurements

For those whose systolic pressure is over 150 mmHg

1. Fasten the cuff around your wrist.
Raise your arm so that the cuff is at the same level as your heart.
2. Press and hold the START button until the memory data appears on the display.
3. Release the START button.
4. Press and hold the START button again until a number about 30 to 40 mmHg higher than your expected systolic pressure appears.
5. When the desired number is reached, release the START button. Measurement will begin.

Note)

The operation above is required only when you measure your blood pressure for the first time. Next time, the UB-401 will provide an appropriate inflation. So, measure your blood pressure as described on page 6.


Recalling Memory Data

Note)

This device stores the last thirty measurements in memory automatically. The data is retained as long as the batteries are in the device. When you remove the batteries or the batteries are drained, the data will be erased.

1. Press and hold the START button until the memory data appears on the display.
2. The newest data appears. First, the Memory Number, then its data.
3. The memorized data appears automatically one after the other, from the most recent to oldest (i.e. 30, 29, 28 ... 3, 2, 1).

After the last data appears, the display goes blank.

With nothing displayed, press and hold the START button.

Note)

Pressing the START button while the memory data is displayed will cause the UB-401 to start

taking a measurement.

About Blood Pressure

1. What is blood pressure?

Blood pressure is the force exerted by blood against the walls of the arteries. Systolic pressure occurs when the heart contracts. Diastolic pressure occurs when the heart expands.

Blood pressure is measured in millimeters of mercury (mmHg). One's natural blood pressure is represented by the fundamental pressure, which is measured first thing in the morning while one is still at rest and before taking food or drink.

2. What is hypertension and how is it controlled?

Hypertension, an abnormally high arterial blood pressure, if left unattended, can cause many health problems including stroke and heart attack.

Hypertension can be controlled by altering lifestyle, avoiding stress and with medication under a doctor's supervision.

To prevent hypertension or keep it under control:

- Do not smoke.
- Reduce salt and fat intake.
- Maintain proper weight.
- Exercise regularly.
- Have regular physical checkups.

3. Why measure blood pressure at home?

Blood pressure measured at a clinic or doctor's office may cause apprehension and can produce an elevated reading, 25 to 30 mmHg higher than that measured at home. Home measurement reduces the effects of outside influences on blood pressure readings, supplements the doctor's readings and provides a more accurate, complete blood pressure history.

4. Blood pressure measured on your wrist may differ from that measured on your upper arm.

Blood pressure measured on your wrist and on your upper arm may be different. This results from the difference in widths of the blood vessels in your wrist and arm. The average difference in values for a healthy adult is about ± 10 mmHg.

The arteries in the upper arm lead to the wrist. Blood pressure variations in the arteries are reflected in the blood pressure in the wrist. So, measuring blood pressure on your wrist is a good way to check your blood pressure.

Note)

With those who suffer from hypertension, arteriosclerosis, diabetes, kidney


failure, and other disorders of the circulatory system, the difference may be greater.

About Blood Pressure

5. WHO blood pressure classification

Standards for assessment of high blood pressure, without regard to age, have been established by the World Health Organization (WHO), as shown in the chart below.

Reference Material:
Journal of Hypertension
1999, Vol 17 No.2


6. Blood pressure variations

An individual's blood pressure varies greatly on daily and seasonal bases. It may vary by 30 to 50 mmHg due to various conditions in a day. In hypertensive individuals, variations are even more pronounced.

Normally, the blood pressure rises while at work or play and falls to its lowest levels during sleep. So, do not be overly concerned by the results of one measurement.

Take measurements at the same time every day using the procedure described in this manual, and know your normal blood pressure. Many readings give a more comprehensive blood pressure history.

Be sure to note date and time when recording your blood pressure. Consult your doctor to interpret your blood pressure data.

Troubleshooting

Problem	Probable cause	Corrective action
Nothing appears in the display, even when the power is turned on.	Batteries are drained. Battery terminals are not in the correct position.	Replace all batteries with new ones. Place the batteries with negative and positive terminals matching those indicated on the cuff.
The cuff will not inflate.	Battery voltage is too low. ( , LOW BATTERY mark, appears. If the batteries are drained, the mark does not appear.)	Replace all batteries with new ones.
The device will not measure. Readings are too high or too low.	The cuff is not fastened properly. You moved your arm or body during the measurement. The cuff position is not correct. _____	Fasten the cuff correctly. Make sure you remain very still and quiet during the measurement. Sit comfortably and still. Raise your hand so that the cuff is at the same level as your heart. If you have a very weak or irregular heart beat, the device may have difficulty in determining your blood pressure.
Others	The value is different from that measured at a clinic or doctor's office. _____	See "3. Why measure blood pressure at home" on page 8. Remove the batteries. Place them back in correctly, and try the measurement again.

Note)

If the actions described above do not solve the problem, contact the dealer.
Do not attempt to repair the device yourself.

Maintenance

Do not open the device. It uses delicate electrical components and an intricate air unit that could be damaged. If you can not fix the problem using the troubleshooting instructions, request service from your supplier or from the A&D service group. The A&D service group will provide authorized suppliers technical information, spare parts and units.

The technical testing procedures, which should be done at least every 3 years, can be performed either by the manufacturer or by an authorized repair service in accordance with the regulations governing manufacturing of medical products.

Technical Data

Type	UB-401
Measurement method	Oscillometric measurement
Measurement range	Pressure : 20 - 280 mmHg Pulse : 40 - 200 beats / minute
Measurement accuracy	Pressure : ± 3 mmHg or 2%, whichever is greater Pulse : $\pm 5\%$
Power supply	2 x 1.5V alkaline batteries (LR03/Micro/AAA)
Wrist circumference	13.5 - 21.5 cm ($5\frac{1}{4}$ - $8\frac{1}{2}$ in.)
Classification	Type BF 
Clinical test	According to ANSI / AAMI SP-10 1987
Operating condition	?10 °C to ?40 °C / 30%RH to 85 %RH
Storage condition	?10 °C to ?60 °C / 30%RH to 85 %RH
Dimensions	Approx. 74 [W] x 64 [H] x 23 [D] mm, excluding cuff
Weight	Approx. 93 g, excluding batteries

 0366

Note)

Specifications are subject to change without prior notice.

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A&D Co., LTD.
A&D INSTRUMENTS LTD.
A&D MERCURY PTY. LTD.