Methods of connecting multiple balances/scales to a PC

Connection method	Number of connectable units *1	Connection equipment *2		Software	Remarks
Connect Serial/USB converter to PC USB port	4	Serial/USB converter	AX-USB-9P AX-USB-25P AX-USB-DIN		WinCT-Plus can display weighing values from multiple
Connect Bluetooth converter to PC USB port	4	AD8529PC-W Bluetooth converter		winCT-Plus or WinCT (RS-Key does not work)	balances/scales in one screen. For WinCT, launch separate WinCT windows for the number of connections and set the Com port for each window.
Connect USB hub to PC USB port and connect Bluetooth converter to hub *3	6	AD8529PC-W Bluetooth converter + USB hub (commercially available)			
Connect serial communication device to PC USB port *4	4	Serial communication device (commercially available)			
Connect via Ethernet	100 *5	Optional Ethernet interface or AD-8526 Serial/Ethernet converter		WinCT-Plus	Ethernet cable and ethernet hub are required

^{*1} Number of units with which operation was confirmed

In any of the connection methods, when a large amount of data is transmitted, such as through stream mode, there is a possibility that the data may be lost, so connection is limited to command mode, auto print mode, etc. (This is a PC buffer problem that depends on PC specs)

^{*2} We do not guarantee operation for non-A&D products / Depending on the specifications of the balance/scale, optional items such as connection cable and AC adapter are required.

^{*3} Operation confirmed with Buffalo BSH7 AE03 Bus power / Self-power (Operation possible but not guaranteed).

^{*4} Operation confirmed with CONTEC COM4CX-USB (Operation possible but not guaranteed).

^{*5} AD-8526 Serial/Ethernet Converter specifications allow up to 100.