

# UA-1000 SERIES

UA-1010 / UA-1020 / UA-1030T



## Blood Pressure Monitors



Advanced Cuff Design  
Advanced Inflation  
Advanced Deflation



JAPAN  
Quality & Design



Technology You Can Trust

A&D Medical  
<http://www.aandd.jp>

# UA-1000 SERIES

## UA-1010 / UA-1020 / UA-1030T

A&D, a pioneer in home blood pressure monitoring, has once again come up with a new series of world-class monitors with 3A Technology that sets new industry standards and enhances value for users.

SlimFit cuff



### UA-1010



**Simple monitoring with more accuracy** It offers easy, one-button operation with new functions that will guide user to accurate measurements. UA-1010 uses the new latex-free and metal-free SlimFit cuff which is more friendly to sensitive skin.

### UA-1020



**Individualized monitoring with more comfort** UA-1020 is equipped with newly designed SmoothFit cuff, which enables painless, comfortable measurement. Also it has two unique modes. MyPressure mode provides individualized inflation and shorten the measurement time. The TriCheck™ mode automatically takes three consecutive readings with one-minute intervals for most accurate results.



SmoothFit cuff  
\*Patent pending



SmoothFit cuff  
\*Patent pending

### UA-1030T



**Talking monitor for everyone** UA-1030T has come with the latest talking technology with all features of UA-1020. It announces the results in 5 languages (English, French, Spanish, Italian, and Chinese). You can increase volume through six levels or switch on to silent mode if required.



Cuff Fit Error Indicator



Movement Error Indicator



Irregular Heartbeat Indicator (IHB)



Latex-free and metal-free, SlimFit Cuff



WHO Blood Pressure Classification Indicator



M 90 90 memories + Average Reading (All/AM/PM)



Clock & Calendar



Newly designed SmoothFit Cuff for painless and comfortable measurement (Latex-free and metal-free)



My Pressure Setting mode (AUTO/180/210/240)



TriCheck™ mode for consecutive 3 measurements and average reading



Clear voice announcement of measurement results



**Advanced Cuff Design**  
**Advanced Inflation**  
**Advanced Deflation**

Latex-free and metal-free plus wrinkle-free cuff designs.

User selection made possible between A&D's fuzzy logic inflation and anxiety-free "MyPressure" setting.

Linear-motor driven active deflation valve makes it measure accurately for all cuff sizes.



#### A&D Company, Limited

3-23-14 Higashi-Ikebukuro, Toshima-ku, Tokyo 170-0013 Japan  
Telephone: [81](3)5391-6132 Fax: [81](3)5391-6148  
1-243, Asahi, Kitamoto-shi, Saitama, 364-8585 Japan  
Telephone: [81](0485)93-1111 Fax: [81](0485)93-1119  
(ISO13485 Certified)  
<http://www.aandd.jp>

#### A&D INSTRUMENTS LTD.

Unit 24/26 Blacklands Way Abingdon Business Park,  
Abingdon, Oxon OX14 1DY UNITED KINGDOM  
Telephone: [44](1235) 550420 Fax: [44](1235) 550485  
(Authorized Representative Established in the European Community)