



AUTOMATIC Blood Pressure Monitor

- Irregular HeartBeat (IHB) Indicator
- Blood Pressure Classification Indicator
- SlimFit Comfortable Cuff
- Fully Automatic Measurement
- Large 30-Memory + Average Reading
- Standard AC adapter Inlet

UA-621

* Asia/EC model



WRIST Blood Pressure Monitor

- Irregular HeartBeat (IHB) Indicator
- Blood Pressure Classification Indicator
- Comfortable Wrist Cuff
- Fast Measurement along with Inflation
- Large 30-Memory + Average Reading

UB-510

* Asia/EC model



A&D
A&D Medical
Technology You Can Trust
<http://www.aand.jp>



A&D
A&D Medical
Technology You Can Trust
<http://www.aand.jp>

UA-621 Upper Arm Automatic Blood Pressure Monitor

UB-510 Wrist Blood Pressure Monitor

*Asia/EC model



Irregular HeartBeat Indicator

The indicator illuminates if an IHB is detected during the measurement. The indicator alerts users of an irregular heartbeat, which may be associated with cardiovascular disease. (Note: This is for warning purposes only. A&D recommends contacting a physician if the user observes this indicator frequently.)

Blood Pressure Classification Indicator

The World Health Organization (WHO) has designated 6 blood pressure classifications. This monitor automatically indicates your classification (Optimal to Severe hypertension) after each measurement is completed.

Systolic Pressure	Diastolic Pressure	
180~	110~	Grade 3 hypertension (severe)
160~179	100~109	Grade 2 hypertension (moderate)
140~159	90~99	Grade 1 hypertension (mild)
130~139	85~89	High-normal
120~129	80~84	Normal
Less than 120	Less than 80	Optimal

Reference: Journal of Hypertension 1999, Vol 17 No.2

Multilingual Operation manual

Asia model		EC model		
English	Chinese	English	Malayalam	Arabic
Korean	Thai	Spanish	French	Chinese
Malay	Vietnamese	Italian	Portuguese	
Hindi	Indonesian	Turkey	German	

SlimFit Comfortable Cuff (UA-621)

During a measurement, our new innovative edge free internal bladder "EFB" distributes pressure evenly around the user's arm. This provides a higher level of comfort while the arm is under pressure.

Fully Automatic Measurement

Pressing the START button is the only step necessary to perform a complete and accurate measurement.

Large 30-Memory + Average Reading

Both the UA-621 and the UB-510 have a 30-reading memory capacity for monitoring your measurement trends. An average of all the stored data is shown at the start of each measurement, offering the user their comprehensive measurement history at a glance.

Standard AC adapter Inlet (UA-621)

The monitor is equipped with an AC inlet for use with an external power source. Since there is no need to worry about battery life, the user can economize while enjoying the use of the monitor everyday.

(The AC adapter is optional and should be bought separately.)

A&D
A&D Medical
Tokyo JAPAN
<http://www.aandd.jp>

Technology You Can Trust