

*Auto Inflation  
Blood Pressure Monitor  
Designed for Women*

UA-782



CE 0123

**A&D**  
A&D Medical

Technology You Can Trust

<http://www.aandd.jp>

# Blood Pressure Monitor

## *Designed for Women*

### Model UA-782

#### *Be feminine...Don't dismiss your curves*

Are you surprised by the number of medical products that dismiss your femininity? Women enjoy their femininity as a part of staying vibrant & healthy. The right diet, clothes, cosmetics, exercise, etc. are just a few areas that women monitor closely.

However, as you are aware, paying attention to your health is one of the most significant factors in maintaining your femininity at any age. We designed a blood pressure monitor to help women do just that!

With the UA-782 we placed latest technology in the uniquely shaped body of a digital BP monitor. Please take a moment to learn more...

#### *Smooth and fast measurements*

You may feel a little uncomfortable during blood pressure measurements because of the pressure applied to the arm and the noise while inflating the cuff. With our newest technology, those discomforts are a thing of the past. The UA-782 requires a shorter measurement time and its newly developed semi-hard shell cuff also caters to a wider range of arm sizes. The new whisper-silent pump further aids smooth and fast measurements. Arm range: 17cm-32cm

#### *Irregular Heart Beat (IHB) indicator*

This indicator illuminates if an IHB is detected during a measurement. The indicator alerts you of an irregular heartbeat, sometimes a symptom of cardiovascular disease. (Note: This is for warning purposes only. A&D recommends contacting a physician if the user observes this indicator frequently).

#### *Blood pressure classification indicator*

The World Health Organization (WHO) distinguishes 6 different blood pressure stages. This monitor automatically indicates your classification (from Optimal to Severe hypertension) after each measurement is completed.

#### *90-reading memory*

The UA-782 has a 90-reading memory capacity for monitoring your measurement trends. An average of all the stored data is shown when the device is started, offering the user their comprehensive measurement history at a glance.

#### *Newly designed e-Cuff*

A newly designed semi-hard shell cuff promises easy fastening around your arm as well as comfortable blood pressure measurements for patients in a wider arm range (17 to 32cm).

#### *Unique heart shaped cuff holder*

The UA-782 is equipped with a convenient cuff holder. Simply place the cuff in the cuff holder located on the backside of the unit for storage after the measurement is completed.

**A&D**  
A&D Medical

*Technology You Can Trust*

<http://www.aandd.jp>

\*UA782CA-ADCC-02-PR2-143\*\*  
REF\*CARD-UA782