

## Verified operation of our software with Windows Vista

Updated: May 2008

Description	Verified operation results with Windows Vista
WinCT Ver.3.00 <sup>i</sup>	<b>Failed:</b> RsKey <sup>ii</sup>
WinCT-Moisture	Verified: RsTemp, RsFig
WinCT-Viscosity	Verified: RSVisco
WinCT-Plus	Verified: RSMulti, DeviceInstaller <b>Failed:</b> Opening Excel from RSMulti
Quick USB	Verified: Automatic installation and data input to Excel 2007 (can be used as the alternative to RsKey)

---

<sup>i</sup> Windows Vista-compliant, Ver.4.01 (not including RsWeight) is now available for download.

<sup>ii</sup> RsKey only. The operations of RsCom and RsWeight were verified.