


# Update to accessory for HR-AZ / A and FX-CT series

Thank you for purchasing the HR-AZ/A or FX-CT series. This is to notify customers of a change to an accessory for the series and to the dimensions inside the breeze break.

1. Both a weighing pan and a pan support are now provided together as a single weighing pan unit.



Weighing pan and pan support (separate components)


Weighing pan unit (one component)

**Previous design**

**New design**

See 'Unpacking the Balance' and 'Precautions' in the instruction manual.

2. Accordingly, the height and inside dimension above the pan have been changed as shown below.



**Previous dimensions**

**New dimensions**

See 'External Dimensions' in the instruction manual.