

DIGITAL BLOOD PRESSURE MONITOR

UA-705

Clinically
Validated

IHB



Compact Desktop Monitor

- Multi-feature Results
- Irregular Heart Beat (IHB) Indicator
- Blood Pressure Classification Indicator
- 30 Memory + Average Reading
- SlimFit Comfortable Cuff
- Uses One AA Battery

CE 0123

A&D
A&D Medical

Technology You Can Trust

<http://www.aandd.jp>

DIGITAL BLOOD PRESSURE MONITOR UA-705

Multi-feature Results

This convenient compact desktop style houses fast measurement, 3-line display for simultaneous reading, IHB indicator, WHO classification indicator, comfortable cuff and long battery life – just press START and squeeze the rubber bulb a few times!

Irregular Heart Beat (IHB) Indicator

This lights up if IHB is detected during the measurement. Helps users to find an irregular heartbeat which may be associated with cardiovascular disease. (Note: This is for warning purposes only. A&D recommends contacting a physician if the user sees this indicator frequently.)

Blood Pressure Classification Indicator

World Health Organization (WHO) classifies blood pressure ranges into 6 stages. This monitor indicates your position (Optimal to Severe Hypertension) automatically after each measurement is complete.

30 Memory + Average Reading

The UA-705 has a large memory of 30 readings for monitoring the trend of measurements. An average of all the stored data is shown at each start up, giving a condensed reading history at a glance.

SlimFit Comfortable Cuff


A new innovative edge free internal bladder “EFB” promises even pressure on the arm during the measurement. This provides the user with a higher level of comfort while the arm is under pressure.

Uses One AA Battery

The monitor is designed to conserve energy, promising a long unit service life to keep your measurements going.

WHO blood pressure classification

Standards for assessment of high blood pressure, without regard to age, have been established by the World Health Organization (WHO), as shown in the chart on the right.



A&D
A&D Medical

Technology You Can Trust

<http://www.aandd.jp>

*UA705CA-ADCC-03-PR3-143**
REF*CARD-UA705