


Disposable arm covers for the blood pressure monitor. These single-use hygienic arm covers are easy to wear. Disposable arm covers offer cost effective protection to arms and garment sleeves.

Features

- ▶ For infection control
- ▶ Made of disposable plastic
- ▶ 400 pcs per carton (200 pcs / batch x 2)
- ▶ Accommodates arm circumferences up to 44 cm
- ▶ Translucent for visibility
- ▶ Each piece is perforated, making them easy to detach
- ▶ Attached strings make each batch simple to hang anywhere



How To Use



Specifications

Item	Description
Product Code	AS-134010367
Material	High density polyethylene
Quantity	400 pcs per carton (200 pcs x 2 batches)
Size	750(W) x 0.015(H) x 220(D) mm per piece
Thickness	15 µm