

Blood Pressure Monitor

Model : UA-1100NFC

General Specifications

- Easy to operate-Just place your arm through the cuff and press the button
- The upper-arm type enables measurement in a stable posture and level with heart height
- The rubber band attached to the cuff keeps the unit compact, making it easy to carry and store
- Easy-to-read screen despite its compact size
- Easy data management through the App and NFC communication function
- No pairing is required, so a single unit can be shared by multiple people, such as family members
- 3 color variations(White/Gray/Purple)



× 3
1.5V (LR03)

Easy to Use



1. Secure the cuff around the upper arm and press the START button.



2. Data can be recorded in the App after measurement.

Color Variations



White



Gray



Purple


Compact

The rubber band attached to the cuff keeps the unit compact.



Easy-to-read Screen

Displays systolic and diastolic blood pressure value and pulse rate simultaneously.



Specifications

Measurement Method	Oscillometric measurement
Display type	Digital LCD, Pressure / Pulse displayed simultaneously
Measurement range	Pressure: 0 - 299 mmHg Systolic pressure: 60 mmHg - 279 mmHg Diastolic pressure: 40 mmHg - 200 mmHg Pulse: 40 - 180 beats / minute
Measurement accuracy	Pressure: ± 3 mmHg / Pulse: ± 5 %
Power supply	3 x 1.5 V batteries (AAA)
Number of measurements	Approx. 400 times AAA (alkaline batteries) With pressure value 170 mmHg, room temperature 23 °C.
Memory Storage	Previous measurement
Cuff	Arm circumference: 22 cm - 32 cm
Operating Conditions	+10 °C to +40 °C / 15 %RH to 85 %RH / 800 hPa to 1060 hPa
Transport / Storage Conditions	-20 °C to +60 °C / 10 %RH to 95 %RH / 700 hPa to 1060 hPa
Dimensions (W) x (H) x (D)	Approx. 36 x 144 x 31 mm
Weight	Approx. 240 g, excluding batteries
Individual box dimensions (W) x (H) x (D)	71 x 180 x 71 mm
Individual box weight	Approx. 373 g, excluding local instruction manual
Master carton dimensions (W) x (H) x (D)	150 x 200 x 370 mm
Master carton weight	Approx. 3994 g, excluding local instruction manual

Download A&D Connect Healthcare



A&D Connect
Healthcare



https://link.aandd.jp/app_downloadlink_1



https://link.aandd.jp/app_downloadlink_2



Guide video of the App



https://link.aandd.jp/app_healthcare_asia

What is In The Box

- 1 Digital Blood pressure monitor
- 1 Instruction manual
- 3 Alkaline Batteries
- 1 Quick Start Guide

- App Store is a trademark of Apple Inc., registered in the U.S. and other countries.
- Google Play and the Google Play logo are trademarks of Google LLC.
- The N-Mark is a trademark or registered trademark of NFC Forum, Inc. in the United States and in other countries.



A&D Company Limited.(JAPAN)
URL:aandd.jp

A&D ENGINEERING, INC.(USA)
URL:andonline.com

A&D Instruments Canada Inc.(CANADA)
URL:autocontrol.com

A&D INSTRUMENTS LIMITED(United Kingdom)
URL:andprecision.com

A&D AUSTRALASIA PTY LTD(Australia)
URL:andastralasia.com.au

A&D RUS CO., LTD.(Russia)
URL:and-rus.ru

A&D Technology Trading(Shanghai)Co., Ltd(China)
URL:aanddtech.cn

A&D INSTRUMENTS INDIA PRIVATE LIMITED(India)
URL:aanddindia.in

A&D Instruments (Thailand) Limited
URL:thai.andprecision.com/

A&D KOREA Limited
URL:andk.co.kr/main/