


DIGITAL BLOOD PRESSURE MONITOR

UA-767 Automatic



Advanced ONE BUTTON OPERATION

- One button operation : automatic power on cuff inflation blood pressure measurement and cuff deflation.
- Easy to read, liquid crystal display : Systolic, diastolic and heart rate shown simultaneously.
- Advanced Oscillometric method for hospital accuracy.



A&D
A&D Medical
Technology You Can Trust
<http://www.aand.jp>

DIGITAL BLOOD PRESSURE MONITOR UA-767 Automatic

One Button Operation

The UA-767 has one button operation : Power-on, start, pressurization at pre-set level, automatic measurement, display of systolic and diastolic blood pressure readings and heart rate, automatic cuff deflation and unit shut off are all accomplished with the push of only one button.

Easy to Read, Liquid Crystal Display

Like its brother model, UA-777, this unique display simultaneously displays systolic /diastolic blood pressure and heart rate, giving the user instantaneous readings of all critical data.

Advanced Oscillometric Measurement Technique


The same Oscillometric technique used by physicians is incorporated into the UA-767. This technique, along with advanced software and an extra sensitive pressure sensor, assures reliable and accurate blood pressure measurement.

Ultra Modern Design

The UA-767 is the most attractive digital blood pressure monitor on the market. It can be used with pride in either home or office.

WHO blood pressure classification

Standards for assessment of high blood pressure, without regard to age, have been established by the World Health Organization (WHO), as shown in the chart on the right.



A&D
A&D Medical
Tokyo JAPAN
Technology You Can Trust
<http://www.aand.jp>