

Connection and Function List of AD-8922A Remote Controller

Series/Model	RS232C output interface (Option name if required)	Required Cable	ON:OFF	CAL	SAMPLE	PRINT	MODE	RE-ZERO	Command set (AD-8922A function settings)
BM	Standard (D-Sub9P)	AX-KO2741-180 or AX-KO2466-200	Yes	Yes	Yes	Yes	Yes	Yes	CSet 1
BA-T	Standard (D-Sub9P)	AX-KO2741-180 or AX-KO2466-200	Yes	Yes	-	Yes	-	Yes	CSet 1
BA	Standard (D-Sub9P)	AX-KO2741-180 or AX-KO2466-200	Yes	Yes	Yes	Yes	Yes	Yes	CSet 1
BH-T	Standard (D-Sub9P)	AX-KO2741-180 or AX-KO2466-200	Yes	Yes	-	Yes	-	Yes	CSet 1
BH	Standard (D-Sub9P)	AX-KO2741-180 or AX-KO2466-200	Yes	Yes	Yes	Yes	Yes	Yes	CSet 1
GH	Standard (D-Sub9P)	AX-KO2741-180 or AX-KO2466-200	Yes	Yes	Yes	Yes	Yes	Yes	CSet 1
GR	Standard (D-Sub9P)	AX-KO1710-200	Yes	Yes	Yes	Yes	Yes	Yes	CSet 2
HR	Standard (D-Sub9P)	AX-KO1710-200	Yes	-	Yes	Yes	Yes	Yes	CSet 3
HR-i	Standard (D-Sub9P)	AX-KO2741-180 or AX-KO2466-200	Yes	-	Yes	Yes	Yes	Yes	CSet 3
HR-A	Standard (D-Sub9P)	AX-KO2741-180 or AX-KO2466-200	Yes	-	Yes	Yes	Yes	Yes	CSet 3
HR-AZ	Standard (D-Sub9P)	AX-KO2741-180 or AX-KO2466-200	Yes	Yes	Yes	Yes	Yes	Yes	CSet 1
GF-A	Standard (D-Sub9P) or GX-A-03	AX-KO2741-180 or AX-KO2466-200	Yes	-	Yes	Yes	Yes	Yes	CSet 3
GX-A	Standard (D-Sub9P) or GX-A-03	AX-KO2741-180 or AX-KO2466-200	Yes	Yes	Yes	Yes	Yes	Yes	CSet 1
GF-AWP	Standard (D-Sub9P) or GX-A-03	AX-KO2737-500EX (Balance-side Dust/Waterproof) AX-KO2741-180 or AX-KO2466-200	Yes	-	Yes	Yes	Yes	Yes	CSet 3
GX-AWP	Standard (D-Sub9P) or GX-A-03	AX-KO2737-500EX (Balance-side Dust/Waterproof)	Yes	Yes	Yes	Yes	Yes	Yes	CSet 1
GF-M	Standard (D-Sub9P)	AX-KO2737-500EX (Balance-side Dust/Waterproof) AX-KO2741-180 or AX-KO2466-200	Yes	-	Yes	Yes	Yes	Yes	CSet 3
GX-M	Standard (D-Sub9P)	AX-KO2737-500EX (Balance-side Dust/Waterproof) AX-KO2741-180 or AX-KO2466-200	Yes	Yes	Yes	Yes	Yes	Yes	CSet 1
GF-L	Standard (D-Sub9P)	AX-KO2737-500EX (Balance-side Dust/Waterproof) AX-KO2741-180 or AX-KO2466-200	Yes	-	Yes	Yes	Yes	Yes	CSet 3
GX-L	Standard (D-Sub9P)	AX-KO2737-500EX (Balance-side Dust/Waterproof) AX-KO2741-180 or AX-KO2466-200	Yes	Yes	Yes	Yes	Yes	Yes	CSet 1
GF	Standard (D-Sub25P) or GX-04 (D-Sub25P)	AX-KO1710-200	Yes	-	Yes	Yes	Yes	Yes	CSet 3
GX	Standard (D-Sub25P) or GX-04 (D-Sub25P)	AX-KO1710-200	Yes	Yes	Yes	Yes	Yes	Yes	CSet 1
FX-i	Standard (D-Sub9P)	AX-KO2741-180 or AX-KO2466-200	Yes	-	Yes	Yes	Yes	Yes	CSet 3
FZ-I	Standard (D-Sub9P)	AX-KO2741-180 or AX-KO2466-200	Yes	Yes	Yes	Yes	Yes	Yes	CSet 1
FX-iWP	Standard (D-Sub9P)	AX-KO2737-500EX (Balance-side Dust/Waterproof) AX-KO2741-180 or AX-KO2466-200	Yes	-	Yes	Yes	Yes	Yes	CSet 3
FZ-iWP	Standard (D-Sub9P)	AX-KO2737-500EX (Balance-side Dust/Waterproof) AX-KO2741-180 or AX-KO2466-200	Yes	Yes	Yes	Yes	Yes	Yes	CSet 1
FX	Standard (D-Sub9P)	AX-KO2741-180 or AX-KO2466-200	Yes	-	Yes	Yes	Yes	Yes	CSet 3
FZ	Standard (D-Sub9P)	AX-KO2741-180 or AX-KO2466-200	Yes	Yes	Yes	Yes	Yes	Yes	CSet 1
FX-WP	Standard (D-Sub9P)	AX-KO2737-500EX (Balance-side Dust/Waterproof) AX-KO2741-180 or AX-KO2466-200	Yes	-	Yes	Yes	Yes	Yes	CSet 3
FZ-WP	Standard (D-Sub9P)	AX-KO2737-500EX (Balance-side Dust/Waterproof) AX-KO2741-180 or AX-KO2466-200	Yes	Yes	Yes	Yes	Yes	Yes	CSet 1
GF-K	Standard (D-Sub25P) or GX-04K (D-Sub25P)	GX-07K (Balance-side Dust/Waterproof) AX-KO1710-200	Yes	-	Yes	Yes	Yes	Yes	CSet 3
GX-K	Standard (D-Sub25P) or GX-04K (D-Sub25P)	GX-07K (Balance-side Dust/Waterproof) AX-KO1710-200	Yes	Yes	Yes	Yes	Yes	Yes	CSet 1
GP	Standard (D-Sub25P) or GP-04 (D-Sub25P)	AX-KO1710-200	Yes	Yes	Yes	Yes	Yes	Yes	CSet 1
MC	Standard (D-Sub9P), GP-04 (-100KS), GX-04K(-10K/30K), or GX-04(-1000/6100)	AX-KO1710-200	Yes	Yes	Yes	Yes	Yes	Yes	CSet 1
AD4212A	Standard (D-Sub25P) or AD4212A-04 (D-Sub25P)	AX-KO1710-200	Yes	-	Yes	Yes	Yes	Yes	CSet 3
AD4212B	Standard (D-Sub25P)	AX-KO1710-200	Yes	-	Yes	Yes	Yes	Yes	CSet 3
AD4212C	Standard (D-Sub9P)	Cable included with AD4212C	Yes	Yes	Yes	Yes	Yes	Yes	CSet 1
AD4212D	Standard (D-Sub9P)	Cable included with AD4212D	Yes	Yes	Yes	Yes	Yes	Yes	CSet 1
AD4212F	Standard (D-Sub9P)	AX-KO3590-200/500/1000	Yes	Yes	Yes	Yes	Yes	Yes	Cset 1
EJ	EJ-03 (D-Sub9P)	AX-KO2741-180 or AX-KO2466-200	-	-	-	Yes	Yes	Yes	CSet 4
EK-i /EW-i	Standard (D-Sub9P)	AX-KO2741-180 or AX-KO2466-200	-	-	-	Yes	Yes	Yes	CSet 4
EK-L	EJ-03 (D-Sub9P)	AX-KO2741-180 or AX-KO2466-200	-	-	-	Yes	Yes	Yes	CSet 4
FG	FG-23 (Din8P)	AX-KO1786-200	-	-	-	Yes	-	-	CSet 0
FS	FS-03 (Din7P)	AX-KO1786-200	-	-	-	Yes	-	-	CSet 0
FS-i	FS-03i (Terminal block)	AX-KO3285-320	-	-	-	Yes	-	-	CSet 0
HV-G / HW-G	Standard (Din8P)	AX-KO1786-200	-	-	-	Yes	-	-	CSet 0
HV-C / HW-C	HVW-03C (D-Sub9P)	AX-KO2741-180 or AX-KO2466-200	-	-	-	Yes	-	-	CSet 0
HV-CWP / HW-CWP	HVW-03CWP (Terminal block)	AX-KO3285-320	-	-	-	Yes	-	-	CSet 0
SC/SE	SCE-03 (Terminal block)	AX-KO3285-320	-	-	-	Yes	-	Yes	CSet 4
SW	SW-03 (Terminal block)	AX-KO3285-320	-	-	-	Yes	-	-	CSet 0
FC-i/Si	Standard (D-Sub9P)	AX-KO2741-180 or AX-KO2466-200	-	-	-	Yes	-	-	CSet 0
HC-i	HC-03i (Din8P)	AX-KO1786-200	-	-	-	Yes	-	-	CSet 0
HD	HD-03i (Din7P)	AX-KO1786-200	-	-	-	Yes	-	-	CSet 0
GC	Standard (D-Sub9P)	AX-KO2466-200	-	-	-	Yes	Yes	-	CSet 0
HT/HT-CL	Connection unavailable								
HT-120	Connection unavailable								
FG-CWP	Connection unavailable								
HL-i	Connection unavailable								
SJ-HS	Connection unavailable								
SK-WP/SK	Connection unavailable								
HV-CEP / HW-CEP	Connection unavailable								