

Windows 10*1 Supported Versions of Communication Tools Software (as of November 29, 2023)

*1 Verified with 64-bit versions only.

For Balances/Scales (Free Software)

		Balances/Scales			Balances with FRD Function	Balances with FRD Function	Weigh Modules AD-4212C Series
A&D WinCT	Software Program	RsCom	RsKey	RsWeight	WinCT-FRD	WinCT-Plus	WinCT-AD4212C
	Version	6.00	6.00	6.00		RsMulti	
Verified OS	Win10 1511	Unverified*2	Unverified*2	Unverified*2	Unverified	Yes	Yes
	Win10 1703	Unverified*2	Unverified*2	Unverified*2	Yes	Yes	Yes
	Win10 1709	Unverified*2	Unverified*2	Unverified*2	Yes	Yes	Yes
	Win10 1803	Unverified*2	Unverified*2	Unverified*2	Yes	Yes	Yes
	Win10 1903	Unverified*2	Unverified*2	Unverified*2	Yes	Yes	Yes
	Win10 22H2	Yes	Yes	Yes	Yes	Yes	Yes

*2 Verified on version 5.40

For Moisture Analyzers (Provided as standard for the MX-50 and MS-70. Also available as a separately-sold software.e)

		Moisture Analyzers			
A&D WinCT	Software Program	RsFig	RsTemp	RsCom	RsKey
	Version	2.41M	2.31M	2.50M	2.50M
Verified OS	Win10 1511	Yes	Yes	Yes	No*3
	Win10 1703	Yes	Yes	Yes	No*3
	Win10 1709	Yes	Yes	Yes	No*3
	Win10 1803	Yes	Yes	Yes	No*3
	Win10 1903	Yes	Yes	Yes	No*3
	Win10 22H2	Yes	Yes	Yes	No*3

*3 RsKey is supported by Windows XP and below. When using Windows Vista and above, please use the free version of RsKey in WinCT for Balances/Scales.

For Viscometers (Included with SV-A/SV Series)

		Viscometers		
A&D WinCT	Software Program	RsVisco	RsCom	RsKey
	Version	1.13V	2.50V	2.50V
Verified OS	Win10 1511	Yes	Yes	No*3
	Win10 1703	Yes	Yes	No*3
	Win10 1709	Yes	Yes	No*3
	Win10 1803	Yes	Yes	No*3
	Win10 1903	Yes	Yes	No*3
	Win10 22H2	Yes	Yes	No*3

*3 RsKey is supported by Windows XP and below. When using Windows Vista and above, please use the free version of RsKey in WinCT for Balances/Scales.

For Pipette Accuracy Testers (Included with Pipette Accuracy Testers)

		Pipette Accuracy Testers
A&D WinCT	Software Program	WinCT-Pipette
	Version	1.25
Verified OS	Win10 1511	Yes
	Win10 1703	Yes
	Win10 1709	Yes
	Win10 1803	Yes
	Win10 1903	Yes
	Win10 22H2	Yes