

# Windows 10<sup>\*1</sup> Supported Versions of Communication Tools Software (as of March 2, 2020)

\*1 Verified with 64-bit versions only.

Free Software						
Ver.	Balance/Scale			WinCT-FRD	Multiple Balances and/or Scales	AD-4212C Series
	WinCT		WinCT-Plus		WinCT-AD4212C	
	RsCom	RsKey	RsWeight	RsMulti		
Ver.	5.40	5.40	5.40	1.00	3.00	1.13
Win10 1511	Yes	Yes	Yes	Not Yet	Yes	Yes
Win10 1703	Yes	Yes	Yes	Yes	Yes	Yes
Win10 1709	Yes	Yes	Yes	Yes	Yes	Yes
Win10 1803	Yes	Yes	Yes	Yes	Yes	Yes
Win10 1903	Yes	Yes	Yes	Yes	Yes	Yes

Paid Software								
Ver.	Moisture Analyzer				Viscometer/Rheometer			Pipette Accuracy Tester
	WinCT-Moisture <sup>*2</sup>				WinCT-Viscosity			WinCT-Pipette <sup>*4</sup>
	RsFig	RsTemp	RsCom	RsKey	RsVisco	RsCom	RsKey	
Ver.	2.41M	2.31M	2.50M	2.50M	1.13V	2.50V	2.50V	1.26
Win10 1511	Yes	Yes	Yes	No <sup>*3</sup>	Yes	Yes	No <sup>*3</sup>	Yes
Win10 1703	Yes	Yes	Yes	No <sup>*3</sup>	Yes	Yes	No <sup>*3</sup>	Yes
Win10 1709	Yes	Yes	Yes	No <sup>*3</sup>	Yes	Yes	No <sup>*3</sup>	Yes
Win10 1803	Yes	Yes	Yes	No <sup>*3</sup>	Yes	Yes	No <sup>*3</sup>	Yes
Win10 1903	Yes	Yes	Yes	No <sup>*3</sup>	Yes	Yes	No <sup>*3</sup>	Yes

\*2 Provided as standard for the MX-50 and MS-70

\*3 WinCT-Moisture/Viscosity RsKey is supported by Windows XP and below. When using Windows Vista and above, please use the free version of RsKey in WinCT for Balance/Scale.

\*4 Included with Pipette Accuracy Testers