

Windows 11*¹ Supported Versions of Communication Tools Software (as of October 20, 2021)

*1 Verified with 64-bit versions only.

For Balances/Scales (Free Software)

		Balances/Scales			Balances with FRD Function	Multiple Balances and/or Scales	Weigh Modules AD-4212C Series	Balances/Scales with UFC Function
A&D WinCT	Software	WinCT			WinCT-FRD	WinCT-Plus	WinCT-AD4212C	WinCT-UFC
	Program	RsCom	RsKey	RsWeight				
	Version	5.40	5.40	5.40				
Verified OS	Win11 21H2	Yes	Yes	Yes	Yes	Yes	Yes	Yes

For Moisture Analyzers (Provided as standard for the MX-50 and MS-70. Also available as a separately-sold software)

		Moisture Analyzers			
A&D WinCT	Software	WinCT-Moisture			
	Program	RsFig	RsTemp	RsCom	RsKey
	Version	2.41M	2.31M	2.50M	2.50M
Verified OS	Win11 21H2	Yes	Yes	Yes	No* ²

*2 RsKey is supported by Windows XP and below. When using Windows Vista and above, please use the free version of RsKey in WinCT for Balances/Scales.

For Viscometers (Included with SV-A/SV Series)

		Viscometers		
A&D WinCT	Software	WinCT-Viscosity		
	Program	RsVisco	RsCom	RsKey
	Version	1.13V	2.50V	2.50V
Verified OS	Win11 21H2	Yes	Yes	No* ²

*2 RsKey is supported by Windows XP and below. When using Windows Vista and above, please use the free version of RsKey in WinCT for Balances/Scales.

For Pipette Accuracy Testers (Included with Pipette Accuracy Testers)

		Pipette Accuracy Testers
A&D WinCT	Software	WinCT-Pipette
	Version	1.26
Verified OS	Win11 21H2	Yes