


A&D SeCUR stands for A&D 'Safe, Electronic Controlling Utility for Recording'.

Secure Your Weighing Data — Simplify Compliance with A&D SeCUR

What is A&D SeCUR?

A&D SeCUR is a lightweight, Windows-based application that enables secure, compliant, and traceable weighing data management from an A&D electronic balance. Designed for regulated industries, it helps you meet FDA 21 CFR Part 11 and ALCOA+ requirements with minimal setup.



Note) Supports one A&D balance per PC running A&D SeCUR.

What You Can Achieve with A&D SeCUR

- ✓ **Compliant Weighing:** Capture and store weighing data with full metadata
- ✓ **Audit-Ready Records:** Output uneditable PDF reports and complete audit trails
- ✓ **User Accountability:** Control user access and permissions with IDs and passwords
- ✓ **Tamper-Proof Data:** Protect authenticity with MD5 checksums and encryption
- ✓ **Traceable Corrections:** Retain invalidated weighing results with reasons

(Continued on the other page)

Home Screen

Weighing Screen

Audit Trail

PDF Report

Supported A&D Balance Models (Verified Compatibility)

- BA-T/BA series
- BH-T/BH series
- BM series
- GH series
- GR series
- GX-A/GF-A series
- GX-AWP/GF-AWP series
- GX-M/GF-M series
- GX-L/GF-L series
- FZ/FX series
- FZ-WP/FX-WP series
- EK-i series

Recommended PC Environment

- Windows 10 Pro (1903+) or Windows 11 Pro
- Display resolution: 1920×1080 (FHD), scaling 100%
- RS-232C, USB, or Ethernet port

Try A&D SeCUR Today!

Free 7-day trial — no license key needed.
To continue use, purchase a license key (USB dongle).

Download: link.aandd.jp/SeCUR



Discover Precision

A&D Company, Ltd. (JAPAN)
URL: aandd.jp

A&D Engineering, Inc. (USA)
URL: andonline.com

A&D Australasia Pty Ltd. (Australia)
URL: andaustralasia.com.au

A&D Instruments Ltd. (United Kingdom)
URL: andprecision.com

A&D Korea Ltd. (South Korea)
URL: andk.co.kr

A&D Rus Co., Ltd. (Russia)
URL: and-rus.ru

A&D Instruments India (P) Ltd. (India)
URL: aanddindia.in

A&D Sciencetech Taiwan Ltd. (Taiwan)
URL: aandd.com.tw

A&D Instruments Thailand Ltd. (Thailand)
URL: thai.andprecision.com

A&D Technology Trading (Shanghai) Co., Ltd. (China)
URL: aanddtech.cn