
Previous and New Icons Table



| New Icons | Descriptions | Previous Icons |
|----------------------------|---|-------------------|
| Comparator Function | Comparator Function Compares weighing results to a reference weight | COMPARATOR |
| -FOND- Auto Power ON | Auto Power ON Automatically turns on the display when plugged in | - €ON ▶ |
| -EOFF | Auto Power OFF Automatically turns off the display after a period of inactivity | ◆OFF |
| Underhook | Underhook Weighing of thermic or magnetic samples | under hook |
| Rechargeable Battery | Rechargeable Battery Works with rechargeable batteries | SLA |
| Normal Battery | Normal Battery Works with normal batteries | |
| Smart Range | Smart Range Automatically switches between the precision range and standard range | SMART RANGE |
| Universal Flexi Coms | Universal Flexi Coms Enables free editing of serial data output format or printout format | |
| EC Type Approval | EC Type Approval Approved in accordacne with the directives for "use in trade" of the EC | [M] |
| IP65 | IP65 IP65-level protection against water and dust | IP65 |
| NSF Compliant | NSF Compliant Meets sanitary guidelines as established by NSF International | NSF |
| Software | Software Special software available | WinCT |

