

## A&D at MEDICA 2009

**MEDICA 2009**  
Düsseldorf Nov. 18 - 21, 2009  
International Trade Fair with Congress

Hall 11  
Stand D78



As in previous years, A&D will be at booth 11D78. Our main theme is "Bridging Home and Hospital" and we will present three product categories (Home Health Care, Hospital, and Telemedicine) to show off our complete line of products and technology.

### Telemedicine

This year A&D will focus on Telemedicine and hold a large demonstration of mobile applications. We'll present [our Bluetooth® devices](#), [the first Continua Certified™ blood pressure monitor and weight scale](#), which will transfer measurement results to 4 mobile phones, each using a different application. Next, the phones will forward the data to a healthcare platform like Google Health or HealthVault and you will see the results on a PC screen. We'll also demonstrate transferring data to a Panasonic Toughbook, a PC designed for healthcare professionals.



### Home Health Care

A&D will exhibit a NEW product that will help with your proactive health care.

### Hospital

A&D is the leader of 24-hour ambulatory blood pressure monitoring. We are proud to present our [ABPM](#), which is used for hypertension research both in Japan and abroad. You can also try our [environmentally- and user-friendly sphygmomanometer](#) here.

Please visit us at booth 11D78 for the latest information.