

**CeBIT 2010**  
**(March 2-6, 2010, in Hannover, Germany)**

A&D Medical's Telemedicine teams from the UK, Germany, the USA, and Japan showed off their Bluetooth and USB/Continua certified products at CeBIT (March 2-6, 2010, in Hannover, Germany) at booth A08 in Hall 8. Displays included the UA-767PBT-C and UA-767PC blood pressure monitors, and the UC-321PBT-C, UC-321PC, and AD-6121APBT weight scales.

For more details about CeBit, please visit the following site.

<http://www.cebit.de/7590>

We demonstrated the coupling of different sets of linking devices with repositories such as Google Health and HealthVault, with the Continua Certified UA-767PBT-C and UC-321PBT-C communicating with a PC, along with an Android phone supported by LNI, USA. The EC-type approved AD-6121APBT medical scale prototype wirelessly connected with a Sony Vaio PC running Bluetooth Stack for Windows by Toshiba. The transferred weight data was displayed using Digital Note2, an application from Intel Japan.

We invited Dr. Gerold Porenta from Medmobile, Austria, for his Nokia E50 mobile phone application with UA-767PBT. Also invited was Dr. Günter Schreier of the Austrian Institute of Technology, who gave CeBIT visitors a presentation about his work using the NFC version of the UA-767PBT for CHD.

During the exhibition week, about 250 people measured their blood pressure with the UA-767PBT and saw the result on the Web at [www.medmobile.at](http://www.medmobile.at).

At a near-by exhibition stand, Robert Bosch Medical demonstrated a home health console with the UA-767PC.

CeBIT 2011 will be held from Tuesday, March 1 to Saturday, March 5, 2011. See you there!

A&D Medical UK team will be at WoHIT2010 (World of Health IT 2010 Conference & Exhibition) in Barcelona, Spain from March 15 to 18, 2010.

